

### Appendix M. MAUDE and Medical Device Recall Reports

This appendix gives results from the US Food and Drug Administration (FDA) Manufacturer and User Facility Device Experience (MAUDE)<sup>2</sup> and Medical Device Recalls<sup>1</sup> databases, in Table M1 through Table M3.

Table M1. US Food and Drug Administration Medical Device Recall Reports<sup>1</sup>

| Recall Number | Product Description   | Trade Name                                      | Recall Class | Center | Center Classification Date | Termination | Firm Name               | Manufacturer Recall Reason  |
|---------------|---|---|--------------|--------|----------------------------|-------------|-------------------------|---|
| Z-1805-2020   | Obalon Balloon System with Model 4300 Touch Dispenser Touchscreen; Software version 02.00.00.0076 | Obalon Balloon System with 4300 Touch Dispenser | 2            | CDRH   | 4/28/2020                  | NR          | Obalon Therapeutics Inc | During initial set-up and/or replacement of the dispenser batteries, the touch dispenser touchscreen can inadvertently lose calibration. This issue can also occur if the touchscreen is pressed when the device is powered on in preparation for a balloon administration. As a result, the touchscreen can become unresponsive and balloon inflation cannot be initiated. |

Table M2. Reports on Intra-gastric Balloon Devices from the FDA MAUDE Database<sup>2</sup>

Note: Key information is presented in this table, but detailed event descriptions are available at the [FDA MAUDE database](#).

| Report Number                 | Event Date | Event Type  | Manufacturer                  | Date Received | Product Code | Brand Name                                      | Device Problem  | Patient Problem   |
|-------------------------------|------------|-------------|-------------------------------|---------------|--------------|---|---|---|
| 30050998<br>03-2024-<br>00169 | 1/3/2024   | Injury      | Apollo<br>Endosurgery,<br>Inc | 1/31/2024     | LTI          | Orbera365<br>Intra-gastric<br>Balloon<br>System | Use of Device<br>Problem;<br>Migration  | Abdominal Pain;<br>Vomiting;<br>Obstruction/<br>Occlusion   |
| 30050998<br>03-2024-<br>00078 | 1/5/2024   | Injury      | Apollo<br>Endosurgery,<br>Inc | 1/31/2024     | LTI          | Orbera365<br>Intra-gastric<br>Balloon<br>System | Deflation<br>Problem; Use of<br>Device Problem  | No Clinical Signs,<br>Symptoms or<br>Conditions   |
| 30050998<br>03-2024-<br>00251 | 1/7/2024   | Injury      | Apollo<br>Endosurgery,<br>Inc | 1/31/2024     | LTI          | Orbera365<br>Intra-gastric<br>Balloon<br>System | Use of Device<br>Problem; Adverse<br>Event Without<br>Identified Device<br>or Use Problem | Abdominal Pain;<br>Cardiac Arrest;<br>Perforation;<br>Pneumonia; Sepsis;<br>Post Operative<br>Wound Infection |
| 30050998<br>03-2024-<br>00170 | 10/23/2023 | Injury      | Apollo<br>Endosurgery,<br>Inc | 1/26/2024     | LTI          | BIB<br>Intra-gastric<br>Balloon<br>System       | Inflation Problem;<br>Difficult to<br>Remove; Use of<br>Device Problem                    | Pain  |
| 30050998<br>03-2024-<br>00175 | 11/9/2023  | Malfunction | Apollo<br>Endosurgery         | 1/26/2024     | LTI          | BIB<br>Intra-gastric<br>Balloon<br>System       | Material<br>Separation;<br>Difficult to<br>Advance  | No Clinical Signs,<br>Symptoms or<br>Conditions   |
| 30126389<br>28-2024-<br>02991 | 12/6/2023  | Malfunction | Spatz Fgia<br>Inc.            | 1/22/2024     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System       | Patient-Device<br>Incompatibility;<br>Air/ Gas in<br>Device                               | Abdominal Pain;<br>Vomiting   |
| 30126389<br>28-2024-<br>02989 | 1/3/2024   | Malfunction | Spatz Fgia<br>Inc.            | 1/22/2024     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System       | Patient-Device<br>Incompatibility;<br>Air/ Gas in<br>Device                               | Insufficient<br>Information   |

| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                     | Device Problem   | Patient Problem                                 |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|--|--|---|
| 30050998<br>03-2023-<br>07189 | 12/26/2023 | Injury      | Apollo<br>Endosurgery,<br>Inc. | 1/19/2024     | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Burst Container<br>or Vessel; Use of<br>Device Problem | Vomiting  |
| 30050998<br>03-2024-<br>00038 | 8/15/2023  | Injury      | Apollo<br>Endosurgery,<br>Inc  | 1/17/2024     | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem;<br>Migration                     | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30050998<br>03-2024-<br>00062 | 11/20/2023 | Injury      | Apollo<br>Endosurgery,<br>Inc  | 1/17/2024     | LTI          | BIB<br>Intragastric<br>Balloon<br>System       | Inflation Problem                                      | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30126389<br>28-2024-<br>02987 | 12/5/2023  | Malfunction | Spatz Fgia<br>Inc.             | 1/17/2024     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem                                   | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30050998<br>03-2023-<br>07171 | 7/31/2023  | Injury      | Apollo<br>Endosurgery,<br>Inc  | 1/16/2024     | LTI          | BIB<br>Intragastric<br>Balloon<br>System       | Inflation Problem                                      | Nausea; Pain                                    |
| 30050998<br>03-2024-<br>00037 | 9/4/2023   | Injury      | Apollo<br>Endosurgery,<br>Inc  | 1/16/2024     | LTI          | BIB<br>Intragastric<br>Balloon<br>System       | Deflation<br>Problem; Use of<br>Device Problem         | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30050998<br>03-2024-<br>00043 | 10/28/2023 | Injury      | Apollo<br>Endosurgery,<br>Inc  | 1/16/2024     | LTI          | BIB<br>Intragastric<br>Balloon<br>System       | Inflation Problem                                      | Nausea; Pain;<br>Vomiting                       |
| 30050998<br>03-2024-<br>00061 | 10/31/2023 | Injury      | Apollo<br>Endosurgery,<br>Inc  | 1/16/2024     | LTI          | BIB<br>Intragastric<br>Balloon<br>System       | Inflation Problem                                      | Abdominal Pain;<br>Nausea                       |
| 30050998<br>03-2024-<br>00028 | 11/7/2023  | Injury      | Apollo<br>Endosurgery,<br>Inc. | 1/16/2024     | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem; Use of<br>Device Problem         | No Clinical Signs,<br>Symptoms or<br>Conditions |

| Report Number                 | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                               |
|-------------------------------|------------|-------------|-------------------------|---------------|--------------|--|--|---|
| 30050998<br>03-2023-<br>07168 | 12/9/2023  | Injury      | Apollo Endosurgery      | 1/16/2024     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Use of Device Problem               | No Clinical Signs, Symptoms or Conditions     |
| 30050998<br>03-2023-<br>07169 | 12/17/2023 | Malfunction | Apollo Endosurgery      | 1/16/2024     | LTI          | BIB Intra-gastric Balloon System       | Material Integrity Problem                             | No Clinical Signs, Symptoms or Conditions     |
| 30126389<br>28-2024-<br>02990 | 1/3/2024   | Malfunction | Spatz Fgia Inc.         | 1/16/2024     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions     |
| 30126389<br>28-2023-<br>02983 | 12/4/2023  | Malfunction | Spatz Fgia Inc.         | 1/15/2024     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions     |
| 30050998<br>03-2023-<br>07151 | 8/1/2023   | Injury      | Apollo Endosurgery      | 1/12/2024     | LTI          | BIB Intra-gastric Balloon System       | Inflation Problem; Use of Device Problem               | Abdominal Pain; Vomiting                      |
| 30126389<br>28-2023-<br>02977 | 10/5/2023  | Malfunction | Spatz Fgia Inc.         | 1/11/2024     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions     |
| 30126389<br>28-2023-<br>02980 | 11/29/2023 | Malfunction | Spatz Fgia Inc.         | 1/11/2024     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions     |
| 30050998<br>03-2023-<br>07131 | 9/30/2023  | Injury      | Apollo Endosurgery, Inc | 1/9/2024      | LTI          | Orbera365 Intra-gastric Balloon System | Adverse Event Without Identified Device or Use Problem | Erosion; Pyrosis/ Heartburn; Nausea; Vomiting |
| 30050998<br>03-2023-<br>07068 | 12/6/2023  | Malfunction | Apollo Endosurgery, Inc | 1/5/2024      | LTI          | Orbera365 Intra-gastric Balloon System | Material Separation; Difficult to Advance              | No Clinical Signs, Symptoms or Conditions     |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem                           | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30050998<br>03-2023-07056 | 12/9/2023  | Injury      | Apollo Endosurgery Inc., | 1/5/2024      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Use of Device Problem | No Clinical Signs, Symptoms or Conditions |
| 30050998<br>03-2023-07076 | 12/11/2023 | Malfunction | Apollo Endosurgery Inc., | 1/5/2024      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem                        | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02970 | 7/28/2023  | Malfunction | Spatz Fgia Inc.          | 1/4/2024      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                        | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02969 | 8/23/2023  | Malfunction | Spatz Fgia Inc.          | 1/4/2024      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                        | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02968 | 9/5/2023   | Malfunction | Spatz Fgia Inc.          | 1/4/2024      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                        | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02971 | 9/13/2023  | Malfunction | Spatz Fgia Inc.          | 1/4/2024      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                        | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02967 | 9/19/2023  | Malfunction | Spatz Fgia Inc.          | 1/4/2024      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                        | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02964 | 10/25/2023 | Malfunction | Spatz Fgia Inc.          | 1/4/2024      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                        | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02966 | 10/25/2023 | Malfunction | Spatz Fgia Inc.          | 1/4/2024      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                        | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer       | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2023-02965 | 11/14/2023 | Malfunction | Spatz Fgia Inc.    | 1/4/2024      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02962 | 11/28/2023 | Injury      | Spatz Fgia Inc.    | 1/4/2024      | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Patient Device Interaction Problem; Air/ Gas in Device | Abdominal Pain                            |
| 30050998<br>03-2023-07041 | 12/10/2023 | Malfunction | Apollo Endosurgery | 1/4/2024      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Migration   | No Clinical Signs, Symptoms or Conditions |
| 30050998<br>03-2023-07042 | 12/11/2023 | Injury      | Apollo Endosurgery | 1/4/2024      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Use of Device Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02961 | 10/25/2023 | Malfunction | Spatz Fgia Inc.    | 1/3/2024      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | Vomiting                                  |
| 30126389<br>28-2023-02959 | 9/27/2023  | Malfunction | Spatz Fgia Inc.    | 1/1/2024      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02956 | 10/9/2023  | Injury      | Spatz Fgia Inc.    | 12/28/2023    | LTI          | Spatz3 Adjustable Balloon System       | Improper or Incorrect Procedure or Method  | Laceration(s)                             |
| 30050998<br>03-2023-06938 | 11/27/2023 | Malfunction | Apollo Endosurgery | 12/28/2023    | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak  | No Clinical Signs, Symptoms or Conditions |

| Report Number                 | Event Date | Event Type  | Manufacturer          | Date Received | Product Code | Brand Name                                | Device Problem   | Patient Problem                                   |
|-------------------------------|------------|-------------|-----------------------|---------------|--------------|---|--|---|
| 30050998<br>03-2023-<br>06814 | 11/10/2023 | Injury      | Apollo<br>Endosurgery | 12/27/2023    | LTI          | BIB<br>Intragastric<br>Balloon<br>System  | Deflation<br>Problem; Use of<br>Device Problem;<br>Material Integrity<br>Problem | Vomiting  |
| 30050998<br>03-2023-<br>06813 | 8/11/2023  | Injury      | Apollo<br>Endosurgery | 12/26/2023    | LTI          | BIB<br>Intragastric<br>Balloon<br>System  | Inflation Problem;<br>Use of Device<br>Problem                                   | Inflammation;<br>Nausea; Vomiting;<br>Paresthesia |
| 30050998<br>03-2023-<br>06894 | 10/24/2023 | Malfunction | Apollo<br>Endosurgery | 12/26/2023    | LTI          | BIB<br>Intragastric<br>Balloon<br>System  | Material<br>Separation;<br>Detachment of<br>Device or Device<br>Component        | No Clinical Signs,<br>Symptoms or<br>Conditions   |
| 30126389<br>28-2023-<br>02955 | 11/23/2023 | Malfunction | Spatz Fgia<br>Inc.    | 12/26/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions   |
| 30126389<br>28-2023-<br>02954 | 11/22/2023 | Malfunction | Spatz Fgia<br>Inc.    | 12/24/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System | Patient-Device<br>Incompatibility;<br>Air/ Gas in<br>Device                      | Abdominal Pain;<br>Nausea; Vomiting               |
| 30126389<br>28-2023-<br>02951 | 11/15/2023 | Malfunction | Spatz Fgia<br>Inc.    | 12/21/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System | Patient-Device<br>Incompatibility;<br>Air/ Gas in<br>Device                      | Insufficient<br>Information                       |
| 30126389<br>28-2023-<br>02952 | 11/21/2023 | Malfunction | Spatz Fgia<br>Inc.    | 12/21/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions   |
| 30126389<br>28-2023-<br>02953 | 11/22/2023 | Malfunction | Spatz Fgia<br>Inc.    | 12/21/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System | Patient-Device<br>Incompatibility;<br>Air/ Gas in<br>Device                      | Insufficient<br>Information                       |

| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                     | Device Problem   | Patient Problem                                 |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2023-<br>02979 | 7/18/2023  | Malfunction | Spatz Fgia<br>Inc.             | 12/17/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30126389<br>28-2023-<br>02946 | 10/25/2023 | Malfunction | Spatz Fgia<br>Inc.             | 12/17/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30126389<br>28-2023-<br>02948 | 11/1/2023  | Malfunction | Spatz Fgia<br>Inc.             | 12/17/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30126389<br>28-2023-<br>02978 | 11/2/2023  | Malfunction | Spatz Fgia<br>Inc.             | 12/17/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30050998<br>03-2023-<br>06691 | 11/15/2023 | Malfunction | Apollo<br>Endosurgery<br>Inc., | 12/15/2023    | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem; Use of<br>Device Problem               | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30126389<br>28-2023-<br>02944 | 11/1/2023  | Malfunction | Spatz Fgia<br>Inc.             | 12/13/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30126389<br>28-2023-<br>02945 | 11/14/2023 | Malfunction | Spatz Fgia<br>Inc.             | 12/13/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30050998<br>03-2023-<br>06549 | 9/27/2023  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 12/11/2023    | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem; Use of<br>Device Problem;<br>Migration | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30126389<br>28-2023-<br>02942 | 10/6/2023  | Malfunction | Spatz Fgia<br>Inc.             | 12/11/2023    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions |



| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem                                     | Patient Problem  |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|--|
| 30126389<br>28-2023-02949 | 11/10/2023 | Malfunction | Spatz Fgia Inc.          | 12/11/2023    | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device | Abdominal Pain; Nausea; Vomiting                         |
| 30126389<br>28-2023-02950 | 11/15/2023 | Malfunction | Spatz Fgia Inc.          | 12/11/2023    | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device | Abdominal Pain; Eructate; Nausea                         |
| 30050998<br>03-2023-06590 | 10/15/2023 | Injury      | Apollo Endosurgery       | 12/8/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Use of Device Problem           | No Clinical Signs, Symptoms or Conditions                |
| 30050998<br>03-2023-06569 | 10/14/2023 | Malfunction | Apollo Endosurgery, Inc. | 12/7/2023     | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak; Inflation Problem               | No Clinical Signs, Symptoms or Conditions                |
| 30050998<br>03-2023-06589 | 11/6/2023  | Injury      | Apollo Endosurgery, Inc  | 12/7/2023     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem                                  | Pain; Discomfort   |
| 30050998<br>03-2023-06558 | 7/31/2023  | Injury      | Apollo Endosurgery, Inc  | 12/6/2023     | LTI          | BIB Intra gastric Balloon System       | Inflation Problem                                  | Nausea; Pain   |
| 30126389<br>28-2023-02938 | 11/8/2023  | Malfunction | Spatz Fgia Inc.          | 12/4/2023     | LTI          | Spatz3 Adjustable Balloon System       | Inflation Problem                                  | No Clinical Signs, Symptoms or Conditions                |
| 30067221<br>12-0230-00228 | NR         | Malfunction | Apollo Endosurgery, Inc  | 12/1/2023     | LTI          | Bib; Intra gastric Balloon System      | Inflation Problem                                  | Abdominal Pain; Nausea; Vomiting; Obstruction/ Occlusion |
| 30067221<br>12-2023-00226 | 3/9/2020   | Malfunction | Apollo Endosurgery, Inc  | 11/30/2023    | LTI          | Orbera® Intra gastric Balloon System   | Inflation Problem                                  | Abdominal Pain; Vomiting                                 |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|---|
| 30067221<br>12-2023-00227 | NR         | Malfunction | Apollo Endosurgery, Inc  | 11/30/2023    | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/Blood Leak; Unintended Deflation                                   | Failure of Implant  |
| 30126389<br>28-2023-02947 | 10/31/2023 | Malfunction | Spatz Fgia Inc.          | 11/28/2023    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions                 |
| 30067221<br>12-2023-00223 | 7/13/2023  | Malfunction | Apollo Endosurgery, Inc. | 11/27/2023    | LTI          | Orbera Intra-gastric Balloon           | Patient-Device Incompatibility  | Dehydration; Fatigue; Nausea; Vomiting                    |
| 30067221<br>12-2023-00224 | 10/25/2023 | Malfunction | Apollo Endosurgery, Inc. | 11/27/2023    | LTI          | Orbera Intra-gastric Balloon           | Patient-Device Incompatibility  | Fatigue; Vomiting   |
| 30067221<br>12-2023-00222 | 10/23/2023 | Injury      | Apollo Endosurgery, Inc. | 11/17/2023    | LTI          | Orbera 365 Intra-gastric Balloon       | Deflation Problem; Fluid/Blood Leak; Migration or Expulsion of Device; Unintended Deflation | Abdominal Pain; Failure of Implant; Obstruction/Occlusion |
| 30126389<br>28-2023-02941 | 11/5/2023  | Malfunction | Spatz Fgia Inc.          | 11/16/2023    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions                 |
| 30126389<br>28-2023-02935 | NR         | Malfunction | Spatz Fgia Inc.          | 11/16/2023    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration or Expulsion of Device   | Vomiting  |
| 30126389<br>28-2023-02931 | 10/11/2023 | Malfunction | Spatz Fgia Inc.          | 11/15/2023    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions                 |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                       | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|----------------------------------|--|---|
| 30126389<br>28-2023-02939 | 11/5/2023  | Malfunction | Spatz Fgia Inc.          | 11/12/2023    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02940 | 11/8/2023  | Malfunction | Spatz Fgia Inc.          | 11/12/2023    | LTI          | Spatz3 Adjustable Balloon System | Patient-Device Incompatibility; Air/ Gas in Device   | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00216 | 9/13/2023  | Malfunction | Apollo Endosurgery, Inc. | 11/10/2023    | LTI          | Orbera 365 Intra gastric Balloon | Deflation Problem; Fluid/ Blood Leak; Migration or Expulsion of Device; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2023-00221 | 10/10/2023 | Malfunction | Apollo Endosurgery, Inc. | 11/10/2023    | LTI          | Orbera 365 Intra gastric Balloon | Deflation Problem; Fluid/ Blood Leak; Migration or Expulsion of Device; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2023-00218 | 10/20/2023 | Malfunction | Apollo Endosurgery, Inc. | 11/10/2023    | LTI          | Orbera 365 Intra gastric Balloon | Deflation Problem; Fluid/ Blood Leak; Migration or Expulsion of Device; Unintended Deflation | Failure of Implant                        |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                       | Device Problem  | Patient Problem                                |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|----------------------------------|---|--|
| 30067221<br>12-2023-00215 | 9/21/2023  | Malfunction | Apollo Endosurgery, Inc. | 11/8/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Deflation Problem; Fluid/Blood Leak; Migration or Expulsion of Device; Unintended Deflation | Failure of Implant                             |
| 30126389<br>28-2023-02934 | 9/29/2023  | Malfunction | Spatz Fgia Inc.          | 11/8/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions      |
| 30067221<br>12-2023-00214 | 10/16/2023 | Malfunction | Apollo Endosurgery, Inc. | 11/8/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Deflation Problem; Fluid/Blood Leak; Migration or Expulsion of Device; Unintended Deflation | Failure of Implant                             |
| 30126389<br>28-2023-02936 | 10/25/2023 | Malfunction | Spatz Fgia Inc.          | 11/8/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions      |
| 30126389<br>28-2023-02932 | 10/20/2023 | Malfunction | Spatz Fgia Inc.          | 11/7/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions      |
| 30126389<br>28-2023-02933 | 10/20/2023 | Malfunction | Spatz Fgia Inc.          | 11/7/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions      |
| 30126389<br>28-2023-02930 | 9/17/2023  | Malfunction | Spatz Fgia Inc.          | 11/6/2023     | LTI          | Spatz3 Adjustable Balloon System | Patient-Device Incompatibility  | Abdominal Pain; Nausea; Vomiting; Constipation |

| Report Number                 | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem   |
|-------------------------------|------------|-------------|--------------------------|---------------|--------------|-------------------------------------|--|---|
| 30126389<br>28-2023-<br>02927 | 10/1/2023  | Injury      | Spatz Fgia Inc.          | 11/5/2023     | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility                             | Abdominal Pain; Pyrosis/ Heartburn; Vomiting                    |
| 30067221<br>12-2023-<br>00212 | 9/30/2023  | Malfunction | Apollo Endosurgery, Inc. | 11/3/2023     | LTI          | Orbera Intra gastric Balloon        | Patient-Device Incompatibility                             | Abdominal Pain; Dehydration; Nausea; Vomiting; Abdominal Cramps |
| 30126389<br>28-2023-<br>02924 | 9/15/2023  | Malfunction | Spatz Fgia Inc.          | 10/29/2023    | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | No Clinical Signs, Symptoms or Conditions                       |
| 30067221<br>12-2023-<br>00208 | 8/30/2023  | Malfunction | Apollo Endosurgery, Inc. | 10/27/2023    | LTI          | Orbera Intra gastric Balloon        | Inflation Problem; Air/ Gas in Device                      | Failure of Implant  |
| 30067221<br>12-2023-<br>00219 | 3/5/2023   | Malfunction | Apollo Endosurgery, Inc  | 10/26/2023    | LTI          | Orbera Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant  |
| 30067221<br>12-2023-<br>00220 | 7/25/2023  | Malfunction | Apollo Endosurgery, Inc. | 10/26/2023    | LTI          | Orbera 365 Intra gastric Balloon    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant  |
| 30126389<br>28-2023-<br>02923 | 9/27/2023  | Malfunction | Spatz Fgia Inc.          | 10/26/2023    | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | No Clinical Signs, Symptoms or Conditions                       |
| 30067221<br>12-2023-<br>00202 | 9/18/2023  | Malfunction | Apollo Endosurgery, Inc. | 10/25/2023    | LTI          | Orbera 365 Intra gastric Balloon    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant  |

| Report Number             | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|-------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2023-02919 | 9/19/2023  | Injury      | Spatz Fgia Inc.         | 10/25/2023    | LTI          | Spatz3 Adjustable Balloon System       | Improper or Incorrect Procedure or Method                  | Perforation of Esophagus                  |
| 30067221<br>12-2023-00205 | 9/18/2023  | Malfunction | Apollo Endosurgery, Inc | 10/23/2023    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2023-00210 | NR         | Malfunction | Apollo Endosurgery, Inc | 10/19/2023    | LTI          | Orbera Intra gastric Balloon System    | Insufficient Information                                   | Nausea; Vomiting; Obstruction/ Occlusion  |
| 30067221<br>12-2023-00209 | NR         | Malfunction | Apollo Endosurgery, Inc | 10/19/2023    | LTI          | Orbera365 Intra gastric Balloon System | Inflation Problem  | Dehydration; Nausea; Vomiting             |
| 30126389<br>28-2023-02929 | 9/14/2023  | Malfunction | Spatz Fgia Inc.         | 10/16/2023    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration or Expulsion of Device        | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00204 | 9/12/2023  | Malfunction | Apollo Endosurgery, Inc | 10/12/2023    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Leak/ Splash; Unintended Deflation      | Failure of Implant                        |
| 30067221<br>12-2023-00203 | 9/8/2023   | Malfunction | Apollo Endosurgery, Inc | 10/10/2023    | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2023-02917 | 9/14/2023  | Malfunction | Spatz Fgia Inc.         | 10/10/2023    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem  |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|-------------------------------------|--|--|
| 30067221<br>12-2023-00200 | 7/13/2023  | Malfunction | Apollo Endosurgery, Inc  | 10/9/2023     | LTI          | Orbera Intra-gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant   |
| 30067221<br>12-2023-00199 | 8/14/2023  | Malfunction | Apollo Endosurgery, Inc  | 10/9/2023     | LTI          | Orbera Intra-gastric Balloon System | Inflation Problem  | Failure of Implant   |
| 30067221<br>12-2023-00196 | 9/28/2023  | Malfunction | Apollo Endosurgery, Inc  | 10/9/2023     | LTI          | Orbera Intra-gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Bowel Perforation  |
| 30067221<br>12-2023-00201 | NR         | Malfunction | Apollo Endosurgery, Inc  | 10/9/2023     | LTI          | Orbera Intra-gastric Balloon System | Adverse Event Without Identified Device or Use Problem     | Cardiac Arrest; Nausea; Vomiting; Syncope/ Fainting                    |
| 30067221<br>12-2023-00193 | 9/2/2023   | Malfunction | Apollo Endosurgery, Inc. | 10/6/2023     | LTI          | Orbera Intra-gastric Balloon        | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant   |
| 30067221<br>12-2023-00189 | NR         | Malfunction | Apollo Endosurgery, Inc. | 10/6/2023     | LTI          | Orbera 365 Intra-gastric Balloon    | Adverse Event Without Identified Device or Use Problem     | Pyrosis/ Heartburn; Failure of Implant; Nausea; Vomiting               |
| 30126389<br>28-2023-02915 | 8/25/2023  | Malfunction | Spatz Fgia Inc.          | 10/5/2023     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | No Clinical Signs, Symptoms or Conditions                              |
| 30067221<br>12-2023-00186 | NR         | Malfunction | Apollo Endosurgery, Inc. | 10/3/2023     | LTI          | Orbera Intra-gastric Balloon        | Adverse Event Without Identified Device or Use Problem     | Pyrosis/ Heartburn; Failure of Implant; Inflammation; Nausea; Vomiting |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|-------------------------------------|--|---|
| 30067221<br>12-2023-00194 | 8/24/2023  | Malfunction | Apollo Endosurgery, Inc  | 10/2/2023     | LTI          | Orbera Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                             | Failure of Implant  |
| 30067221<br>12-2023-00190 | 9/5/2023   | Malfunction | Apollo Endosurgery, Inc. | 10/2/2023     | LTI          | Orbera Intra gastric Balloon        | Compatibility Problem  | Anemia; Gastritis; Failure of Implant; Low Oxygen Saturation; Gastrointestinal Hemorrhage |
| 30067221<br>12-2023-00188 | NR         | Malfunction | Apollo Endosurgery, Inc  | 10/2/2023     | LTI          | Orbera Intra gastric Balloon System | Insufficient Information   | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available                 |
| 30067221<br>12-2023-00187 | NR         | Malfunction | Apollo Endosurgery, Inc  | 10/2/2023     | LTI          | Orbera Intra gastric Balloon System | Migration  | Abdominal Pain; Obstruction/ Occlusion  |
| 30067221<br>12-2023-00198 | 4/24/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/29/2023     | LTI          | Orbera Intra gastric Balloon        | Deflation Problem; Unintended Deflation  | Failure of Implant  |
| 30126389<br>28-2023-02909 | 4/12/2023  | Injury      | Spatz Fgia Inc.          | 9/28/2023     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem; Migration or Expulsion of Device                                    | Obstruction/ Occlusion  |
| 30126389<br>28-2023-02913 | 7/27/2023  | Injury      | Spatz Fgia Inc.          | 9/28/2023     | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility; Patient Device Interaction Problem; Air/ Gas in Device | Stomach Ulceration  |



| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                       | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|----------------------------------|--|---|
| 30126389<br>28-2023-02912 | 8/4/2023   | Malfunction | Spatz Fgia Inc.          | 9/28/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02914 | 8/27/2023  | Malfunction | Spatz Fgia Inc.          | 9/27/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02911 | 3/17/2023  | Malfunction | Spatz Fgia Inc.          | 9/26/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02910 | 6/29/2023  | Malfunction | Spatz Fgia Inc.          | 9/26/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00163 | 7/17/2023  | Malfunction | Apollo Endosurgery, Inc  | 9/21/2023     | LTI          | Orbera Intra gastric Balloon     | Inflation Problem  | Failure of Implant                        |
| 30067221<br>12-2023-00184 | 4/25/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/19/2023     | LTI          | Orbera Intra gastric Balloon     | Inflation Problem; Air/ Gas in Device                      | Failure of Implant                        |
| 30067221<br>12-2023-00182 | 6/28/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/19/2023     | LTI          | Orbera Intra gastric Balloon     | Inflation Problem; Air/ Gas in Device                      | Failure of Implant                        |
| 30067221<br>12-2023-00183 | 7/14/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/19/2023     | LTI          | Orbera Intra gastric Balloon     | Inflation Problem; Air/ Gas in Device                      | Failure of Implant                        |
| 30067221<br>12-2023-00185 | 7/16/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/19/2023     | LTI          | Orbera Intra gastric Balloon     | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2023-00192 | NR         | Malfunction | Apollo Endosurgery, Inc  | 9/18/2023     | LTI          | Orbera Intra-gastric Balloon System    | Migration  | Perforation of Esophagus                  |
| 30067221<br>12-2023-00191 | 6/17/2023  | Malfunction | Apollo Endosurgery, Inc  | 9/15/2023     | LTI          | Orbera365 Intra-gastric Balloon System | Inflation Problem; Insufficient Information  | Failure of Implant                        |
| 30067221<br>12-2023-00180 | 7/21/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/15/2023     | LTI          | Orbera 365 Intra-gastric Balloon       | Deflation Problem; Fluid/ Blood Leak; Migration or Expulsion of Device; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2023-00181 | 8/22/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/15/2023     | LTI          | Orbera 365 Intra-gastric Balloon       | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                                   | Failure of Implant                        |
| 30126389<br>28-2023-02907 | 8/7/2023   | Malfunction | Spatz Fgia Inc.          | 9/14/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02908 | 8/21/2023  | Malfunction | Spatz Fgia Inc.          | 9/14/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02903 | 8/14/2023  | Malfunction | Spatz Fgia Inc.          | 9/13/2023     | LTI          | Spatz3 Adjustable Balloon System       | Inflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                       | Device Problem   | Patient Problem                            |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|----------------------------------|--|--|
| 30126389<br>28-2023-02904 | 7/31/2023  | Malfunction | Spatz Fgia Inc.          | 9/12/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions  |
| 30126389<br>28-2023-02906 | 8/9/2023   | Malfunction | Spatz Fgia Inc.          | 9/12/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device        | No Clinical Signs, Symptoms or Conditions  |
| 30126389<br>28-2023-02902 | 8/11/2023  | Malfunction | Spatz Fgia Inc.          | 9/11/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions  |
| 30126389<br>28-2023-02905 | 8/16/2023  | Malfunction | Spatz Fgia Inc.          | 9/10/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device        | Abdominal Pain; Diarrhea; Nausea; Vomiting |
| 30067221<br>12-2023-00170 | 7/19/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/8/2023      | LTI          | Orbera 365 Intra gastric Balloon | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant; Vomiting               |
| 30126389<br>28-2023-02901 | 7/19/2023  | Malfunction | Spatz Fgia Inc.          | 9/8/2023      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions  |
| 30067221<br>12-2023-00172 | 7/31/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/8/2023      | LTI          | Orbera 365 Intra gastric Balloon | Insufficient Information                                   | Insufficient Information                   |
| 30126389<br>28-2023-02896 | 8/1/2023   | Malfunction | Spatz Fgia Inc.          | 9/7/2023      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions  |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem  |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|--|
| 30126389<br>28-2023-02900 | 8/9/2023   | Malfunction | Spatz Fgia Inc.          | 9/7/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions  |
| 30126389<br>28-2023-02897 | 8/1/2023   | Malfunction | Spatz Fgia Inc.          | 9/6/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions  |
| 30126389<br>28-2023-02895 | 8/7/2023   | Malfunction | Spatz Fgia Inc.          | 9/4/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions  |
| 30067221<br>12-2023-00168 | 7/19/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/1/2023      | LTI          | Orbera 365 Intra gastric Balloon       | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                                   | Failure of Implant   |
| 30067221<br>12-2023-00169 | 8/4/2023   | Malfunction | Apollo Endosurgery, Inc. | 9/1/2023      | LTI          | Orbera 365 Intra gastric Balloon       | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                                   | Failure of Implant   |
| 30067221<br>12-2023-00175 | 7/1/2023   | Malfunction | Apollo Endosurgery, Inc  | 8/30/2023     | LTI          | Orbera365 Intra gastric Balloon System | Increase in Pressure   | Dehydration; Diarrhea; Vomiting; Pancreatitis; Unspecified Kidney or Urinary Problem |
| 30067221<br>22-2023-00167 | 7/17/2023  | Malfunction | Apollo Endosurgery, Inc. | 8/30/2023     | LTI          | Orbera 365 Intra gastric Balloon       | Deflation Problem; Fluid/ Blood Leak; Migration or Expulsion of Device; Unintended Deflation | Failure of Implant   |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem  |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|--|
| 30067221<br>12-2023-00166 | 8/1/2023   | Malfunction | Apollo Endosurgery, Inc. | 8/30/2023     | LTI          | Orbera 365 Intra gastric Balloon       | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation         | Failure of Implant   |
| 30067221<br>22-2023-00178 | 9/1/2022   | Malfunction | Apollo Endosurgery, Inc. | 8/28/2023     | LTI          | Orbera 365 Intra gastric Balloon       | Patient-Device Incompatibility; Patient Device Interaction Problem | Dehydration; Failure of Implant  |
| 30067221<br>12-2023-00177 | 2/2/2023   | Malfunction | Apollo Endosurgery, Inc. | 8/28/2023     | LTI          | Orbera 365 Intra gastric Balloon       | Patient-Device Incompatibility; Patient Device Interaction Problem | Dehydration; Failure of Implant; Vomiting  |
| 30067221<br>12-2023-00173 | 4/20/2023  | Malfunction | Apollo Endosurgery, Inc  | 8/28/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                            | Failure of Implant   |
| 30126389<br>28-2023-02899 | 6/9/2023   | Malfunction | Spatz Fgia Inc.          | 8/28/2023     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device                 | No Clinical Signs, Symptoms or Conditions  |
| 30067221<br>12-2023-00174 | 7/3/2023   | Malfunction | Apollo Endosurgery, Inc  | 8/28/2023     | LTI          | Orbera365 Intra gastric Balloon System | Partial Blockage   | Abdominal Pain; Pyrosis/ Heartburn; Failure of Implant; Nausea; Obstruction/ Occlusion |
| 30067221<br>12-2023-00176 | 7/6/2023   | Malfunction | Apollo Endosurgery, Inc  | 8/28/2023     | LTI          | Orbera365 Intra gastric Balloon System | Increase in Pressure   | Abdominal Pain; Nausea; Pancreatitis   |

| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                | Device Problem  | Patient Problem                                 |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|---|---|---|
| 30126389<br>28-2023-<br>02898 | 7/6/2023   | Malfunction | Spatz Fgia<br>Inc.             | 8/28/2023     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System | Patient-Device<br>Incompatibility;<br>Air/ Gas in<br>Device   | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30067221<br>22-2023-<br>00165 | 7/31/2023  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 8/28/2023     | LTI          | Orbera 365<br>Intragastric<br>Balloon     | Inflation Problem   | Fever; Failure of<br>Implant                    |
| 30067221<br>12-2023-<br>00161 | 3/8/2023   | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 8/25/2023     | LTI          | Orbera 365<br>Intragastric<br>Balloon     | Deflation<br>Problem; Fluid/<br>Blood Leak;<br>Migration or<br>Expulsion of<br>Device;<br>Unintended<br>Deflation | Failure of Implant;<br>Vomiting                 |
| 30067221<br>12-2023-<br>00162 | 5/25/2023  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 8/25/2023     | LTI          | Orbera 365<br>Intragastric<br>Balloon     | Inflation Problem   | Failure of Implant                              |
| 30067221<br>22-2023-<br>00163 | 7/17/2023  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 8/25/2023     | LTI          | Orbera<br>Intragastric<br>Balloon         | Inflation Problem   | Failure of Implant                              |
| 30067221<br>12-2023-<br>00164 | 7/19/2023  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 8/25/2023     | LTI          | Orbera<br>Intragastric<br>Balloon         | Inflation Problem   | Failure of Implant                              |
| 30126389<br>28-2023-<br>02894 | 7/25/2023  | Malfunction | Spatz Fgia<br>Inc.             | 8/23/2023     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System | Deflation<br>Problem  | Nausea  |
| 30126389<br>28-2023-<br>02887 | NR         | Injury      | Spatz Fgia<br>Inc.             | 8/23/2023     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System | Adverse Event<br>Without<br>Identified Device<br>or Use Problem   | Coma  |
| 30126389<br>28-2023-<br>02892 | 5/2/2023   | Malfunction | Spatz Fgia<br>Inc.             | 8/22/2023     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System | Deflation<br>Problem  | No Clinical Signs,<br>Symptoms or<br>Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                       | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|----------------------------------|--|---|
| 30126389<br>28-2023-02890 | 6/6/2023   | Malfunction | Spatz Fgia Inc.          | 8/22/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02893 | 7/25/2023  | Malfunction | Spatz Fgia Inc.          | 8/22/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02888 | 5/2/2023   | Malfunction | Spatz Fgia Inc.          | 8/21/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | Discomfort                                |
| 30067221<br>12-2023-00158 | 5/26/2023  | Malfunction | Apollo Endosurgery, Inc. | 8/21/2023     | LTI          | Orbera 365 Intra gastric Balloon | Deflation Problem; Fluid/ Blood Leak; Migration or Expulsion of Device; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2023-00157 | 7/26/2023  | Malfunction | Apollo Endosurgery, Inc. | 8/21/2023     | LTI          | Orbera 365 Intra gastric Balloon | Deflation Problem; Fluid/ Blood Leak; Migration or Expulsion of Device; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2023-02883 | 7/21/2023  | Malfunction | Spatz Fgia Inc.          | 8/20/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02884 | 7/21/2023  | Malfunction | Spatz Fgia Inc.          | 8/20/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                       | Device Problem  | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|----------------------------------|---|---|
| 30067221<br>22-2023-00156 | 6/21/2023  | Malfunction | Apollo Endosurgery, Inc. | 8/18/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>22-2023-00155 | 7/9/2023   | Malfunction | Apollo Endosurgery, Inc. | 8/18/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2023-02829 | 4/26/2023  | Malfunction | Spatz Fgia Inc.          | 8/17/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02881 | 6/6/2023   | Malfunction | Spatz Fgia Inc.          | 8/16/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device       | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>22-2023-00153 | 7/14/2023  | Malfunction | Apollo Endosurgery, Inc. | 8/16/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>22-2023-00159 | 7/13/2023  | Malfunction | Apollo Endosurgery, Inc. | 8/15/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Separation Problem; Premature Separation                  | Failure of Implant                        |
| 30067221<br>12-2023-00150 | 6/14/2023  | Malfunction | Apollo Endosurgery, Inc. | 8/14/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant                        |



| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                  | Device Problem  | Patient Problem                                 |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|---|---|---|
| 30126389<br>28-2023-<br>02880 | 7/13/2023  | Malfunction | Spatz Fgia<br>Inc.             | 8/13/2023     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System   | Deflation<br>Problem  | No Clinical Signs,<br>Symptoms or<br>Conditions |
| MW5135<br>363                 | 2/23/1962  | Death       | Apollo<br>Endosurgery,<br>Inc  | 8/12/2023     | LTI          | Orbera<br>Gastric<br>Balloon                | Adverse Event<br>Without<br>Identified Device<br>or Use Problem                             | Unspecified<br>Infection; Pain;<br>Perforation  |
| 30067221<br>12-2023-<br>00151 | 5/30/2023  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 8/11/2023     | LTI          | Orbera<br>Intragastric<br>Balloon<br>System | Inflation Problem   | Failure of Implant                              |
| 30067221<br>12-2023-<br>00152 | 7/2/2023   | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 8/11/2023     | LTI          | Orbera<br>Intragastric<br>Balloon<br>System | Inflation Problem   | Failure of Implant                              |
| 30067221<br>12-2023-<br>00148 | 7/9/2023   | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 8/11/2023     | LTI          | Orbera 365<br>Intragastric<br>Balloon       | Deflation<br>Problem;<br>Migration or<br>Expulsion of<br>Device;<br>Unintended<br>Deflation | Failure of Implant                              |
| 30126389<br>28-2023-<br>02876 | 6/21/2023  | Malfunction | Spatz Fgia<br>Inc.             | 8/10/2023     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System   | Air/ Gas in<br>Device   | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30126389<br>28-2023-<br>02878 | 7/3/2023   | Malfunction | Spatz Fgia<br>Inc.             | 8/10/2023     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System   | Deflation<br>Problem  | No Clinical Signs,<br>Symptoms or<br>Conditions |

| Report Number                 | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem                           |
|-------------------------------|------------|-------------|-------------------------|---------------|--------------|--|---|---|
| 30126389<br>28-2023-<br>02875 | 1/11/2023  | Injury      | Spatz Fgia Inc.         | 8/9/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Patient-Device Incompatibility; Patient Device Interaction Problem | Abdominal Pain                            |
| 30126389<br>28-2023-<br>02873 | 4/18/2023  | Malfunction | Spatz Fgia Inc.         | 8/9/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-<br>00149 | 6/13/2023  | Malfunction | Apollo Endosurgery, Inc | 8/9/2023      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                            | Failure of Implant                        |
| 30126389<br>28-2023-<br>02877 | 7/5/2023   | Malfunction | Spatz Fgia Inc.         | 8/9/2023      | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device                                    | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-<br>00147 | 7/11/2023  | Malfunction | Apollo Endosurgery, Inc | 8/9/2023      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                            | Failure of Implant                        |
| 30126389<br>28-2023-<br>02874 | 5/24/2023  | Malfunction | Spatz Fgia Inc.         | 8/8/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-<br>02871 | NR         | Malfunction | Spatz Fgia Inc.         | 8/8/2023      | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device                                    | Abdominal Pain; Discomfort                |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                       | Device Problem  | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|----------------------------------|---|---|
| 30126389<br>28-2023-02869 | 5/25/2023  | Malfunction | Spatz Fgia Inc.          | 8/6/2023      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02872 | 7/7/2023   | Malfunction | Spatz Fgia Inc.          | 8/6/2023      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02867 | 5/2/2023   | Malfunction | Spatz Fgia Inc.          | 8/3/2023      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02863 | 5/11/2023  | Injury      | Spatz Fgia Inc.          | 8/3/2023      | LTI          | Spatz3 Adjustable Balloon System | Patient-Device Incompatibility; Air/ Gas in Device  | Abdominal Pain; Nausea; Pancreatitis      |
| 30126389<br>28-2023-02868 | 5/26/2023  | Malfunction | Spatz Fgia Inc.          | 8/3/2023      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00146 | 3/21/2023  | Malfunction | Apollo Endosurgery, Inc. | 8/1/2023      | LTI          | Orbera 365 Intragastic Balloon   | Deflation Problem; Fluid/ Blood Leak; Migration or Expulsion of Device; Migration; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2023-02865 | 7/3/2023   | Malfunction | Spatz Fgia Inc.          | 7/31/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                       | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|----------------------------------|--|---|
| 30067221<br>12-2023-00143 | 6/19/2023  | Malfunction | Apollo Endosurgery, Inc. | 7/28/2023     | LTI          | Orbera 365 Intra gastric Balloon | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation         | Failure of Implant  |
| 30126389<br>28-2023-02858 | 6/22/2023  | Malfunction | Spatz Fgia Inc.          | 7/26/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions   |
| 30067221<br>12-2023-00144 | 6/30/2023  | Injury      | Apollo Endosurgery, Inc. | 7/26/2023     | LTI          | Orbera Intra gastric Balloon     | Patient-Device Incompatibility; Patient Device Interaction Problem | Abdominal Pain; Failure of Implant; Nausea; Vomiting; Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30126389<br>28-2023-02855 | 6/21/2023  | Malfunction | Spatz Fgia Inc.          | 7/24/2023     | LTI          | Spatz3 Adjustable Balloon System | Inflation Problem  | No Clinical Signs, Symptoms or Conditions   |
| 30126389<br>28-2023-02860 | 6/29/2023  | Malfunction | Spatz Fgia Inc.          | 7/24/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions   |
| 30067221<br>12-2023-00141 | 2/6/2023   | Malfunction | Apollo Endosurgery, Inc. | 7/21/2023     | LTI          | Orbera 365 Intra gastric Balloon | Migration or Expulsion of Device; Expulsion                        | Failure of Implant  |
| 30126389<br>28-2023-02861 | 6/16/2023  | Malfunction | Spatz Fgia Inc.          | 7/20/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions   |
| 30126389<br>28-2023-02859 | 6/27/2023  | Malfunction | Spatz Fgia Inc.          | 7/20/2023     | LTI          | Spatz3 Adjustable Balloon System | Patient-Device Incompatibility; Air/ Gas in Device                 | No Clinical Signs, Symptoms or Conditions   |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                       | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|----------------------------------|--|---|
| 30126389<br>28-2023-02854 | 6/1/2023   | Injury      | Spatz Fgia Inc.          | 7/19/2023     | LTI          | Spatz3 Adjustable Balloon System | Inadequate or Insufficient Training; Improper or Incorrect Procedure or Method; Misassembly by Users | Laceration(s); Perforation of Esophagus   |
| 30067221<br>12-2023-00140 | 6/15/2023  | Malfunction | Apollo Endosurgery, Inc. | 7/19/2023     | LTI          | Orbera Intra gastric Balloon     | Separation Problem; Premature Separation   | Failure of Implant                        |
| 30067221<br>12-2023-00142 | 6/22/2023  | Malfunction | Apollo Endosurgery, Inc. | 7/18/2023     | LTI          | Orbera 365 Intra gastric Balloon | Deflation Problem; Fluid/ Blood Leak; Migration or Expulsion of Device; Unintended Deflation         | Failure of Implant                        |
| 30126389<br>28-2023-02853 | 2/17/2023  | Malfunction | Spatz Fgia Inc.          | 7/16/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02850 | 6/14/2023  | Malfunction | Spatz Fgia Inc.          | 7/13/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02852 | 4/24/2023  | Malfunction | Spatz Fgia Inc.          | 7/11/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number                 | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem                                      |
|-------------------------------|------------|-------------|--------------------------|---------------|--------------|-------------------------------------|--|--|
| 30126389<br>28-2023-<br>02845 | 5/12/2023  | Malfunction | Spatz Fgia Inc.          | 7/11/2023     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem; Migration or Expulsion of Device        | Constipation   |
| 30067221<br>12-2023-<br>00139 | 6/12/2023  | Malfunction | Apollo Endosurgery, Inc. | 7/10/2023     | LTI          | Orbera 365 Intra gastric Balloon    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                                   |
| 30126389<br>28-2023-<br>02849 | 6/5/2023   | Malfunction | Spatz Fgia Inc.          | 7/5/2023      | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | No Clinical Signs, Symptoms or Conditions            |
| 30126389<br>28-2023-<br>02848 | 6/13/2023  | Malfunction | Spatz Fgia Inc.          | 7/5/2023      | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | No Clinical Signs, Symptoms or Conditions            |
| 30126389<br>28-2023-<br>02846 | 5/30/2023  | Malfunction | Spatz Fgia Inc.          | 7/3/2023      | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility; Air/ Gas in Device         | Abdominal Pain; Dehydration; Nausea; Vomiting        |
| 30126389<br>28-2023-<br>02841 | 5/23/2023  | Injury      | Spatz Fgia Inc.          | 6/29/2023     | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility                             | Nausea; Vomiting                                     |
| 30067221<br>12-2023-<br>00132 | 5/25/2023  | Injury      | Apollo Endosurgery, Inc. | 6/23/2023     | LTI          | Orbera 365 Intra gastric Balloon    | Difficult to Remove  | Laceration(s) of Esophagus; Perforation of Esophagus |
| 30067221<br>12-2023-<br>00136 | NR         | Malfunction | Apollo Endosurgery, Inc  | 6/21/2023     | LTI          | Orbera Intra gastric Balloon System | Inflation Problem  | Nausea; Vomiting; Pancreatitis                       |

| Report Number                 | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem                                     | Patient Problem  |
|-------------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|--|
| 30126389<br>28-2023-<br>02847 | 5/30/2023  | Malfunction | Spatz Fgia Inc.          | 6/20/2023     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device | Abdominal Pain   |
| 30126389<br>28-2023-<br>02844 | 6/7/2023   | Malfunction | Spatz Fgia Inc.          | 6/20/2023     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device | Abdominal Pain; Nausea   |
| 30126389<br>28-2023-<br>02839 | 5/25/2023  | Malfunction | Spatz Fgia Inc.          | 6/19/2023     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device | Abdominal Pain; Nausea   |
| 30067221<br>12-2023-<br>00128 | 7/1/2022   | Malfunction | Apollo Endosurgery, Inc  | 6/16/2023     | LTI          | Orbera Intra gastric Balloon System    | Migration or Expulsion of Device; Migration        | Perforation  |
| 30067221<br>12-2023-<br>00131 | 3/23/2023  | Malfunction | Apollo Endosurgery, Inc  | 6/16/2023     | LTI          | Orbera Intra gastric Balloon System    | Migration  | Abdominal Pain; Failure of Implant; Sepsis; Hernia; Obstruction/ Occlusion |
| 30067221<br>12-2023-<br>00130 | 5/12/2023  | Malfunction | Apollo Endosurgery, Inc  | 6/16/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation            | Failure of Implant   |
| 30067221<br>12-2023-<br>00111 | 5/21/2023  | Malfunction | Apollo Endosurgery, Inc. | 6/16/2023     | LTI          | Orbera Intra gastric Balloon           | Deflation Problem; Unintended Deflation            | Abdominal Pain; Failure of Implant   |
| 30126389<br>28-2023-<br>02838 | 5/18/2023  | Malfunction | Spatz Fgia Inc.          | 6/15/2023     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device | Nausea; Vomiting   |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|---|
| 30126389<br>28-2023-02836 | 4/23/2023  | Malfunction | Spatz Fgia Inc.          | 6/13/2023     | LTI          | Spatz3 Adjustable Balloon System       | Insufficient Information                              | No Clinical Signs, Symptoms or Conditions   |
| 30126389<br>28-2023-02837 | 5/20/2023  | Malfunction | Spatz Fgia Inc.          | 6/13/2023     | LTI          | Spatz3 Adjustable Balloon System       | Insufficient Information                              | No Clinical Signs, Symptoms or Conditions   |
| 30067221<br>12-2023-00126 | 5/16/2023  | Malfunction | Apollo Endosurgery, Inc  | 6/6/2023      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation               | Failure of Implant  |
| 30067221<br>12-2023-00124 | 5/6/2023   | Malfunction | Apollo Endosurgery, Inc  | 6/5/2023      | LTI          | Orbera365 Intra gastric Balloon System | Migration   | Abdominal Pain; Dehydration; Nausea; Vomiting; Obstruction/ Occlusion; Kidney Infection |
| 30067221<br>12-2023-00118 | NR         | Malfunction | Apollo Endosurgery, Inc. | 6/5/2023      | LTI          | Orbera 365 Intra gastric Balloon       | Deflation Problem; Unintended Deflation               | Abdominal Pain; Failure of Implant  |
| 30067221<br>12-2023-00112 | 2/10/2023  | Malfunction | Apollo Endosurgery, Inc  | 6/2/2023      | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Leak/ Splash; Unintended Deflation | Failure of Implant  |
| 30067221<br>12-2023-00116 | 3/10/2023  | Malfunction | Apollo Endosurgery, Inc  | 6/2/2023      | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem                                     | Failure of Implant  |
| 30067221<br>12-2023-00117 | 4/10/2023  | Malfunction | Apollo Endosurgery, Inc  | 6/2/2023      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Leak/ Splash; Unintended Deflation | Failure of Implant  |



| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem                    |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|------------------------------------|
| 30067221<br>12-2023-00115 | 4/20/2023  | Malfunction | Apollo Endosurgery, Inc  | 6/2/2023      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                     | Failure of Implant                 |
| 30067221<br>12-2023-00120 | 4/28/2023  | Malfunction | Apollo Endosurgery, Inc  | 6/2/2023      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                     | Failure of Implant                 |
| 30067221<br>12-2023-00119 | 5/4/2023   | Malfunction | Apollo Endosurgery, Inc  | 6/2/2023      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation  | Failure of Implant                 |
| 30067221<br>12-2023-00108 | 3/29/2023  | Malfunction | Apollo Endosurgery, Inc. | 5/31/2023     | LTI          | Orbera Intra gastric Balloon           | Deflation Problem; Unintended Deflation                     | Failure of Implant                 |
| 30067221<br>12-2023-00114 | 4/22/2023  | Malfunction | Apollo Endosurgery, Inc. | 5/30/2023     | LTI          | Orbera 365 Intra gastric Balloon       | Deflation Problem; Unintended Deflation                     | Failure of Implant; Vomiting       |
| 30067221<br>12-2023-00109 | 5/2/2023   | Malfunction | Apollo Endosurgery, Inc. | 5/26/2023     | LTI          | Orbera Intra gastric Balloon           | Deflation Problem; Mechanical Problem; Unintended Deflation | Failure of Implant                 |
| 30126389<br>28-2023-02831 | 4/14/2023  | Malfunction | Spatz Fgia Inc.          | 5/24/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | Abdominal Pain                     |
| 30067221<br>12-2023-00107 | 4/4/2023   | Malfunction | Apollo Endosurgery, Inc. | 5/23/2023     | LTI          | Obera Intra gastric Balloon            | Patient Device Interaction Problem                          | Vomiting; Perforation of Esophagus |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem                                     | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2023-00105 | 4/4/2023   | Malfunction | Apollo Endosurgery, Inc. | 5/19/2023     | LTI          | Orbera Intra gastric Balloon           | Deflation Problem; Unintended Deflation            | Failure of Implant                        |
| 30067221<br>12-2023-00102 | 4/8/2023   | Malfunction | Apollo Endosurgery, Inc. | 5/19/2023     | LTI          | Orbera 365 Intra gastric Balloon       | Deflation Problem; Unintended Deflation            | Failure of Implant                        |
| 30067221<br>12-2023-00103 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/17/2023     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem                                  | Hypersensitivity/ Allergic reaction       |
| 30126389<br>28-2023-02821 | 3/7/2023   | Malfunction | Spatz Fgia Inc.          | 5/16/2023     | LTI          | Spatz3 Adjustable Balloon System       | Insufficient Information                           | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00085 | 3/20/2023  | Malfunction | Apollo Endosurgery, Inc  | 5/15/2023     | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Unintended Deflation            | Failure of Implant                        |
| 30067221<br>12-2023-00093 | 4/10/2023  | Malfunction | Apollo Endosurgery, Inc  | 5/15/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation            | Failure of Implant                        |
| 30067221<br>12-2023-00089 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/15/2023     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem; Air/ Gas in Device              | Failure of Implant                        |
| 30067221<br>12-2023-00084 | 3/26/2023  | Malfunction | Apollo Endosurgery, Inc. | 5/12/2023     | LTI          | Orbera 365 Intra gastric Balloon       | Deflation Problem; Unintended Deflation            | Failure of Implant                        |
| 30126389<br>28-2023-02833 | 4/19/2023  | Malfunction | Spatz Fgia Inc.          | 5/10/2023     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device | Abdominal Pain; Nausea; Vomiting          |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem                                  |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|-------------------------------------|--|--|
| 30126389<br>28-2023-02834 | 5/7/2023   | Malfunction | Spatz Fgia Inc.          | 5/9/2023      | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | No Clinical Signs, Symptoms or Conditions        |
| 30126389<br>28-2023-02830 | 4/6/2023   | Malfunction | Spatz Fgia Inc.          | 5/7/2023      | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility; Air/ Gas in Device         | No Clinical Signs, Symptoms or Conditions        |
| 30067221<br>12-2023-00083 | 4/13/2023  | Malfunction | Apollo Endosurgery, Inc. | 5/5/2023      | LTI          | Orbera Intra gastric Balloon System | Adverse Event Without Identified Device or Use Problem     | Abdominal Pain; Vomiting; Aspiration Pneumonitis |
| 30126389<br>28-2023-02820 | 4/14/2023  | Malfunction | Spatz Fgia Inc.          | 5/1/2023      | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | No Clinical Signs, Symptoms or Conditions        |
| 30067221<br>12-2023-00080 | 2/18/2023  | Malfunction | Apollo Endosurgery, Inc. | 4/28/2023     | LTI          | Orbera Intra gastric Balloon        | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                               |
| 30067221<br>12-2023-00081 | 4/1/2023   | Malfunction | Apollo Endosurgery, Inc. | 4/28/2023     | LTI          | Orbera 365 Intra gastric Balloon    | Deflation Problem; Unintended Deflation                    | Failure of Implant                               |
| 30126389<br>28-2023-02828 | 1/18/2023  | Malfunction | Spatz Fgia Inc.          | 4/27/2023     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | No Clinical Signs, Symptoms or Conditions        |
| 30126389<br>28-2023-02825 | 1/18/2023  | Malfunction | Spatz Fgia Inc.          | 4/27/2023     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | No Clinical Signs, Symptoms or Conditions        |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem                                     | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2023-02817 | 3/11/2023  | Malfunction | Spatz Fgia Inc.          | 4/27/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00091 | 3/24/2023  | Malfunction | Apollo Endosurgery, Inc  | 4/26/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation            | Failure of Implant                        |
| 30067221<br>12-2023-00092 | 3/24/2023  | Malfunction | Apollo Endosurgery, Inc  | 4/26/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation            | Failure of Implant                        |
| 30126389<br>28-2023-02813 | 6/15/2022  | Malfunction | Spatz Fgia Inc.          | 4/24/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02815 | 1/15/2023  | Malfunction | Spatz Fgia Inc.          | 4/24/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | Discomfort                                |
| 30126389<br>28-2023-02824 | 1/18/2023  | Malfunction | Spatz Fgia Inc.          | 4/24/2023     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device | Discomfort                                |
| 30126389<br>28-2023-02814 | 3/29/2023  | Malfunction | Spatz Fgia Inc.          | 4/24/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00079 | 1/19/2023  | Malfunction | Apollo Endosurgery, Inc. | 4/21/2023     | LTI          | Orbera 365 Intra gastric Balloon       | Deflation Problem; Unintended Deflation            | Abdominal Pain; Failure of Implant        |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                       | Device Problem  | Patient Problem                                 |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|----------------------------------|---|---|
| 30067221<br>12-2023-00078 | 2/13/2023  | Malfunction | Apollo Endosurgery, Inc. | 4/21/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Adverse Event Without Identified Device or Use Problem; Migration | Abdominal Pain; Failure of Implant; Perforation |
| 30126389<br>28-2023-02822 | 1/29/2023  | Malfunction | Spatz Fgia Inc.          | 4/20/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions       |
| 30126389<br>28-2023-02823 | 3/10/2023  | Malfunction | Spatz Fgia Inc.          | 4/20/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem   | No Clinical Signs, Symptoms or Conditions       |
| 30126389<br>28-2023-02812 | 3/27/2023  | Malfunction | Spatz Fgia Inc.          | 4/19/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device               | Abdominal Pain; Diarrhea; Nausea; Vomiting      |
| 30126389<br>28-2023-02816 | NR         | Malfunction | Spatz Fgia Inc.          | 4/19/2023     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device               | No Clinical Signs, Symptoms or Conditions       |
| 30067221<br>12-2023-00073 | 3/2/2023   | Malfunction | Apollo Endosurgery, Inc. | 4/18/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Deflation Problem; Unintended Deflation                           | Failure of Implant                              |
| 30067221<br>12-2023-00066 | 1/8/2023   | Malfunction | Apollo Endosurgery, Inc. | 4/17/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Insufficient Information  | Necrosis  |
| 30067221<br>12-2023-00069 | 3/5/2023   | Malfunction | Apollo Endosurgery, Inc. | 4/17/2023     | LTI          | Orbera 365 Intra-gastric Balloon | Deflation Problem; Fluid/Blood Leak; Unintended Deflation         | Failure of Implant                              |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                            |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|--|
| 30067221<br>12-2023-00070 | 3/21/2023  | Malfunction | Apollo Endosurgery, Inc. | 4/17/2023     | LTI          | Orbera 365 Intra gastric Balloon       | Unintended Deflation                                       | Failure of Implant                         |
| 30067221<br>12-2023-00071 | NR         | Malfunction | Apollo Endosurgery, Inc. | 4/17/2023     | LTI          | BIB System Intra gastric Balloon       | Air/ Gas in Device   | Failure of Implant                         |
| 30126389<br>28-2023-02808 | 9/13/2022  | Malfunction | Spatz Fgia Inc.          | 4/13/2023     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device         | No Clinical Signs, Symptoms or Conditions  |
| 30067221<br>12-2023-00074 | NR         | Malfunction | Apollo Endosurgery, Inc  | 4/12/2023     | LTI          | Orbera® Intra gastric Balloon System   | Obstruction of Flow  | Obstruction/ Occlusion                     |
| 30067221<br>12-2023-00067 | 3/3/2023   | Malfunction | Apollo Endosurgery, Inc  | 4/4/2023      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                    | Failure of Implant                         |
| 30067221<br>12-2023-00065 | 3/14/2023  | Malfunction | Apollo Endosurgery, Inc  | 4/4/2023      | LTI          | Orbera Intra gastric Balloon System    | No Apparent Adverse Event                                  | Abdominal Pain; Flatus; Pyrosis/ Heartburn |
| 30126389<br>28-2023-02806 | 3/10/2023  | Malfunction | Spatz Fgia Inc.          | 4/3/2023      | LTI          | Spatz3 Adjustable Balloon System       | Insufficient Information                                   | No Clinical Signs, Symptoms or Conditions  |
| 30067221<br>12-2023-00068 | 12/1/2022  | Malfunction | Apollo Endosurgery, Inc. | 3/31/2023     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem  | Abdominal Pain; Failure of Implant         |
| 30067221<br>12-2023-00063 | 2/6/2023   | Malfunction | Apollo Endosurgery, Inc  | 3/31/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                         |

| Report Number             | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|-------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2023-00057 | 2/13/2023  | Malfunction | Apollo Endosurgery, Inc | 3/31/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2023-02810 | 3/23/2023  | Malfunction | Spatz Fgia Inc.         | 3/28/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02811 | 3/23/2023  | Malfunction | Spatz Fgia Inc.         | 3/28/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02798 | 2/3/2023   | Malfunction | Spatz Fgia Inc.         | 3/26/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00055 | 2/20/2023  | Malfunction | Apollo Endosurgery, Inc | 3/24/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                    | Failure of Implant                        |
| 30126389<br>28-2023-02801 | 2/18/2023  | Malfunction | Spatz Fgia Inc.         | 3/20/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02805 | NR         | Injury      | Spatz Fgia Inc.         | 3/20/2023     | LTI          | Spatz3 Adjustable Balloon System       | Obstruction of Flow  | Abdominal Pain; Nausea; Constipation      |
| 30067221<br>12-2023-00060 | 10/25/2022 | Malfunction | Apollo Endosurgery, Inc | 3/17/2023     | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |

| Report Number             | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|-------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2023-00052 | 1/31/2023  | Malfunction | Apollo Endosurgery, Inc | 3/17/2023     | LTI          | BIB Intra gastric Balloon System       | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-0230-00053 | 3/9/2023   | Malfunction | Apollo Endosurgery, Inc | 3/17/2023     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem  | Failure of Implant                        |
| 30067221<br>12-2023-00059 | 12/18/2022 | Malfunction | Apollo Endosurgery, Inc | 3/16/2023     | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2023-02803 | 2/24/2023  | Malfunction | Spatz Fgia Inc.         | 3/16/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00050 | 2/25/2023  | Malfunction | Apollo Endosurgery, Inc | 3/16/2023     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem  | Abdominal Pain; Pyrosis/ Heartburn        |
| 30067221<br>12-2023-00058 | NR         | Malfunction | Apollo Endosurgery, Inc | 3/16/2023     | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Gastritis; Failure of Implant             |
| 30126389<br>28-2023-02800 | 1/26/2023  | Malfunction | Spatz Fgia Inc.         | 3/15/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00048 | 12/25/2023 | Malfunction | Apollo Endosurgery, Inc | 3/15/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |



| Report Number                 | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|-------------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2023-<br>02795 | 10/29/2022 | Malfunction | Spatz Fgia Inc.          | 3/14/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration or Expulsion of Device                              | No Clinical Signs, Symptoms or Conditions |
| 16528829                      | 1/13/2023  | Malfunction | Apollo Endosurgery, Inc. | 3/13/2023     | LTI          | Orbera Intra gastric Balloon System    | Inadequacy of Device Shape and/or Size   | Abdominal Pain                            |
| 30126389<br>28-2023-<br>02792 | 1/13/2023  | Malfunction | Spatz Fgia Inc.          | 3/9/2023      | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device                               | Nausea; Vomiting                          |
| 30126389<br>28-2023-<br>02790 | 2/3/2023   | Malfunction | Spatz Fgia Inc.          | 3/9/2023      | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device                               | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-<br>02789 | 2/5/2023   | Malfunction | Spatz Fgia Inc.          | 3/6/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration or Expulsion of Device                              | Abdominal Pain                            |
| 30126389<br>28-2023-<br>02802 | 2/27/2023  | Malfunction | Spatz Fgia Inc.          | 3/2/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-<br>02785 | 1/31/2023  | Malfunction | Spatz Fgia Inc.          | 2/28/2023     | LTI          | Spatz3 Adjustable Balloon System       | Insufficient Information   | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-<br>00047 | 2/2/2023   | Malfunction | Apollo Endosurgery, Inc  | 2/24/2023     | LTI          | Orbera365 Intra gastric Balloon System | Adverse Event Without Identified Device or Use Problem; Insufficient Information | Renal Failure; Vomiting                   |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                    |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|------------------------------------|
| 30067221<br>12-2023-00049 | NR         | Malfunction | Apollo Endosurgery, Inc  | 2/24/2023     | LTI          | Orbera Intra gastric Balloon System    | Adverse Event Without Identified Device or Use Problem     | Fever; Pain                        |
| 30067221<br>12-2023-00046 | 3/5/2022   | Malfunction | Apollo Endosurgery, Inc  | 2/23/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                 |
| 30067221<br>12-2023-00042 | 2/11/2023  | Malfunction | Apollo Endosurgery, Inc  | 2/23/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                 |
| 30067221<br>12-2023-00041 | 2/16/2023  | Malfunction | Apollo Endosurgery, Inc. | 2/23/2023     | LTI          | Orbera Intra gastric Balloon Systems   | Fluid/ Blood Leak; Migration                               | Abdominal Pain; Failure of Implant |
| 30067221<br>12-2023-00043 | 2/20/2023  | Malfunction | Apollo Endosurgery, Inc  | 2/23/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                 |
| 30067221<br>12-2023-00044 | 1/30/2023  | Malfunction | Apollo Endosurgery, Inc  | 2/22/2023     | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak  | Perforation of Esophagus           |
| 30067221<br>12-2023-00037 | 2/9/2023   | Malfunction | Apollo Endosurgery, Inc. | 2/22/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                 |

| Report Number             | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem  |
|---------------------------|------------|-------------|-------------------------|---------------|--------------|--|--|--|
| 30067221<br>12-2023-00036 | 2/10/2023  | Malfunction | Apollo Endosurgery, Inc | 2/22/2023     | LTI          | Orbera Intra gastric Balloon System    | Partial Blockage; Obstruction of Flow                      | Dehydration; Failure of Implant; Vomiting; Obstruction/ Occlusion    |
| 30067221<br>12-2023-00035 | 12/1/2023  | Malfunction | Apollo Endosurgery, Inc | 2/22/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant   |
| 30067221<br>12-2023-00021 | 1/15/2023  | Malfunction | Apollo Endosurgery, Inc | 2/21/2023     | LTI          | Bib; Intra gastric Balloon System      | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant   |
| 30067221<br>12-2023-00026 | 1/24/2023  | Malfunction | Apollo Endosurgery, Inc | 2/21/2023     | LTI          | Orbera365 Intra gastric Balloon System | No Apparent Adverse Event                                  | Dehydration; Vomiting  |
| 30067221<br>12-2023-00022 | 1/26/2023  | Malfunction | Apollo Endosurgery, Inc | 2/21/2023     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem  | Failure of Implant   |
| 30067221<br>12-2023-00031 | NR         | Malfunction | Apollo Endosurgery, Inc | 2/21/2023     | LTI          | Bib; Intra gastric Balloon System      | Deflation Problem; Unintended Deflation                    | Abdominal Pain; Pyrosis/ Heartburn; Perforation of Vessels; Vomiting |
| 30067221<br>12-2023-00038 | NR         | Death       | Apollo Endosurgery, Inc | 2/21/2023     | LTI          | Orbera Intra gastric Balloon System    | Adverse Event Without Identified Device or Use Problem     | Perforation  |
| 30126389<br>28-2023-02783 | 12/21/2022 | Malfunction | Spatz Fgia Inc.         | 2/20/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | Abdominal Pain   |

| Report Number             | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem  |
|---------------------------|------------|-------------|-------------------------|---------------|--------------|--|--|--|
| 30126389<br>28-2023-02784 | 1/19/2023  | Malfunction | Spatz Fgia Inc.         | 2/20/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions                          |
| 30126389<br>28-2023-02766 | 10/18/2022 | Malfunction | Spatz Fgia Inc.         | 2/16/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | Abdominal Pain   |
| 30067221<br>12-2023-00030 | NR         | Malfunction | Apollo Endosurgery, Inc | 2/16/2023     | LTI          | Orbera Intra gastric Balloon System    | Partial Blockage   | Failure of Implant; Nausea; Pain; Vomiting; Obstruction/ Occlusion |
| 30067221<br>12-2023-00029 | NR         | Malfunction | Apollo Endosurgery, Inc | 2/16/2023     | LTI          | Orbera® Intra gastric Balloon System   | Inflation Problem  | Failure of Implant   |
| 30126389<br>28-2023-02776 | 12/25/2022 | Malfunction | Spatz Fgia Inc.         | 2/15/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions                          |
| 30067221<br>12-2023-00023 | 2/10/2022  | Malfunction | Apollo Endosurgery, Inc | 2/14/2023     | LTI          | Orbera Intra gastric Balloon System    | Patient-Device Incompatibility                             | Failure of Implant; Pain; Implant Pain                             |
| 30067221<br>12-2023-00034 | 11/1/2022  | Malfunction | Apollo Endosurgery, Inc | 2/14/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant   |
| 30067221<br>12-2023-00028 | 1/18/2023  | Malfunction | Apollo Endosurgery, Inc | 2/14/2023     | LTI          | Orbera365 Intra gastric Balloon System | Infusion or Flow Problem; Unintended Deflation             | Failure of Implant   |

| Report Number             | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                             | Device Problem                                     | Patient Problem                           |
|---------------------------|------------|-------------|-------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2023-00032 | 1/26/2023  | Malfunction | Apollo Endosurgery, Inc | 2/14/2023     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation            | Failure of Implant                        |
| 30126389<br>28-2023-02779 | 2/26/2022  | Malfunction | Spatz Fgia Inc.         | 2/13/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02778 | 11/20/2022 | Malfunction | Spatz Fgia Inc.         | 2/13/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | Abdominal Pain                            |
| 30126389<br>28-2023-02786 | 2/3/2023   | Malfunction | Spatz Fgia Inc.         | 2/9/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02775 | 9/27/2022  | Malfunction | Spatz Fgia Inc.         | 2/8/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02773 | 11/27/2022 | Malfunction | Spatz Fgia Inc.         | 2/8/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02777 | 12/26/2022 | Malfunction | Spatz Fgia Inc.         | 2/8/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02780 | 2/10/2022  | Malfunction | Spatz Fgia Inc.         | 2/7/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02768 | 1/10/2023  | Malfunction | Spatz Fgia Inc.         | 2/7/2023      | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device | Abdominal Pain; Nausea; Vomiting          |

| Report Number         | Event Date | Event Type  | Manufacturer                         | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem   |
|-----------------------|------------|-------------|--------------------------------------|---------------|--------------|--|--|---|
| MW5114745             | NR         | Injury      | Apollo Endosurgery Costa Rica S.R.L. | 2/6/2023      | LTI          | Orbera Intra Gastric Balloon           | Improper or Incorrect Procedure or Method                  | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 3006722112-2023-00024 | NR         | Malfunction | Apollo Endosurgery, Inc              | 2/3/2023      | LTI          | Orbera Intra Gastric Balloon System    | Inflation Problem  | Abdominal Pain; Failure of Implant; Perforation of Vessels; Pancreatitis  |
| 3006722112-2023-00025 | NR         | Malfunction | Apollo Endosurgery, Inc              | 2/3/2023      | LTI          | Orbera Intra Gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant  |
| 3006722112-2023-00020 | 7/28/2022  | Malfunction | Apollo Endosurgery, Inc              | 2/2/2023      | LTI          | Orbera365 Intra Gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant  |
| 3006722112-2023-00019 | 9/26/2022  | Malfunction | Apollo Endosurgery, Inc              | 2/2/2023      | LTI          | Orbera365 Intra Gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant  |
| 3012638928-2023-02769 | 11/2/2022  | Malfunction | Spatz Fgia Inc.                      | 2/2/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration                               | No Clinical Signs, Symptoms or Conditions                                 |
| 3006722112-2023-00005 | 11/8/2022  | Malfunction | Apollo Endosurgery, Inc              | 2/2/2023      | LTI          | Orbera365 Intra Gastric Balloon System | Pressure Problem   | Neuralgia   |
| 3006722112-2023-00017 | 12/22/2022 | Malfunction | Apollo Endosurgery, Inc              | 2/2/2023      | LTI          | Orbera365 Intra Gastric Balloon System | Fluid/ Blood Leak; Reflux within Device                    | Gastritis; Nausea; Obstruction/ Occlusion                                 |

| Report Number             | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|-------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2023-00007 | 1/13/2023  | Malfunction | Apollo Endosurgery, Inc | 2/2/2023      | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem; Air/ Gas in Device                      | Abdominal Pain; Failure of Implant        |
| 30126389<br>28-2023-02764 | 8/15/2022  | Malfunction | Spatz Fgia Inc.         | 2/1/2023      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2023-00001 | 11/4/2022  | Malfunction | Apollo Endosurgery, Inc | 2/1/2023      | LTI          | Orbera365 Intra gastric Balloon System | Adverse Event Without Identified Device or Use Problem     | Abdominal Pain; Perforation; Vomiting     |
| 30067221<br>12-2023-00015 | 11/21/2022 | Malfunction | Apollo Endosurgery, Inc | 2/1/2023      | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem  | Failure of Implant                        |
| 30067221<br>12-2023-00016 | 12/10/2022 | Malfunction | Apollo Endosurgery, Inc | 2/1/2023      | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2023-00003 | 1/12/2023  | Malfunction | Apollo Endosurgery, Inc | 2/1/2023      | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem  | Nausea; Vomiting                          |
| 30067221<br>12-2023-00018 | NR         | Malfunction | Apollo Endosurgery, Inc | 1/28/2023     | LTI          | Orbera Intra gastric Balloon System    | Appropriate Term/ Code Not Available                       | Vomiting; Stomach Ulceration              |
| 30067221<br>12-2023-00014 | NR         | Malfunction | Apollo Endosurgery, Inc | 1/25/2023     | LTI          | Orbera Intra gastric Balloon System    | Appropriate Term/ Code Not Available                       | Abdominal Pain; Bowel Perforation         |

| Report Number             | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem                                  |
|---------------------------|------------|-------------|-------------------------|---------------|--------------|-------------------------------------|--|--|
| 30067221<br>12-2023-00006 | NR         | Malfunction | Apollo Endosurgery, Inc | 1/24/2023     | LTI          | Orbera Intra gastric Balloon System | Deflation Problem; Migration or Expulsion of Device; Migration; Unintended Deflation | Abdominal Pain; Vomiting; Obstruction/ Occlusion |
| 30067221<br>12-2023-00012 | NR         | Malfunction | Apollo Endosurgery, Inc | 1/23/2023     | LTI          | Orbera Intra gastric Balloon System | Appropriate Term/ Code Not Available   | Failure of Implant; Pain; Vomiting; Constipation |
| 30067221<br>12-2023-00013 | NR         | Malfunction | Apollo Endosurgery, Inc | 1/23/2023     | LTI          | Orbera Intra gastric Balloon System | No Apparent Adverse Event  | Abdominal Pain; Diarrhea; Nausea; Constipation   |
| 30126389<br>28-2023-02762 | 12/15/2022 | Malfunction | Spatz Fgia Inc.         | 1/18/2023     | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility; Air/ Gas in Device                                   | No Clinical Signs, Symptoms or Conditions        |
| 30126389<br>28-2023-02761 | 12/20/2022 | Malfunction | Spatz Fgia Inc.         | 1/18/2023     | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility; Air/ Gas in Device                                   | Abdominal Pain                                   |
| 30126389<br>28-2023-02765 | 1/6/2023   | Malfunction | Spatz Fgia Inc.         | 1/18/2023     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | No Clinical Signs, Symptoms or Conditions        |
| 30126389<br>28-2023-02758 | 12/15/2022 | Injury      | Spatz Fgia Inc.         | 1/17/2023     | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility   | Abdominal Pain; Nausea; Vomiting                 |



| Report Number             | Event Date | Event Type  | Manufacturer            | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                    |
|---------------------------|------------|-------------|-------------------------|---------------|--------------|--|--|------------------------------------|
| 30067221<br>12-2023-00002 | 10/26/2022 | Malfunction | Apollo Endosurgery, Inc | 1/13/2023     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                           | Abdominal Pain; Failure of Implant |
| 30067221<br>12-2022-00134 | 11/22/2022 | Malfunction | Apollo Endosurgery, Inc | 1/13/2023     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                           | Failure of Implant                 |
| 30067221<br>12-2022-00136 | 12/13/2022 | Malfunction | Apollo Endosurgery, Inc | 1/13/2023     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                           | Failure of Implant                 |
| 30067221<br>12-2022-00135 | 12/20/2022 | Malfunction | Apollo Endosurgery, Inc | 1/13/2023     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Migration or Expulsion of Device; Migration; Unintended Deflation | Failure of Implant                 |
| 30067221<br>12-2022-00137 | 12/26/2022 | Malfunction | Apollo Endosurgery, Inc | 1/13/2023     | LTI          | Orbera Intra-gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                           | Failure of Implant                 |
| 30067221<br>12-2022-00133 | 10/7/2022  | Malfunction | Apollo Endosurgery, Inc | 1/11/2023     | LTI          | Orbera Intra-gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                           | Failure of Implant                 |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|---|
| 30067221<br>12-2022-00132 | 10/26/2022 | Malfunction | Apollo Endosurgery, Inc  | 1/11/2023     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2023-02760 | 11/17/2022 | Malfunction | Spatz Fgia Inc.          | 1/11/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2023-02757 | 12/6/2022  | Malfunction | Spatz Fgia Inc.          | 1/11/2023     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00131 | 12/13/2022 | Malfunction | Apollo Endosurgery, Inc  | 1/11/2023     | LTI          | Orbera Intra-gastric Balloon System    | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2022-00126 | 2/8/2022   | Malfunction | Apollo Endosurgery, Inc  | 12/29/2022    | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Unintended Deflation                   | Failure of Implant                        |
| 30067221<br>12-2022-00127 | 9/7/2022   | Malfunction | Apollo Endosurgery, Inc. | 12/29/2022    | LTI          | Orbera365 Intra-gastric Balloon System | Insufficient Information                                  | Obstruction/Occlusion                     |
| 30126389<br>28-2022-02751 | 12/5/2022  | Injury      | Spatz Fgia Inc.          | 12/29/2022    | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility                            | Abdominal Pain; Nausea; Vomiting          |
| 30126389<br>28-2022-02755 | 12/7/2022  | Malfunction | Spatz Fgia Inc.          | 12/21/2022    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                    |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|------------------------------------|
| 30126389<br>28-2022-02750 | 11/24/2022 | Malfunction | Spatz Fgia Inc.          | 12/20/2022    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | Abdominal Pain                     |
| 30126389<br>28-2022-02749 | 11/12/2022 | Injury      | Spatz Fgia Inc.          | 12/18/2022    | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility                             | Abdominal Pain; Stomach Ulceration |
| 30126389<br>28-2022-02753 | 11/4/2022  | Malfunction | Spatz Fgia Inc.          | 12/15/2022    | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device         | Abdominal Pain; Vomiting           |
| 30126389<br>28-2022-02754 | 12/3/2022  | Malfunction | Spatz Fgia Inc.          | 12/15/2022    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | Abdominal Pain                     |
| 30067221<br>12-2022-00130 | NR         | Malfunction | Apollo Endosurgery, Inc. | 12/15/2022    | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Abdominal Pain; Failure of Implant |
| 30067221<br>12-2022-00116 | 10/10/2022 | Malfunction | Apollo Endosurgery, Inc  | 12/14/2022    | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                 |
| 30067221<br>12-2022-00128 | 10/19/2022 | Malfunction | Apollo Endosurgery, Inc  | 12/14/2022    | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak  | Failure of Implant                 |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                                |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|--|
| 30067221<br>12-2022-00117 | 10/20/2022 | Malfunction | Apollo Endosurgery, Inc  | 12/14/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Mechanical Problem; Unintended Deflation | Failure of Implant                             |
| 30067221<br>12-2022-00118 | 10/23/2022 | Malfunction | Apollo Endosurgery, Inc  | 12/14/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                     | Failure of Implant                             |
| 30067221<br>12-2022-00125 | 10/30/2022 | Malfunction | Apollo Endosurgery, Inc  | 12/14/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                     | Failure of Implant                             |
| 30067221<br>12-2022-00101 | 7/27/2022  | Malfunction | Apollo Endosurgery, Inc. | 12/13/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                     | Failure of Implant                             |
| 30067221<br>12-2022-00121 | 9/9/2022   | Malfunction | Apollo Endosurgery, Inc. | 12/13/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                     | Failure of Implant; Laceration(s) of Esophagus |
| 30067221<br>12-2022-00120 | 9/13/2022  | Malfunction | Apollo Endosurgery, Inc. | 12/13/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                     | Failure of Implant                             |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2022-00123 | 9/19/2022  | Malfunction | Apollo Endosurgery, Inc. | 12/13/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2022-00100 | 9/23/2022  | Malfunction | Apollo Endosurgery, Inc. | 12/13/2022    | LTI          | Orbera365 Intra gastric Balloon System | Inflation Problem  | Failure of Implant                        |
| 30067221<br>12-2022-00122 | 10/7/2022  | Malfunction | Apollo Endosurgery, Inc. | 12/13/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2022-00119 | 10/19/2022 | Malfunction | Apollo Endosurgery, Inc. | 12/13/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2022-00124 | 10/22/2022 | Malfunction | Apollo Endosurgery, Inc. | 12/13/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                    | Failure of Implant                        |
| 30067221<br>12-2022-00114 | 10/24/2022 | Malfunction | Apollo Endosurgery, Inc. | 12/13/2022    | LTI          | Orbera Intra gastric Balloon System    | Infusion or Flow Problem                                   | Abdominal Pain; Vomiting                  |
| 30126389<br>28-2022-02748 | 11/7/2022  | Malfunction | Spatz Fgia Inc.          | 12/13/2022    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | Abdominal Pain                            |
| 30126389<br>28-2022-02744 | 11/13/2022 | Malfunction | Spatz Fgia Inc.          | 12/13/2022    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|-------------------------------------|--|---|
| 30067221<br>12-2022-00129 | 12/11/2022 | Malfunction | Apollo Endosurgery, Inc. | 12/13/2022    | LTI          | Orbera Intra-gastric Balloon System | Inflation Problem                                      | Abdominal Pain; Failure of Implant        |
| 30126389<br>28-2022-02440 | 10/10/2021 | Malfunction | Spatz Fgia Inc.          | 12/8/2022     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02502 | 11/4/2021  | Malfunction | Spatz Fgia Inc.          | 12/8/2022     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02521 | 2/1/2022   | Malfunction | Spatz Fgia Inc.          | 12/8/2022     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem; Migration or Expulsion of Device    | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02498 | 2/12/2022  | Injury      | Spatz Fgia Inc.          | 12/8/2022     | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility                         | Abdominal Pain                            |
| 30126389<br>28-2022-02514 | 2/13/2022  | Injury      | Spatz Fgia Inc.          | 12/8/2022     | LTI          | Spatz3 Adjustable Balloon System    | Adverse Event Without Identified Device or Use Problem | Pancreatitis                              |
| 30126389<br>28-2022-02620 | 6/24/2022  | Malfunction | Spatz Fgia Inc.          | 12/8/2022     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02739 | 5/10/2022  | Malfunction | Spatz Fgia Inc.          | 12/7/2022     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2022-02738 | 6/15/2022  | Malfunction | Spatz Fgia Inc.          | 12/7/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02752 | 12/1/2022  | Malfunction | Spatz Fgia Inc.          | 12/7/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02454 | NR         | Malfunction | Spatz Fgia Inc.          | 12/7/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration or Expulsion of Device    | Abdominal Pain                            |
| 30126389<br>28-2022-02519 | NR         | Malfunction | Spatz Fgia Inc.          | 12/7/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02735 | NR         | Malfunction | Spatz Fgia Inc.          | 12/7/2022     | LTI          | Spatz3 Adjustable Balloon System       | Air/ Gas in Device                                     | Abdominal Pain                            |
| 30126389<br>28-2021-02442 | 12/10/2021 | Malfunction | Spatz Fgia Inc.          | 12/6/2022     | LTI          | Spatz3 Adjustable Balloon System       | Leak/ Splash; Therapeutic or Diagnostic Output Failure | Failure of Implant; Vomiting              |
| 30067221<br>12-2022-00096 | 7/1/2022   | Malfunction | Apollo Endosurgery, Inc. | 12/6/2022     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem                                      | Failure of Implant                        |
| 30067221<br>12-2022-00094 | 8/18/2022  | Malfunction | Apollo Endosurgery, Inc. | 12/6/2022     | LTI          | Orbera365 Intra gastric Balloon System | Inflation Problem                                      | Failure of Implant                        |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem              |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|------------------------------|
| 30067221<br>12-2022-00113 | 9/13/2022  | Malfunction | Apollo Endosurgery, Inc. | 12/6/2022     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem   | Failure of Implant           |
| 30067221<br>12-2022-00099 | 9/22/2022  | Malfunction | Apollo Endosurgery, Inc. | 12/6/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                    | Failure of Implant           |
| 30067221<br>12-2022-00111 | 9/30/2022  | Malfunction | Apollo Endosurgery, Inc. | 12/6/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                    | Failure of Implant           |
| 30067221<br>12-2022-00108 | 10/17/2022 | Malfunction | Apollo Endosurgery, Inc. | 12/6/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                    | Failure of Implant           |
| 30067221<br>12-2022-00112 | 10/18/2022 | Malfunction | Apollo Endosurgery, Inc. | 12/6/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                                       | Failure of Implant; Vomiting |
| 30067221<br>12-2022-00110 | 10/18/2022 | Malfunction | Apollo Endosurgery, Inc. | 12/6/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation                    | Failure of Implant           |
| 30067221<br>12-2022-00115 | 10/20/2022 | Malfunction | Apollo Endosurgery, Inc. | 12/6/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Inflation Problem; Unintended Deflation | Failure of Implant           |



| Report Number             | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                          | Device Problem                 | Patient Problem                                    |
|---------------------------|------------|-------------|--------------------------------|---------------|--------------|-------------------------------------|--------------------------------|--|
| 30067221<br>12-2022-00107 | 10/21/2022 | Malfunction | Apollo Endosurgery, Inc.       | 12/6/2022     | LTI          | Orbera Intra-gastric Balloon System | Inflation Problem              | Abdominal Pain; Flatus; Nausea; Vomiting           |
| 30126389<br>28-2022-02737 | 5/11/2022  | Malfunction | Spatz Fgia Inc.                | 12/5/2022     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem              | No Clinical Signs, Symptoms or Conditions          |
| 30126389<br>28-2022-02736 | 10/19/2022 | Malfunction | Spatz Fgia Inc.                | 12/1/2022     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem              | No Clinical Signs, Symptoms or Conditions          |
| 30126389<br>28-2022-02745 | 9/6/2022   | Injury      | Journeylite Physicians         | 11/29/2022    | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility | Abdominal Pain; Dehydration; Vomiting              |
| 30126389<br>28-2022-02472 | 9/28/2022  | Injury      | True You Weight Loss (Atlanta) | 11/29/2022    | LTI          | Spatz3 Adjustable Balloon System    | Patient-Device Incompatibility | Abdominal Pain; Dehydration; Electrolyte Imbalance |
| 30126389<br>28-2022-02476 | 10/18/2021 | Malfunction | Spatz Fgia Inc.                | 11/24/2022    | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem              | No Clinical Signs, Symptoms or Conditions          |
| 30126389<br>28-2022-02529 | 3/12/2022  | Malfunction | Spatz Fgia Inc.                | 11/24/2022    | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem              | No Clinical Signs, Symptoms or Conditions          |
| 30126389<br>28-2022-02531 | 3/21/2022  | Malfunction | Spatz Fgia Inc.                | 11/24/2022    | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem              | No Clinical Signs, Symptoms or Conditions          |
| 30126389<br>28-2022-02533 | 3/26/2022  | Malfunction | Spatz Fgia Inc.                | 11/24/2022    | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem              | No Clinical Signs, Symptoms or Conditions          |

| Report Number             | Event Date | Event Type  | Manufacturer    | Date Received | Product Code | Brand Name                       | Device Problem                                      | Patient Problem   |
|---------------------------|------------|-------------|-----------------|---------------|--------------|----------------------------------|---|---|
| 30126389<br>28-2022-02660 | 7/20/2022  | Malfunction | Spatz Fgia Inc. | 11/24/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | Vomiting  |
| 30126389<br>28-2022-02707 | 9/29/2022  | Malfunction | Spatz Fgia Inc. | 11/24/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions                     |
| 30126389<br>28-2022-02733 | 11/6/2022  | Malfunction | Spatz Fgia Inc. | 11/23/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device | Abdominal Pain  |
| 30126389<br>28-2022-02728 | NR         | Malfunction | Spatz Fgia Inc. | 11/22/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions                     |
| 30126389<br>28-2022-02723 | 8/15/2022  | Malfunction | Spatz Fgia Inc. | 11/21/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device | Abdominal Pain  |
| 30126389<br>28-2022-02724 | 10/24/2022 | Malfunction | Spatz Fgia Inc. | 11/21/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device | Abdominal Pain  |
| 30126389<br>28-2022-02693 | 5/31/2022  | Injury      | Spatz Fgia Inc. | 11/17/2022    | LTI          | Spatz3 Adjustable Balloon System | Patient-Device Incompatibility                      | Dehydration; Nausea; Vomiting                                 |
| 30126389<br>28-2022-02690 | 6/16/2022  | Injury      | Spatz Fgia Inc. | 11/17/2022    | LTI          | Spatz3 Adjustable Balloon System | Patient-Device Incompatibility                      | Abdominal Pain; Diarrhea; Nausea; Vomiting; Syncope/ Fainting |

| Report Number             | Event Date | Event Type  | Manufacturer    | Date Received | Product Code | Brand Name                       | Device Problem   | Patient Problem                                     |
|---------------------------|------------|-------------|-----------------|---------------|--------------|----------------------------------|--|---|
| 30126389<br>28-2022-02695 | 8/23/2022  | Injury      | Spatz Fgia Inc. | 11/17/2022    | LTI          | Spatz3 Adjustable Balloon System | Adverse Event Without Identified Device or Use Problem | Abdominal Pain; Pain; Aspiration Pneumonitis; Cough |
| 30126389<br>28-2022-02681 | 8/28/2022  | Injury      | Spatz Fgia Inc. | 11/17/2022    | LTI          | Spatz3 Adjustable Balloon System | Patient-Device Incompatibility                         | Abdominal Pain; Nausea; Vomiting                    |
| 30126389<br>28-2022-02718 | 8/29/2022  | Malfunction | Spatz Fgia Inc. | 11/17/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions           |
| 30126389<br>28-2022-02734 | 9/10/2022  | Malfunction | Spatz Fgia Inc. | 11/17/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions           |
| 30126389<br>28-2022-02720 | 9/14/2022  | Malfunction | Spatz Fgia Inc. | 11/17/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions           |
| 30126389<br>28-2022-02719 | 9/23/2022  | Malfunction | Spatz Fgia Inc. | 11/17/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions           |
| 30126389<br>28-2022-02721 | 10/10/2022 | Malfunction | Spatz Fgia Inc. | 11/17/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions           |
| 30126389<br>28-2022-02715 | 9/19/2022  | Malfunction | Spatz Fgia Inc. | 11/16/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions           |
| 30126389<br>28-2022-02714 | 10/13/2022 | Malfunction | Spatz Fgia Inc. | 11/16/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions           |

| Report Number             | Event Date | Event Type  | Manufacturer    | Date Received | Product Code | Brand Name                       | Device Problem                                      | Patient Problem                                     |
|---------------------------|------------|-------------|-----------------|---------------|--------------|----------------------------------|---|---|
| 30126389<br>28-2022-02732 | 9/30/2022  | Malfunction | Spatz Fgia Inc. | 11/10/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration                        | Vomiting; No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02713 | 10/3/2022  | Injury      | Spatz Fgia Inc. | 11/10/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration                        | Abdominal Pain                                      |
| 30126389<br>28-2022-02740 | 10/31/2022 | Malfunction | Spatz Fgia Inc. | 11/10/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration                        | No Clinical Signs, Symptoms or Conditions           |
| 30126389<br>28-2022-02711 | 8/2/2022   | Malfunction | Spatz Fgia Inc. | 11/4/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions           |
| 30126389<br>28-2022-02705 | 9/15/2022  | Malfunction | Spatz Fgia Inc. | 11/3/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions           |
| 30126389<br>28-2022-02710 | 8/6/2022   | Malfunction | Spatz Fgia Inc. | 11/2/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device | Abdominal Pain                                      |
| 30126389<br>28-2022-02709 | 9/23/2022  | Malfunction | Spatz Fgia Inc. | 11/2/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration                        | No Clinical Signs, Symptoms or Conditions           |
| 30126389<br>28-2022-02704 | 9/19/2022  | Malfunction | Spatz Fgia Inc. | 10/26/2022    | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions           |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2022-00102 | 2/28/2020  | Malfunction | Apollo Endosurgery, Inc. | 10/20/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                    | Increased Appetite                        |
| 30067221<br>12-2022-00103 | 10/1/2020  | Malfunction | Apollo Endosurgery, Inc. | 10/20/2022    | LTI          | Orbera365 Intra gastric Balloon System | Migration or Expulsion of Device; Migration                | Abdominal Pain; Perforation               |
| 30067221<br>12-2022-00106 | 11/20/2020 | Malfunction | Apollo Endosurgery, Inc. | 10/20/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Abdominal Pain; Failure of Implant        |
| 30067221<br>12-0220-00104 | 1/15/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/20/2022    | LTI          | Orbera365 Intra gastric Balloon System | Leak/ Splash   | Abdominal Pain                            |
| 30067221<br>12-2022-00105 | 3/28/2022  | Malfunction | Apollo Endosurgery, Inc. | 10/20/2022    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2022-02697 | 8/3/2022   | Malfunction | Spatz Fgia Inc.          | 10/13/2022    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02698 | 8/12/2022  | Malfunction | Spatz Fgia Inc.          | 10/13/2022    | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device         | Abdominal Pain                            |
| 30126389<br>28-2022-02700 | 8/19/2022  | Malfunction | Spatz Fgia Inc.          | 10/13/2022    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2022-02696 | 8/24/2022  | Malfunction | Spatz Fgia Inc.          | 10/13/2022    | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device     | Abdominal Pain                            |
| 30126389<br>28-2022-02689 | 9/14/2022  | Malfunction | Spatz Fgia Inc.          | 10/13/2022    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02699 | 9/15/2022  | Malfunction | Spatz Fgia Inc.          | 10/13/2022    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02686 | 5/30/2022  | Malfunction | Spatz Fgia Inc.          | 10/6/2022     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device     | Abdominal Pain                            |
| 30126389<br>28-2022-02687 | 8/30/2022  | Malfunction | Spatz Fgia Inc.          | 10/6/2022     | LTI          | Spatz3 Adjustable Balloon System       | Air/ Gas in Device                                     | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02683 | NR         | Malfunction | Spatz Fgia Inc.          | 10/6/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00089 | 8/10/2022  | Malfunction | Apollo Endosurgery, Inc. | 10/5/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Leak/ Splash; Unintended Deflation  | Gastrointestinal Regurgitation            |
| 30067221<br>12-2022-00098 | 6/9/2022   | Malfunction | Apollo Endosurgery, Inc. | 10/4/2022     | LTI          | Orbera365 Intra gastric Balloon System | Adverse Event Without Identified Device or Use Problem | Abdominal Pain                            |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2022-00097 | 8/12/2022  | Malfunction | Apollo Endosurgery, Inc. | 10/4/2022     | LTI          | Orbera Intra gastric Balloon System    | Leak/ Splash   | Failure of Implant                        |
| 30067221<br>12-2022-00093 | 8/20/2022  | Malfunction | Apollo Endosurgery, Inc  | 10/4/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                    | Increased Appetite                        |
| 30067221<br>12-2022-00090 | 8/16/2022  | Malfunction | Apollo Endosurgery, Inc. | 9/30/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2022-02677 | 8/24/2022  | Malfunction | Spatz Fgia Inc.          | 9/29/2022     | LTI          | Spatz3 Adjustable Balloon System       | Air/ Gas in Device   | Abdominal Pain; Vomiting                  |
| 30126389<br>28-2022-02680 | NR         | Malfunction | Spatz Fgia Inc.          | 9/29/2022     | LTI          | Spatz3 Adjustable Balloon System       | Migration or Expulsion of Device                           | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02675 | 7/21/2022  | Malfunction | Spatz Fgia Inc.          | 9/28/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02673 | 8/15/2022  | Malfunction | Spatz Fgia Inc.          | 9/28/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02672 | 8/10/2022  | Malfunction | Spatz Fgia Inc.          | 9/25/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                                     |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2022-02671 | 8/26/2022  | Malfunction | Spatz Fgia Inc.          | 9/23/2022     | LTI          | Spatz3 Adjustable Balloon System       | Air/ Gas in Device   | Abdominal Pain                                      |
| 30126389<br>28-2022-02670 | 8/25/2022  | Malfunction | Spatz Fgia Inc.          | 9/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration                               | Abdominal Pain; Vomiting                            |
| 30126389<br>28-2022-02664 | 7/30/2022  | Injury      | Spatz Fgia Inc.          | 9/12/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration or Expulsion of Device        | Abdominal Pain                                      |
| 30067221<br>12-2022-00091 | 4/2/2022   | Malfunction | Apollo Endosurgery, Inc. | 9/7/2022      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                                  |
| 30067221<br>12-2022-00092 | 8/2/2022   | Malfunction | Apollo Endosurgery, Inc. | 9/7/2022      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                    | Failure of Implant                                  |
| 30126389<br>28-2022-02668 | 8/16/2022  | Injury      | Spatz Fgia Inc.          | 9/7/2022      | LTI          | Spatz3 Adjustable Balloon System       | Adverse Event Without Identified Device or Use Problem     | Abdominal Pain                                      |
| 30067221<br>12-2022-00079 | 6/9/2022   | Malfunction | Apollo Endosurgery, Inc. | 8/30/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                                  |
| 30067221<br>12-2022-00084 | NR         | Malfunction | Apollo Endosurgery, Inc. | 8/30/2022     | LTI          | Orbera Intra gastric Balloon           | Fluid/ Blood Leak  | Fever; Hemorrhage/ Bleeding; Nausea; Pain; Vomiting |



| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2022-00078 | NR         | Malfunction | Apollo Endosurgery, Inc. | 8/30/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2022-02655 | 7/21/2022  | Malfunction | Spatz Fgia Inc.          | 8/28/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00082 | 7/6/2022   | Malfunction | Apollo Endosurgery, Inc. | 8/23/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                    | Failure of Implant                        |
| 30067221<br>12-2022-00076 | 7/15/2022  | Malfunction | Apollo Endosurgery, Inc. | 8/23/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem  | Failure of Implant                        |
| 30067221<br>12-2022-00081 | 7/27/2022  | Malfunction | Apollo Endosurgery, Inc. | 8/23/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                    | Abdominal Pain                            |
| 30126389<br>28-2022-02669 | 8/15/2022  | Malfunction | Spatz Fgia Inc.          | 8/21/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02662 | 7/23/2022  | Malfunction | Spatz Fgia Inc.          | 8/18/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02657 | 7/13/2022  | Malfunction | Spatz Fgia Inc.          | 8/16/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                         | Device Problem    | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|------------------------------------|-------------------|---|
| 30067221<br>12-2022-00080 | 7/20/2022  | Malfunction | Apollo Endosurgery, Inc. | 8/16/2022     | LTI          | Orbera IntraGastric Balloon System | Inflation Problem | Failure of Implant                        |
| 30126389<br>28-2022-02658 | 7/25/2022  | Malfunction | Spatz Fgia Inc.          | 8/15/2022     | LTI          | Spatz3 Adjustable Balloon System   | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02661 | 7/25/2022  | Malfunction | Spatz Fgia Inc.          | 8/15/2022     | LTI          | Spatz3 Adjustable Balloon System   | Migration         | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02652 | 7/13/2022  | Malfunction | Spatz Fgia Inc.          | 8/11/2022     | LTI          | Spatz3 Adjustable Balloon System   | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02641 | 1/27/2022  | Malfunction | Spatz Fgia Inc.          | 8/3/2022      | LTI          | Spatz3 Adjustable Balloon System   | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02644 | 2/21/2022  | Malfunction | Spatz Fgia Inc.          | 8/3/2022      | LTI          | Spatz3 Adjustable Balloon System   | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02645 | 4/8/2022   | Malfunction | Spatz Fgia Inc.          | 8/3/2022      | LTI          | Spatz3 Adjustable Balloon System   | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02643 | 4/14/2022  | Malfunction | Spatz Fgia Inc.          | 8/3/2022      | LTI          | Spatz3 Adjustable Balloon System   | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02642 | 6/20/2022  | Malfunction | Spatz Fgia Inc.          | 8/3/2022      | LTI          | Spatz3 Adjustable Balloon System   | Deflation Problem | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem  |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|--|
| 30067221<br>12-2022-00071 | 6/27/2022  | Malfunction | Apollo Endosurgery, Inc. | 8/3/2022      | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Appropriate Clinical Signs, Symptoms, Conditions Term/Code Not Available |
| 30126389<br>28-2022-02659 | 7/25/2022  | Malfunction | Spatz Fgia Inc.          | 8/2/2022      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions                                |
| 30067221<br>12-2022-00077 | NR         | Malfunction | Apollo Endosurgery, Inc. | 8/1/2022      | LTI          | Orbera Intra-gastric Balloon System    | Inflation Problem   | Abdominal Pain; Nausea; Vomiting   |
| 30126389<br>28-2022-02656 | 6/26/2022  | Malfunction | Spatz Fgia Inc.          | 7/30/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration                              | No Clinical Signs, Symptoms or Conditions                                |
| 30126389<br>28-2022-02634 | 6/17/2022  | Malfunction | Spatz Fgia Inc.          | 7/28/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions                                |
| 30067221<br>12-2022-00075 | 5/10/2022  | Malfunction | Apollo Endosurgery, Inc. | 7/26/2022     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant   |
| 30067221<br>12-2022-00070 | 5/26/2022  | Malfunction | Apollo Endosurgery, Inc. | 7/26/2022     | LTI          | Orbera Intra-gastric Balloon System    | Inflation Problem   | Abdominal Pain; Failure of Implant; Pain                                 |
| 30067221<br>12-2022-00074 | 6/17/2022  | Malfunction | Apollo Endosurgery, Inc. | 7/26/2022     | LTI          | Orbera365 Intra-gastric Balloon System | Expulsion   | Hemorrhage/ Bleeding; Gastrointestinal Hemorrhage                        |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2022-02647 | 7/23/2021  | Malfunction | Spatz Fgia Inc.          | 7/25/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02651 | 6/13/2022  | Malfunction | Spatz Fgia Inc.          | 7/25/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00066 | 5/13/2022  | Malfunction | Apollo Endosurgery, Inc. | 7/21/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Abdominal Pain                            |
| 30126389<br>28-2022-02650 | 6/13/2022  | Malfunction | Spatz Fgia Inc.          | 7/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | Abdominal Pain                            |
| 30126389<br>28-2022-02615 | 6/14/2022  | Malfunction | Spatz Fgia Inc.          | 7/15/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00061 | 4/22/2022  | Malfunction | Apollo Endosurgery, Inc. | 7/13/2022     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem  | Failure of Implant                        |
| 30067221<br>12-2022-00060 | 6/2/2022   | Malfunction | Apollo Endosurgery, Inc. | 7/13/2022     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem  | Abdominal Pain; Pain                      |
| 30067221<br>12-2022-00073 | 6/6/2022   | Malfunction | Apollo Endosurgery, Inc. | 7/13/2022     | LTI          | Orbera Intra gastric Balloon System    | Infusion or Flow Problem                                   | Abdominal Pain; Nausea; Vomiting          |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2022-00063 | 6/10/2022  | Malfunction | Apollo Endosurgery, Inc. | 7/13/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Vomiting  |
| 30126389<br>28-2022-02646 | 6/13/2022  | Malfunction | Spatz Fgia Inc.          | 7/13/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration                               | No Clinical Signs, Symptoms or Conditions   |
| 30126389<br>28-2022-02648 | 7/5/2022   | Malfunction | Spatz Fgia Inc.          | 7/13/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions   |
| 30067221<br>12-2022-00072 | NR         | Malfunction | Apollo Endosurgery, Inc. | 7/8/2022      | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem  | Insufficient Information  |
| 30126389<br>28-2022-02632 | 6/28/2022  | Malfunction | Spatz Fgia Inc.          | 7/7/2022      | LTI          | Spatz3 Adjustable Balloon System       | Migration or Expulsion of Device; Migration                | Abdominal Pain  |
| 30126389<br>28-2022-02639 | 7/4/2022   | Malfunction | Spatz Fgia Inc.          | 7/7/2022      | LTI          | Spatz3 Adjustable Balloon System       | Infusion or Flow Problem; Air/ Gas in Device               | Abdominal Pain  |
| 30067221<br>12-2022-00068 | 6/23/2022  | Malfunction | Apollo Endosurgery, Inc. | 7/6/2022      | LTI          | Orbera365 Intra gastric Balloon System | Inflation Problem; Mechanical Problem                      | Abdominal Pain; Vomiting; Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30126389<br>28-2022-02622 | 6/24/2022  | Malfunction | Spatz Fgia Inc.          | 7/5/2022      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions   |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2022-02633 | 6/29/2022  | Malfunction | Spatz Fgia Inc.          | 7/5/2022      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02621 | 5/28/2021  | Malfunction | Spatz Fgia Inc.          | 7/4/2022      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00067 | 11/18/2021 | Malfunction | Apollo Endosurgery, Inc. | 7/1/2022      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Insufficient Information                  |
| 30126389<br>28-2022-02612 | 2/12/2022  | Malfunction | Spatz Fgia Inc.          | 6/28/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02607 | 6/9/2022   | Malfunction | Spatz Fgia Inc.          | 6/28/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02618 | 6/23/2022  | Malfunction | Spatz Fgia Inc.          | 6/28/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02609 | 5/10/2022  | Malfunction | Spatz Fgia Inc.          | 6/23/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02600 | 5/24/2022  | Malfunction | Spatz Fgia Inc.          | 6/23/2022     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device         | Abdominal Pain; Vomiting                  |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2022-02613 | 5/26/2022  | Malfunction | Spatz Fgia Inc.          | 6/23/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02614 | 6/20/2022  | Malfunction | Spatz Fgia Inc.          | 6/23/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02598 | 4/10/2022  | Malfunction | Spatz Fgia Inc.          | 6/21/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00064 | 4/11/2022  | Malfunction | Apollo Endosurgery, Inc. | 6/21/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2022-02608 | 5/20/2022  | Malfunction | Spatz Fgia Inc.          | 6/21/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02599 | 5/21/2022  | Malfunction | Spatz Fgia Inc.          | 6/21/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02601 | 5/21/2022  | Malfunction | Spatz Fgia Inc.          | 6/21/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02606 | 6/6/2022   | Malfunction | Spatz Fgia Inc.          | 6/21/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem  |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|--|
| 30067221<br>12-2022-00054 | 4/6/2022   | Malfunction | Apollo Endosurgery, Inc. | 6/20/2022     | LTI          | Orbera Intra-gastric Balloon System    | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant   |
| 30067221<br>12-2022-00051 | 5/5/2022   | Malfunction | Apollo Endosurgery, Inc  | 6/20/2022     | LTI          | Orbera365 Intra-gastric Balloon System | Inflation Problem   | Failure of Implant   |
| 30067221<br>12-2022-00055 | 5/16/2022  | Malfunction | Apollo Endosurgery, Inc. | 6/20/2022     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant   |
| 30067221<br>12-2022-00058 | 5/22/2022  | Malfunction | Apollo Endosurgery, Inc. | 6/20/2022     | LTI          | Orbera365 Intra-gastric Balloon System | Inflation Problem   | Failure of Implant; Vomiting   |
| 30067221<br>12-2022-00057 | 5/22/2022  | Malfunction | Apollo Endosurgery, Inc. | 6/20/2022     | LTI          | Orbera365 Intra-gastric Balloon System | Inflation Problem   | Failure of Implant; Vomiting   |
| 30067221<br>12-2022-00059 | 5/24/2022  | Malfunction | Apollo Endosurgery, Inc. | 6/20/2022     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Unintended Deflation                   | Failure of Implant   |
| 30067221<br>12-2022-00062 | 5/31/2022  | Malfunction | Apollo Endosurgery, Inc  | 6/20/2022     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem   | Appropriate Clinical Signs, Symptoms, Conditions Term/Code Not Available |
| 30126389<br>28-2022-02604 | 5/27/2022  | Malfunction | Spatz Fgia Inc.          | 6/17/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration                              | Abdominal Pain   |



| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                      | Device Problem   | Patient Problem   |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|---|--|---|
| 30126389<br>28-2022-<br>02605 | 6/3/2022   | Malfunction | Spatz Fgia<br>Inc.             | 6/17/2022     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System       | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions                 |
| 30067221<br>12-2022-<br>00056 | 6/5/2022   | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/17/2022     | LTI          | Orbera<br>Intragastric<br>Balloon<br>System     | Insufficient<br>Information  | Abdominal Pain;<br>Nausea; Vomiting;<br>Abdominal<br>Distention |
| 30067221<br>12-2022-<br>00046 | 5/10/2022  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/15/2022     | LTI          | BIB System<br>Intragastric<br>Balloon<br>System | Inflation Problem  | Failure of Implant  |
| 30067221<br>12-2022-<br>00049 | 5/18/2022  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/15/2022     | LTI          | Orbera<br>Intragastric<br>Balloon<br>System     | Deflation<br>Problem; Fluid/<br>Blood Leak;<br>Unintended<br>Deflation | Failure of Implant  |
| 30067221<br>12-2022-<br>00053 | NR         | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/15/2022     | LTI          | Orbera<br>Intragastric<br>Balloon<br>System     | Insufficient<br>Information  | Hemorrhage/<br>Bleeding   |
| 30067221<br>12-2022-<br>00045 | 4/22/2022  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/14/2022     | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System  | Deflation<br>Problem; Fluid/<br>Blood Leak;<br>Unintended<br>Deflation | Failure of Implant  |
| 30126389<br>28-2022-<br>02591 | 4/16/2022  | Malfunction | Spatz Fgia<br>Inc.             | 6/12/2022     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System       | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions                 |
| 30126389<br>28-2022-<br>02596 | 4/19/2022  | Malfunction | Spatz Fgia<br>Inc.             | 6/12/2022     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System       | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions                 |

| Report Number                 | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|-------------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2022-<br>02595 | 4/26/2022  | Malfunction | Spatz Fgia Inc.          | 6/12/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-<br>02594 | 5/28/2022  | Malfunction | Spatz Fgia Inc.          | 6/12/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | Abdominal Pain                            |
| 30126389<br>28-2022-<br>02586 | 12/28/2021 | Malfunction | Spatz Fgia Inc.          | 6/9/2022      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-<br>00050 | 4/30/2022  | Malfunction | Apollo Endosurgery, Inc. | 6/9/2022      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2022-<br>00048 | 5/5/2022   | Malfunction | Apollo Endosurgery, Inc. | 6/9/2022      | LTI          | Orbera Intra gastric Balloon System    | Fluid/ Blood Leak  | Failure of Implant                        |
| 30126389<br>28-2022-<br>02585 | 5/11/2022  | Injury      | Spatz Fgia Inc.          | 6/9/2022      | LTI          | Spatz3 Adjustable Balloon System       | Adverse Event Without Identified Device or Use Problem     | Aspiration/ Inhalation                    |
| 30126389<br>28-2022-<br>02580 | 5/10/2022  | Malfunction | Spatz Fgia Inc.          | 6/8/2022      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-<br>02587 | 5/16/2022  | Malfunction | Spatz Fgia Inc.          | 6/2/2022      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                        | Device Problem               | Patient Problem  |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|-----------------------------------|------------------------------|--|
| 30126389<br>28-2022-02582 | 10/8/2021  | Malfunction | Spatz Fgia Inc.          | 5/31/2022     | LTI          | Spatz3 Adjustable Balloon System  | Deflation Problem            | No Clinical Signs, Symptoms or Conditions                                    |
| 30126389<br>28-2022-02583 | 11/12/2021 | Malfunction | Spatz Fgia Inc.          | 5/31/2022     | LTI          | Spatz3 Adjustable Balloon System  | Deflation Problem            | No Clinical Signs, Symptoms or Conditions                                    |
| 30126389<br>28-2022-02584 | 1/15/2022  | Malfunction | Spatz Fgia Inc.          | 5/31/2022     | LTI          | Spatz3 Adjustable Balloon System  | Deflation Problem; Migration | No Clinical Signs, Symptoms or Conditions                                    |
| 30126389<br>28-2022-02579 | 3/22/2022  | Malfunction | Spatz Fgia Inc.          | 5/29/2022     | LTI          | Spatz3 Adjustable Balloon System  | Air/ Gas in Device           | No Clinical Signs, Symptoms or Conditions                                    |
| 30126389<br>28-2022-02581 | 5/5/2022   | Malfunction | Spatz Fgia Inc.          | 5/29/2022     | LTI          | Spatz3 Adjustable Balloon System  | Deflation Problem; Migration | No Clinical Signs, Symptoms or Conditions                                    |
| 30126389<br>28-2022-02570 | 4/26/2022  | Malfunction | Spatz Fgia Inc.          | 5/26/2022     | LTI          | Spatz3 Adjustable Balloon System  | Migration                    | Abdominal Pain   |
| 30126389<br>28-2022-02575 | 5/2/2022   | Malfunction | Spatz Fgia Inc.          | 5/26/2022     | LTI          | Spatz3 Adjustable Balloon System  | Deflation Problem            | No Clinical Signs, Symptoms or Conditions                                    |
| 30067221<br>12-2022-00043 | NR         | Malfunction | Apollo Endosurgery, Inc. | 5/26/2022     | LTI          | Orbera Intragastic Balloon System | Material Perforation         | Abdominal Pain; Pyrosis/ Heartburn; Pain; Vomiting; Perforation of Esophagus |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                                   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2022-00044 | NR         | Malfunction | Apollo Endosurgery, Inc. | 5/26/2022     | LTI          | Orbera Intra gastric Balloon System    | Adverse Event Without Identified Device or Use Problem | Dehydration; Fever; Nausea; Pain; Vomiting; Ulcer |
| 30126389<br>28-2022-02571 | 4/29/2022  | Malfunction | Spatz Fgia Inc.          | 5/24/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions         |
| 30126389<br>28-2022-02573 | NR         | Malfunction | Spatz Fgia Inc.          | 5/24/2022     | LTI          | Spatz3 Adjustable Balloon System       | Migration  | No Clinical Signs, Symptoms or Conditions         |
| 30126389<br>28-2022-02564 | 1/28/2022  | Malfunction | Spatz Fgia Inc.          | 5/21/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions         |
| 30126389<br>28-2022-02561 | 3/14/2022  | Malfunction | Spatz Fgia Inc.          | 5/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions         |
| 30067221<br>12-2022-00040 | 3/21/2022  | Malfunction | Apollo Endosurgery, Inc. | 5/20/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                | Failure of Implant                                |
| 30126389<br>28-2022-02560 | 3/3/2022   | Malfunction | Spatz Fgia Inc.          | 5/19/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions         |
| 30126389<br>28-2022-02558 | 3/16/2022  | Malfunction | Spatz Fgia Inc.          | 5/19/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions         |
| 30126389<br>28-2022-02565 | 12/28/2021 | Malfunction | Spatz Fgia Inc.          | 5/18/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                      | No Clinical Signs, Symptoms or Conditions         |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                           | Device Problem                          | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--------------------------------------|---|---|
| 30126389<br>28-2022-02563 | 4/6/2021   | Malfunction | Spatz Fgia Inc.          | 5/17/2022     | LTI          | Spatz3 Adjustable Balloon System     | Deflation Problem                       | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02559 | 3/23/2022  | Malfunction | Spatz Fgia Inc.          | 5/16/2022     | LTI          | Spatz3 Adjustable Balloon System     | Deflation Problem                       | Abdominal Pain                            |
| 30126389<br>28-2022-02562 | 3/9/2022   | Malfunction | Spatz Fgia Inc.          | 5/15/2022     | LTI          | Spatz3 Adjustable Balloon System     | Deflation Problem                       | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02557 | 3/25/2022  | Malfunction | Spatz Fgia Inc.          | 5/15/2022     | LTI          | Spatz3 Adjustable Balloon System     | Deflation Problem                       | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02554 | 4/14/2022  | Malfunction | Spatz Fgia Inc.          | 5/15/2022     | LTI          | Spatz3 Adjustable Balloon System     | Deflation Problem                       | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02566 | 12/20/2021 | Malfunction | Spatz Fgia Inc.          | 5/12/2022     | LTI          | Spatz3 Adjustable Balloon System     | Deflation Problem                       | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02555 | 4/12/2022  | Malfunction | Spatz Fgia Inc.          | 5/10/2022     | LTI          | Spatz3 Adjustable Balloon System     | Deflation Problem                       | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00041 | 4/13/2022  | Malfunction | Apollo Endosurgery, Inc. | 5/10/2022     | LTI          | Orbera365 Intragastic Balloon System | Deflation Problem; Unintended Deflation | Vomiting                                  |
| 30126389<br>28-2022-02552 | 3/8/2022   | Malfunction | Spatz Fgia Inc.          | 5/9/2022      | LTI          | Spatz3 Adjustable Balloon System     | Deflation Problem                       | No Clinical Signs, Symptoms or Conditions |

| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                     | Device Problem  | Patient Problem                                 |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|--|---|---|
| 30126389<br>28-2022-<br>02553 | 4/6/2022   | Malfunction | Spatz Fgia<br>Inc.             | 5/8/2022      | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem  | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30067221<br>12-2022-<br>00039 | 3/26/2022  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 5/7/2022      | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem;<br>Unintended<br>Deflation                | Failure of Implant                              |
| 30126389<br>28-2022-<br>02551 | 4/9/2022   | Malfunction | Spatz Fgia<br>Inc.             | 5/4/2022      | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem  | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30126389<br>28-2022-<br>02547 | 1/26/2022  | Malfunction | Spatz Fgia<br>Inc.             | 5/3/2022      | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem  | Abdominal Pain                                  |
| 30126389<br>28-2022-<br>02548 | 3/23/2022  | Malfunction | Spatz Fgia<br>Inc.             | 5/3/2022      | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem  | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30126389<br>28-2022-<br>02549 | 4/1/2022   | Malfunction | Spatz Fgia<br>Inc.             | 5/3/2022      | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Air/ Gas in<br>Device   | Abdominal Pain                                  |
| 30126389<br>28-2022-<br>02526 | NR         | Malfunction | Spatz Fgia<br>Inc.             | 5/3/2022      | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem  | No Clinical Signs,<br>Symptoms or<br>Conditions |
| 30067221<br>12-2022-<br>00036 | 2/15/2022  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 4/28/2022     | LTI          | Orbera®<br>Intragastric<br>Balloon<br>System   | Adverse Event<br>Without<br>Identified Device<br>or Use Problem | Abdominal Pain;<br>Vomiting; Ulcer              |
| 30067221<br>12-2022-<br>00035 | 2/24/2022  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 4/28/2022     | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem;<br>Unintended<br>Deflation                | Failure of Implant                              |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                              | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|--|---|
| 30067221<br>22-0220-00038 | 3/30/2022  | Malfunction | Apollo Endosurgery, Inc. | 4/28/2022     | LTI          | BIB System Intra gastric Balloon System | Inflation Problem  | Abdominal Pain                            |
| 30126389<br>28-2022-02545 | 12/27/2021 | Malfunction | Spatz Fgia Inc.          | 4/26/2022     | LTI          | Spatz3 Adjustable Balloon System        | Deflation Problem  | Vomiting                                  |
| 30126389<br>28-2022-02544 | 3/2/2022   | Malfunction | Spatz Fgia Inc.          | 4/26/2022     | LTI          | Spatz3 Adjustable Balloon System        | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02537 | 3/17/2022  | Malfunction | Spatz Fgia Inc.          | 4/26/2022     | LTI          | Spatz3 Adjustable Balloon System        | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02535 | 3/24/2022  | Injury      | Spatz Fgia Inc.          | 4/26/2022     | LTI          | Spatz3 Adjustable Balloon System        | Patient Device Interaction Problem   | Stomach Ulceration                        |
| 30126389<br>28-2022-02536 | 3/25/2022  | Injury      | Spatz Fgia Inc.          | 4/26/2022     | LTI          | Spatz3 Adjustable Balloon System        | Adverse Event Without Identified Device or Use Problem; Patient Device Interaction Problem | Stomach Ulceration                        |
| 30126389<br>28-2022-02541 | 4/4/2022   | Malfunction | Spatz Fgia Inc.          | 4/26/2022     | LTI          | Spatz3 Adjustable Balloon System        | Deflation Problem; Output Problem  | Insufficient Information                  |
| 30126389<br>28-2022-02540 | NR         | Malfunction | Spatz Fgia Inc.          | 4/25/2022     | LTI          | Spatz3 Adjustable Balloon System        | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer    | Date Received | Product Code | Brand Name                       | Device Problem     | Patient Problem                           |
|---------------------------|------------|-------------|-----------------|---------------|--------------|----------------------------------|--------------------|---|
| 30126389<br>28-2022-02534 | 3/24/2022  | Malfunction | Spatz Fgia Inc. | 4/17/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02527 | 1/20/2022  | Malfunction | Spatz Fgia Inc. | 4/13/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02528 | 2/25/2022  | Malfunction | Spatz Fgia Inc. | 4/11/2022     | LTI          | Spatz3 Adjustable Balloon System | Air/ Gas in Device | Abdominal Pain                            |
| 30126389<br>28-2022-02522 | 3/1/2022   | Malfunction | Spatz Fgia Inc. | 4/11/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | Abdominal Pain                            |
| 30126389<br>28-2022-02530 | 3/21/2022  | Malfunction | Spatz Fgia Inc. | 4/11/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02532 | 3/21/2022  | Malfunction | Spatz Fgia Inc. | 4/11/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02525 | 1/7/2022   | Malfunction | Spatz Fgia Inc. | 4/10/2022     | LTI          | Spatz3 Adjustable Balloon System | Air/ Gas in Device | Abdominal Pain                            |
| 30126389<br>28-2022-02523 | 2/1/2022   | Malfunction | Spatz Fgia Inc. | 4/10/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | Abdominal Pain                            |
| 30126389<br>28-2022-02524 | 12/27/2022 | Malfunction | Spatz Fgia Inc. | 4/10/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | Abdominal Pain                            |



| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|-------------------------------------|--|---|
| 30067221<br>12-2022-00034 | 3/9/2022   | Malfunction | Apollo Endosurgery, Inc. | 4/6/2022      | LTI          | BIB System Intra gastric Balloon    | Inflation Problem  | Nausea; Abdominal Distention              |
| 30067221<br>12-2022-00037 | 3/9/2022   | Malfunction | Apollo Endosurgery, Inc. | 4/6/2022      | LTI          | Bib® Intra gastric Balloon System   | Inflation Problem  | Nausea; Abdominal Distention              |
| 30126389<br>28-2022-02520 | 1/19/2022  | Malfunction | Spatz Fgia Inc.          | 4/5/2022      | LTI          | Spatz3 Adjustable Balloon System    | Migration or Expulsion of Device                           | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00032 | 2/20/2022  | Malfunction | Apollo Endosurgery, Inc. | 3/28/2022     | LTI          | Orbera Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2022-02509 | 12/21/2021 | Malfunction | Spatz Fgia Inc.          | 3/23/2022     | LTI          | Spatz3 Adjustable Balloon System    | Deflation Problem  | Abdominal Pain                            |
| 30126389<br>28-2022-02501 | 2/21/2022  | Malfunction | Spatz Fgia Inc.          | 3/23/2022     | LTI          | Spatz3 Adjustable Balloon System    | Fluid/ Blood Leak  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02510 | 1/7/2022   | Malfunction | Spatz Fgia Inc.          | 3/21/2022     | LTI          | Spatz3 Adjustable Balloon System    | Fluid/ Blood Leak  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02511 | 2/4/2022   | Malfunction | Sptaz Fgia Inc.          | 3/21/2022     | LTI          | Spatz3 Adjustable Balloon System    | Fluid/ Blood Leak  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2022-02513 | 2/16/2022  | Malfunction | Spatz Fgia Inc.          | 3/21/2022     | LTI          | Spatz3 Adjsutable Balloon System       | Deflation Problem; Migration or Expulsion of Device        | Abdominal Pain                            |
| 30126389<br>28-2022-02512 | 2/16/2022  | Malfunction | Spatz Fgia Inc.          | 3/21/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02507 | 2/28/2022  | Malfunction | Spatz Fgia Inc.          | 3/21/2022     | LTI          | Spatz3 Adjustable Balloon System       | Tear, Rip or Hole in Device Packaging                      | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00028 | 1/19/2022  | Malfunction | Apollo Endosurgery, Inc. | 3/15/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation |   |
| 30067221<br>12-2022-00029 | 1/21/2022  | Malfunction | Apollo Endosurgery, Inc. | 3/15/2022     | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant; Vomiting              |
| 30067221<br>12-2022-00026 | 2/7/2022   | Malfunction | Apollo Endosurgery, Inc. | 3/8/2022      | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30126389<br>28-2022-24445 | 12/7/2021  | Malfunction | Spatz Fgia Inc.          | 3/7/2022      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02487 | 12/17/2021 | Malfunction | Spatz Fgia Inc.          | 3/7/2022      | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                              | Device Problem   | Patient Problem  |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|--|--|
| 30126389<br>28-2022-02486 | 12/20/2021 | Malfunction | Spatz Fgia Inc.          | 3/7/2022      | LTI          | Spatz3 Adjustable Balloon System        | Deflation Problem  | Abdominal Pain   |
| 30126389<br>28-2022-02506 | 2/23/2022  | Malfunction | Spatz Fgia Inc.          | 3/7/2022      | LTI          | Spatz3 Adjustable Balloon System        | Deflation Problem  | No Clinical Signs, Symptoms or Conditions                                      |
| 30067221<br>12-2022-00025 | 11/3/2021  | Malfunction | Apollo Endosurgery, Inc. | 3/4/2022      | LTI          | BIB System Intra gastric Balloon System | Appropriate Term/ Code Not Available                       | Abdominal Pain; Failure of Implant; Obstruction/ Occlusion                     |
| 30067221<br>12-2022-00022 | 11/10/2021 | Malfunction | Apollo Endosurgery, Inc. | 3/3/2022      | LTI          | Orbera Intra gastric Balloon System     | Insufficient Information                                   | Hemorrhage/ Bleeding; Low Blood Pressure/ Hypotension; Nausea; Vomiting; Ulcer |
| 30067221<br>12-2022-00021 | 1/23/2022  | Malfunction | Apollo Endosurgery, Inc. | 3/3/2022      | LTI          | BIB System Intra gastric Balloon System | Inflation Problem  | Failure of Implant   |
| 30067221<br>12-2022-00024 | 1/26/2022  | Malfunction | Apollo Endosurgery, Inc. | 3/3/2022      | LTI          | Orbera365 Intra gastric Balloon System  | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant; Cramp(s)/ Muscle Spasm(s)                                  |
| 30067221<br>12-2022-00020 | 12/3/2021  | Malfunction | Apollo Endosurgery, Inc. | 3/2/2022      | LTI          | Orbera365 Intra gastric Balloon System  | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant   |
| 30067221<br>12-2022-00019 | 12/7/2021  | Malfunction | Apollo Endosurgery, Inc. | 3/2/2022      | LTI          | Orbera365 Intra gastric Balloon System  | Deflation Problem; Unintended Deflation                    | Failure of Implant; Vomiting   |

| Report Number             | Event Date | Event Type  | Manufacturer                         | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                                |
|---------------------------|------------|-------------|--------------------------------------|---------------|--------------|--|--|--|
| 30067221<br>12-2022-00031 | 1/31/2022  | Malfunction | Apollo Endosurgery, Inc.             | 3/2/2022      | LTI          | Orbera365 Intra gastric Balloon System | Inflation Problem  | Failure of Implant; Vomiting                   |
| 30067221<br>12-2022-00023 | 2/1/2022   | Malfunction | Apollo Endosurgery, Inc.             | 3/2/2022      | LTI          | Orbera Intra gastric Balloon System    | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Dehydration; Failure of Implant; Vomiting      |
| MW5107<br>793             | 11/9/2021  | Injury      | Apollo Endosurgery Costa Rica S.R.L. | 2/28/2022     | LTI          | Orbera Endoscopy/ Orbera               | Unexpected Therapeutic Results                             | Hemorrhage/ Bleeding; Nausea; Shock; Dizziness |
| 30126389<br>28-2021-02421 | 9/14/2021  | Malfunction | Spatz Fgia Inc.                      | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | Abdominal Pain                                 |
| 30126389<br>28-2022-02484 | 11/5/2021  | Malfunction | Spatz Fgia Inc.                      | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | Abdominal Pain                                 |
| 30126389<br>28-2022-02436 | 12/8/2021  | Malfunction | Spatz Fgia Inc.                      | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions      |
| 30126389<br>28-2022-02490 | 12/16/2021 | Malfunction | Spatz Fgia Inc.                      | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions      |
| 30126389<br>28-2022-02489 | 12/20/2021 | Malfunction | Spatz Fgia Inc.                      | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions      |

| Report Number             | Event Date | Event Type  | Manufacturer    | Date Received | Product Code | Brand Name                       | Device Problem                                     | Patient Problem                           |
|---------------------------|------------|-------------|-----------------|---------------|--------------|----------------------------------|--|---|
| 30126389<br>28-2022-02494 | 12/27/2021 | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02493 | 12/28/2021 | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02492 | 1/12/2022  | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02491 | 1/13/2022  | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02495 | 1/31/2022  | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02497 | 2/11/2022  | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                  | Abdominal Pain                            |
| 30126389<br>28-2022-02496 | 2/12/2022  | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02499 | 2/14/2022  | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Patient-Device Incompatibility; Air/ Gas in Device | Abdominal Pain                            |
| 30126389<br>28-2022-02500 | 2/16/2022  | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer    | Date Received | Product Code | Brand Name                       | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|-----------------|---------------|--------------|----------------------------------|--|---|
| 30126389<br>28-2022-02503 | 2/18/2022  | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02485 | 12/20/2022 | Malfunction | Spatz Fgia Inc. | 2/27/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | Abdominal Pain                            |
| 30126389<br>28-2022-02477 | 5/15/2021  | Malfunction | Spatz Fgia Inc. | 2/21/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02478 | 8/23/2021  | Malfunction | Spatz Fgia Inc. | 2/21/2022     | LTI          | Spatz3 Adjustable Balloon System | Patient-Device Incompatibility; Infusion or Flow Problem | Abdominal Pain                            |
| 30126389<br>28-2022-02483 | 12/7/2021  | Malfunction | Spatz Fgia Inc. | 2/21/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02482 | 12/16/2021 | Malfunction | Spatz Fgia Inc. | 2/21/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02481 | 12/27/2021 | Malfunction | Spatz Fgia Inc. | 2/21/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | Abdominal Pain                            |
| 30126389<br>28-2022-02488 | 1/12/2022  | Malfunction | Spatz Fgia Inc. | 2/21/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2022-00018 | 1/24/2022  | Malfunction | Apollo Endosurgery, Inc. | 2/17/2022     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation   | Failure of Implant                        |
| 30126389<br>28-2022-02742 | NR         | Death       | Spatz Fgia Inc.          | 2/16/2022     | LTI          | Spatz3 Adjustable Balloon System       | Insufficient Information; Patient Device Interaction Problem | Abdominal Pain; Perforation               |
| 30126389<br>28-2021-02395 | 5/19/2021  | Malfunction | Spatz Fgia Inc.          | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System       | Migration; Unintended Deflation                              | Failure of Implant                        |
| 30126389<br>28-2021-02396 | 7/30/2021  | Malfunction | Spatz Fgia Inc.          | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System       | Migration or Expulsion of Device; Unintended Deflation       | Vomiting                                  |
| 30126389<br>28-2021-02453 | 10/22/2021 | Malfunction | Spatz Fgia Inc.          | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2021-02397 | 11/10/2021 | Malfunction | Spatz Fgia Inc.          | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System       | Unintended Deflation   | Failure of Implant                        |
| 30126389<br>28-2021-02405 | 11/19/2021 | Malfunction | Spatz Fgia Inc.          | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System       | Unintended Deflation   | Failure of Implant                        |
| 30126389<br>28-2021-02434 | 11/25/2021 | Malfunction | Spatz Fgia Inc.          | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem  | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer    | Date Received | Product Code | Brand Name                       | Device Problem    | Patient Problem                           |
|---------------------------|------------|-------------|-----------------|---------------|--------------|----------------------------------|-------------------|---|
| 30126389<br>28-2021-02430 | 11/29/2021 | Malfunction | Spatz Fgia Inc. | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2021-02427 | 12/14/2021 | Malfunction | Spatz Fgia Inc. | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02475 | 12/21/2021 | Malfunction | Spatz Fgia Inc. | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02463 | 12/27/2021 | Malfunction | Spatz Fgia Inc. | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02461 | 12/29/2021 | Malfunction | Spatz Fgia Inc. | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem | Abdominal Pain                            |
| 30126389<br>28-2021-02457 | 1/8/2022   | Malfunction | Spatz Fgia Inc. | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02467 | 1/16/2022  | Malfunction | Spatz Fgia Inc. | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem | Abdominal Pain                            |
| 30126389<br>28-2022-02470 | 1/28/2022  | Malfunction | Spatz Fgia Inc. | 2/15/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem | Abdominal Pain                            |



| Report Number             | Event Date | Event Type  | Manufacturer    | Date Received | Product Code | Brand Name                       | Device Problem                                      | Patient Problem                           |
|---------------------------|------------|-------------|-----------------|---------------|--------------|----------------------------------|---|---|
| 30126328<br>89-2021-02439 | 4/9/2021   | Malfunction | Spatz Fgia Inc. | 2/6/2022      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device | Abdominal Pain                            |
| 30126389<br>28-2021-02429 | 7/23/2021  | Malfunction | Spatz Fgia Inc. | 2/6/2022      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2021-02414 | 11/12/2021 | Malfunction | Spatz Fgia Inc. | 2/6/2022      | LTI          | Spatz3 Adjustable Balloon System | Material Rupture                                    | Insufficient Information                  |
| 30126389<br>28-2022-02466 | 12/18/2021 | Malfunction | Spatz Fgia Inc. | 2/6/2022      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device | Abdominal Pain                            |
| 30126389<br>28-2022-02464 | 1/10/2022  | Malfunction | Spatz Fgia Inc. | 2/6/2022      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | Insufficient Information                  |
| 30126389<br>28-2022-02471 | 1/20/2022  | Malfunction | Spatz Fgia Inc. | 2/6/2022      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02469 | 1/21/2022  | Malfunction | Spatz Fgia Inc. | 2/6/2022      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02468 | 1/25/2022  | Malfunction | Spatz Fgia Inc. | 2/6/2022      | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | Abdominal Pain                            |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|---|
| 30067221<br>12-2022-00016 | 1/12/2022  | Malfunction | Apollo Endosurgery, Inc. | 1/26/2022     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant                        |
| 30067221<br>12-2022-00017 | 1/17/2022  | Malfunction | Apollo Endosurgery, Inc. | 1/26/2022     | LTI          | Orbera® Intra-gastric Balloon System   | Inflation Problem   | Failure of Implant; Vomiting              |
| 30126389<br>28-2021-02449 | 9/2/2021   | Malfunction | Spatz Fgia Inc.          | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem; Migration or Expulsion of Device       | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2021-02437 | 9/3/2021   | Malfunction | Spatz Fgia Inc.          | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2021-02451 | 9/24/2021  | Malfunction | Spatz Fgia Inc.          | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2021-02435 | 10/25/2021 | Malfunction | Spatz Fgia Inc.          | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2021-02428 | 10/30/2021 | Malfunction | Spatz Fgia Inc.          | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2021-02431 | 11/3/2021  | Malfunction | Spatz Fgia Inc.          | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | No Clinical Signs, Symptoms or Conditions |

| Report Number             | Event Date | Event Type  | Manufacturer    | Date Received | Product Code | Brand Name                       | Device Problem                                      | Patient Problem   |
|---------------------------|------------|-------------|-----------------|---------------|--------------|----------------------------------|---|---|
| 30126389<br>28-2021-02426 | 11/8/2021  | Malfunction | Spatz Fgia Inc. | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions                                 |
| 30126389<br>28-2021-02423 | 11/22/2021 | Malfunction | Spatz Fgia Inc. | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem; Migration or Expulsion of Device | No Clinical Signs, Symptoms or Conditions                                 |
| 30126389<br>28-2021-02425 | 11/24/2021 | Malfunction | Spatz Fgia Inc. | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30126389<br>28-2021-02432 | 11/29/2021 | Malfunction | Spatz Fgia Inc. | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions                                 |
| 30126389<br>28-2021-02433 | 11/29/2021 | Malfunction | Spatz Fgia Inc. | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions                                 |
| 30126389<br>28-2021-02443 | 12/6/2021  | Malfunction | Spatz Fgia Inc. | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions                                 |
| 30126389<br>28-2021-02454 | 12/6/2021  | Malfunction | Spatz Fgia Inc. | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions                                 |
| 30126389<br>28-2022-02460 | 12/13/2021 | Malfunction | Spatz Fgia Inc. | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System | Deflation Problem                                   | No Clinical Signs, Symptoms or Conditions                                 |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem                                     | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2022-02459 | 12/25/2021 | Malfunction | Spatz Fgia Inc.          | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2021-02447 | 1/3/2022   | Malfunction | Spatz Fgia Inc.          | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30126389<br>28-2022-02462 | 1/12/2022  | Malfunction | Spatz Fgia Inc.          | 1/20/2022     | LTI          | Spatz3 Adjustable Balloon System       | Patient-Device Incompatibility; Air/ Gas in Device | Abdominal Pain                            |
| 30067221<br>12-2022-00012 | 10/4/2021  | Malfunction | Apollo Endosurgery, Inc. | 1/19/2022     | LTI          | Orbera365 Intra gastric Balloon System | Insufficient Information                           | Vomiting                                  |
| 30067221<br>12-2022-00015 | 11/23/2021 | Malfunction | Apollo Endosurgery, Inc. | 1/19/2022     | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak; Unintended Deflation            | Failure of Implant                        |
| 30067221<br>12-2022-00014 | 12/16/2021 | Malfunction | Apollo Endosurgery, Inc. | 1/19/2022     | LTI          | Orbera® Intra gastric Balloon System   | Insufficient Information                           | Dehydration; Nausea; Pain; Vomiting       |
| 30067221<br>12-2022-00013 | 1/11/2022  | Malfunction | Apollo Endosurgery, Inc. | 1/19/2022     | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak; Unintended Deflation            | Failure of Implant                        |
| 30126389<br>28-2021-02422 | 11/11/2021 | Malfunction | Spatz Fgia Inc.          | 1/18/2022     | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem                                  | No Clinical Signs, Symptoms or Conditions |
| 30067221<br>12-2022-00011 | 12/30/2021 | Malfunction | Apollo Endosurgery, Inc. | 1/12/2022     | LTI          | Orbera Intra gastric Balloon System    | Insufficient Information                           | Flatus; Vomiting                          |

| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                     | Device Problem   | Patient Problem   |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2021-<br>02424 | 12/4/2021  | Malfunction | Spatz Fgia<br>Inc.             | 1/11/2022     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Deflation<br>Problem   | No Clinical Signs,<br>Symptoms or<br>Conditions                       |
| 30067221<br>12-2022-<br>00010 | 12/27/2021 | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 1/10/2022     | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem; Fluid/<br>Blood Leak;<br>Unintended<br>Deflation | Failure of Implant  |
| 30067221<br>12-2022-<br>00008 | 11/9/2021  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 1/7/2022      | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem; Fluid/<br>Blood Leak;<br>Unintended<br>Deflation | Failure of Implant  |
| 30067221<br>12-2022-<br>00006 | 11/23/2021 | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 1/7/2022      | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem; Fluid/<br>Blood Leak;<br>Unintended<br>Deflation | Failure of Implant  |
| 30067221<br>12-2022-<br>00005 | 12/1/2021  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 1/7/2022      | LTI          | Orbera<br>Intragastric<br>Balloon<br>System    | Inflation Problem  | Failure of Implant;<br>Nausea; Vomiting;<br>Obstruction/<br>Occlusion |
| 30067221<br>12-2022-<br>00003 | 12/17/2021 | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 1/7/2022      | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem; Fluid/<br>Blood Leak;<br>Unintended<br>Deflation | Failure of Implant  |
| 30067221<br>12-2022-<br>00004 | 12/17/2021 | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 1/7/2022      | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem; Fluid/<br>Blood Leak;<br>Unintended<br>Deflation | Failure of Implant  |

| Report Number             | Event Date | Event Type  | Manufacturer                         | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem   |
|---------------------------|------------|-------------|--------------------------------------|---------------|--------------|--|---|---|
| 30067221<br>12-2022-00007 | 12/21/2021 | Malfunction | Apollo Endosurgery, Inc.             | 1/7/2022      | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Abdominal Pain; Failure of Implant; Vomiting; Obstruction/Occlusion |
| 30067221<br>12-2022-00009 | 12/23/2021 | Malfunction | Apollo Endosurgery, Inc.             | 1/7/2022      | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant; Vomiting  |
| 30067221<br>12-2022-00001 | 9/20/2021  | Malfunction | Apollo Endosurgery, Inc.             | 1/6/2022      | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant  |
| 30067221<br>12-2022-00002 | 12/14/2021 | Malfunction | Apollo Endosurgery, Inc.             | 1/6/2022      | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/Blood Leak; Unintended Deflation | Failure of Implant  |
| MW5106<br>299             | 11/9/2021  | Injury      | Apollo Endosurgery Costa Rica S.R.L. | 12/27/2021    | LTI          | Orbera                                 | Adverse Event Without Identified Device or Use Problem    | Cardiac Arrest; Hemorrhage/Bleeding; Nausea; Dizziness              |
| 30126289<br>28-2021-02409 | 10/26/2021 | Malfunction | Spatz Fgia Inc.                      | 12/26/2021    | LTI          | Spatz3 Adjustable Balloon System       | Unintended Deflation                                      | Failure of Implant; Pain  |
| 30126389<br>28-2021-00241 | 10/31/2021 | Malfunction | Spatz Fgia Inc.                      | 12/26/2021    | LTI          | Spatz3 Adjustable Balloon System       | Unintended Deflation                                      | Failure of Implant  |
| 30126389<br>28-2021-02406 | 11/8/2021  | Malfunction | Spatz Fgia Inc.                      | 12/26/2021    | LTI          | Spatz3 Adjustable Balloon System       | Deflation Problem   | Insufficient Information  |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                         | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|------------------------------------|--|---|
| 30126389<br>28-2021-02408 | 11/20/2021 | Malfunction | Spatz Fgia Inc.          | 12/26/2021    | LTI          | Spatz3 Adjustable Balloon System   | Deflation Problem; Leak/ Splash                        | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30126389<br>28-2021-02417 | 11/20/2021 | Malfunction | Spatz Fgia Inc.          | 12/26/2021    | LTI          | Spatz3 Adjustable Balloon System   | Insufficient Information                               | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30126389<br>28-2021-02416 | 11/26/2021 | Malfunction | Spatz Fgia Inc.          | 12/26/2021    | LTI          | Spatz3 Adjustable Balloon System   | Output Problem   | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2021-00129 | 12/20/2021 | Malfunction | Apollo Endosurgery, Inc. | 12/21/2021    | LTI          | Orbera® Intragastic Balloon System | Obstruction of Flow                                    | Ischemia  |
| 30126389<br>28-2021-02399 | 10/1/2021  | Malfunction | Spatz Fgia Inc.          | 12/20/2021    | LTI          | Spatz3 Adjustable Balloon System   | Output Problem   | Abdominal Pain  |
| 30126389<br>28-2021-02401 | 10/4/2021  | Malfunction | Spatz Fgia Inc.          | 12/20/2021    | LTI          | Spatz3 Adjustable Balloon System   | Output Problem   | Insufficient Information  |
| 30067221<br>12-2021-00125 | 11/10/2021 | Malfunction | Apollo Endosurgery, Inc. | 12/20/2021    | LTI          | Orbera® Intragastic Balloon System | Inflation Problem                                      | Abdominal Pain; Failure of Implant  |
| 30126389<br>28-2021-02403 | 11/16/2021 | Malfunction | Spatz Fgia Inc.          | 12/20/2021    | LTI          | Spatz3 Adjustable Balloon System   | Adverse Event Without Identified Device or Use Problem | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30126389<br>28-2021-02402 | 11/17/2021 | Malfunction | Spatz Fgia Inc.          | 12/20/2021    | LTI          | Spatz3 Adjustable Balloon System   | Unintended Deflation                                   | Failure of Implant  |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30126389<br>28-2021-02392 | 10/25/2021 | Malfunction | Spatz Fgia Inc.          | 12/14/2021    | LTI          | Spatz3 Adjustable Balloon System       | Biofilm coating in Device; Deflation Problem             | Failure of Implant  |
| 30126389<br>28-2021-02393 | 10/31/2021 | Malfunction | Spatz Fgia Inc.          | 12/14/2021    | LTI          | Spatz3 Adjustable Balloon System       | Leak/ Splash   | Insufficient Information  |
| 30126389<br>28-2021-02394 | 11/3/2021  | Malfunction | Spatz Fgia Inc.          | 12/14/2021    | LTI          | Spatz3 Adjustable Balloon System       | Leak/ Splash   | Insufficient Information  |
| 30067221<br>12-2021-00123 | 11/19/2021 | Malfunction | Apollo Endosurgery, Inc. | 12/14/2021    | LTI          | Orbera Intra-gastric Balloon System    | Inflation Problem  | Failure of Implant  |
| 30067221<br>12-2021-00124 | 12/2/2021  | Malfunction | Apollo Endosurgery, Inc. | 12/14/2021    | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Failure to Deflate | Failure of Implant; Obstruction/ Occlusion                                |
| 30067221<br>12-2021-00120 | 11/12/2021 | Malfunction | Apollo Endosurgery, Inc. | 12/9/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Fluid/ Blood Leak; Unintended Deflation                  | Failure of Implant  |
| 30067221<br>12-2021-00121 | 11/20/2021 | Malfunction | Apollo Endosurgery, Inc. | 12/9/2021     | LTI          | Orbera Intra-gastric Balloon System    | Deflation Problem  | Erosion; Hemorrhage/ Bleeding; Failure of Implant                         |
| 30067221<br>12-2021-00122 | NR         | Malfunction | Apollo Endosurgery, Inc. | 12/9/2021     | LTI          | Orbera® Intra-gastric Balloon System   | Deflation Problem  | Pyrosis/ Heartburn; Nausea; Vomiting; Pancreatitis                        |
| 30126389<br>28-2021-02388 | 10/16/2020 | Malfunction | Spatz Fgia Inc.          | 12/8/2021     | LTI          | Spatz3 Adjustable Balloon System       | Output Problem   | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |



| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                     | Device Problem   | Patient Problem  |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|--|--|--|
| 30126389<br>28-2021-<br>02382 | 7/31/2021  | Malfunction | Spatz Fgia<br>Inc.             | 12/8/2021     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Insufficient<br>Information  | Appropriate Clinical<br>Signs, Symptoms,<br>Conditions Term/<br>Code Not Available |
| 30067221<br>12-2021-<br>00118 | 10/23/2021 | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 12/8/2021     | LTI          | Orbera<br>Intragastric<br>Balloon<br>System    | Fluid/ Blood<br>Leak; Mechanical<br>Problem                            | Failure of Implant   |
| 30067221<br>12-2021-<br>00119 | 11/14/2021 | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 12/8/2021     | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Fluid/ Blood Leak  | Failure of Implant;<br>Vomiting;<br>Obstruction/<br>Occlusion                      |
| 30126389<br>28-2021-<br>02381 | 9/2/2021   | Malfunction | Spatz Fgia<br>Inc.             | 12/4/2021     | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Output Problem   | Insufficient<br>Information  |
| 30126389<br>28-2021-<br>02380 | 9/27/2021  | Malfunction | Spatz Fgia<br>Inc.             | 11/30/2021    | LTI          | Spatz3<br>Adjustable<br>Balloon<br>System      | Output Problem   | Insufficient<br>Information  |
| 30067221<br>12-2021-<br>00112 | 11/4/2021  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 11/18/2021    | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Deflation<br>Problem; Fluid/<br>Blood Leak;<br>Unintended<br>Deflation | Failure of Implant   |
| 30067221<br>12-2021-<br>00113 | 11/5/2021  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 11/18/2021    | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Inflation Problem  | Failure of Implant   |
| 30067221<br>12-2021-<br>00111 | NR         | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 11/17/2021    | LTI          | Orbera<br>Intragastric<br>Balloon<br>System    | Inflation Problem  | Failure of Implant;<br>Nausea; Vomiting  |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem                    |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|------------------------------------|
| 30067221<br>12-2021-00117 | 11/10/2021 | Malfunction | Apollo Endosurgery, Inc. | 11/16/2021    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Failure of Implant                 |
| 30067221<br>12-2021-00115 | 11/10/2021 | Malfunction | Apollo Endosurgery, Inc. | 11/16/2021    | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Fluid/ Blood Leak; Unintended Deflation | Abdominal Pain; Pyrosis/ Heartburn |
| 30067221<br>12-2021-00107 | 10/4/2021  | Malfunction | Apollo Endosurgery, Inc. | 11/2/2021     | LTI          | Orbera Intra gastric Balloon System    | Adverse Event Without Identified Device or Use Problem     | Dehydration; Vomiting; Malaise     |
| 30067221<br>12-2021-00110 | 10/25/2021 | Malfunction | Apollo Endosurgery, Inc. | 11/2/2021     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                    | Failure of Implant                 |
| 30067221<br>12-2021-00108 | 10/14/2021 | Malfunction | Apollo Endosurgery, Inc. | 11/2/2021     | LTI          | Orbera Intra gastric Balloon System    | Insufficient Information                                   | Pancreatitis                       |
| 30067221<br>12-2021-00109 | NR         | Malfunction | Apollo Endosurgery, Inc. | 11/2/2021     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem; Unintended Deflation                    | Failure of Implant                 |
| 30067221<br>12-2021-00095 | 9/29/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/28/2021    | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak  | Failure of Implant                 |
| 30067221<br>12-2021-00093 | 10/2/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/28/2021    | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem  | Failure of Implant                 |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name  | Device Problem                                   | Patient Problem                                      |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|--|--|
| 30067221<br>12-2021-00098 | 10/2/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/28/2021    | LTI          | Orbera365 Intra gastric Balloon System            | Inflation Problem                                | Failure of Implant                                   |
| 30067221<br>12-2021-00104 | 10/5/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/28/2021    | LTI          | Orbera Intra gastric Balloon System               | Fluid/ Blood Leak                                | Failure of Implant                                   |
| 30067221<br>12-2021-00105 | NR         | Malfunction | Apollo Endosurgery, Inc. | 10/28/2021    | LTI          | Orbera365 Intra gastric Balloon System Directions | Fluid/ Blood Leak                                | Abdominal Pain; Failure of Implant; Nausea; Vomiting |
| 30067221<br>12-2021-00094 | 9/18/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/26/2021    | LTI          | Orbera Intra gastric Balloon System               | Migration  | Failure of Implant; Perforation                      |
| 30067221<br>12-2021-00091 | 7/9/2021   | Malfunction | Apollo Endosurgery, Inc. | 10/22/2021    | LTI          | Orbera Intra gastric Balloon System               | Migration  | Abdominal Pain; Vomiting; Obstruction/ Occlusion     |
| 30067221<br>12-2021-00099 | 9/20/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/22/2021    | LTI          | Orbera Intra gastric Balloon System               | Insufficient Information                         | Malaise  |
| 30067221<br>12-2021-00096 | 9/24/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/22/2021    | LTI          | Orbera365 Intra gastric Balloon System            | Deflation Problem; Migration; Failure to Deflate | Airway Obstruction; Failure of Implant               |
| 30067221<br>12-2021-00097 | 9/24/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/22/2021    | LTI          | Orbera365 Intra gastric Balloon System            | Deflation Problem; Migration                     | Airway Obstruction; Failure of Implant               |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem                                      | Patient Problem                              |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|--|
| 30067221<br>12-2021-00092 | 10/1/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/22/2021    | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem; Unintended Deflation             | Failure of Implant                           |
| 30067221<br>12-2021-00103 | 9/2/2021   | Malfunction | Apollo Endosurgery, Inc. | 10/21/2021    | LTI          | Orbera365 Intra-gastric Balloon System | Inflation Problem                                   | Abdominal Pain; Failure of Implant; Vomiting |
| 30067221<br>12-2021-00090 | 9/25/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/21/2021    | LTI          | Orbera Intra-gastric Balloon System    | Inflation Problem                                   | Insufficient Information                     |
| 30067221<br>12-2021-00084 | 9/3/2021   | Malfunction | Apollo Endosurgery, Inc. | 10/20/2021    | LTI          | Orbera365 Intra-gastric Balloon System | Inflation Problem                                   | Failure of Implant                           |
| 30067221<br>12-2021-00086 | 6/7/2021   | Malfunction | Apollo Endosurgery, Inc. | 10/8/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Inflation Problem                                   | Failure of Implant                           |
| 30067221<br>12-2021-00089 | 8/28/2021  | Malfunction | Apollo Endosurgery, Inc. | 10/8/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Unintended Deflation                                | Failure of Implant; Vomiting                 |
| 30067221<br>12-2021-00088 | 8/9/2021   | Malfunction | Apollo Endosurgery, Inc. | 9/30/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Fluid/ Blood Leak                                   | Failure of Implant                           |
| 30067221<br>12-2021-00085 | 9/1/2021   | Malfunction | Apollo Endosurgery, Inc. | 9/29/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Inflation Problem; Migration or Expulsion of Device | Abdominal Pain; Fatigue; Failure of Implant  |
| 30067221<br>12-2021-00082 | 7/10/2021  | Malfunction | Apollo Endosurgery, Inc. | 9/20/2021     | LTI          | Orbera Intra-gastric Balloon System    | Deflation Problem                                   | Abdominal Pain; Failure of Implant; Vomiting |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem                   | Patient Problem                                      |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|----------------------------------|--|
| 30067221<br>12-2021-00083 | 8/18/2021  | Malfunction | Apollo Endosurgery, Inc. | 9/16/2021     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem                | Abdominal Pain; Failure of Implant; Vomiting         |
| 30067221<br>22-2021-00078 | NR         | Malfunction | Apollo Endosurgery, Inc. | 9/15/2021     | LTI          | Orbera365 Intra gastric Balloon System | Inflation Problem; Migration     | Abdominal Pain; Failure of Implant; Nausea; Vomiting |
| 30067221<br>12-2021-00079 | 7/15/2021  | Malfunction | Apollo Endosurgery, Inc. | 9/2/2021      | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem                | Failure of Implant; Pain                             |
| 30067221<br>12-2021-00077 | 7/22/2021  | Malfunction | Apollo Endosurgery, Inc. | 9/2/2021      | LTI          | Orbera365 Intra gastric Balloon System | Migration or Expulsion of Device | Failure of Implant; Vomiting                         |
| 30067221<br>12-2021-00080 | 8/18/2021  | Malfunction | Apollo Endosurgery, Inc. | 9/2/2021      | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem                | Failure of Implant; Vomiting; Obstruction/ Occlusion |
| 30067221<br>12-2021-00076 | 7/20/2021  | Malfunction | Apollo Endosurgery, Inc. | 8/31/2021     | LTI          | Orbera365 Intra gastric Balloon System | Inflation Problem                | Failure of Implant; Pain; Vomiting                   |
| 30067221<br>12-2021-00075 | 7/26/2021  | Malfunction | Apollo Endosurgery, Inc. | 8/31/2021     | LTI          | Orbera365 Intra gastric Balloon System | Inflation Problem                | Failure of Implant; Pain; Vomiting                   |
| 30067221<br>12-2021-00081 | 8/10/2021  | Malfunction | Apollo Endosurgery, Inc. | 8/26/2021     | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak                | Failure of Implant                                   |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem  | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|---|
| 30067221<br>12-2021-00074 | 2/2/2021   | Malfunction | Apollo Endosurgery, Inc. | 8/25/2021     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem; Adverse Event Without Identified Device or Use Problem | Failure of Implant  |
| 30067221<br>12-2021-00073 | 7/19/2021  | Malfunction | Apollo Endosurgery, Inc. | 8/23/2021     | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak; Unintended Deflation                                   | Failure of Implant  |
| 30067221<br>12-2021-00072 | 6/19/2021  | Malfunction | Apollo Endosurgery, Inc. | 8/23/2021     | LTI          | Orbera Intra gastric Balloon System    | Inflation Problem   | Abdominal Pain; Failure of Implant; Vomiting              |
| 30067221<br>12-2021-00069 | 4/25/2021  | Malfunction | Apollo Endosurgery, Inc. | 8/18/2021     | LTI          | Orbera Intra gastric Balloon System    | Insufficient Information  | Abdominal Pain; Failure of Implant                        |
| 30067221<br>12-2021-00067 | 7/8/2021   | Malfunction | Apollo Endosurgery, Inc. | 8/10/2021     | LTI          | Orbera365 Intra gastric Balloon System | Adverse Event Without Identified Device or Use Problem                    | Perforation of Esophagus                                  |
| 30067221<br>12-2021-00070 | 7/16/2021  | Malfunction | Apollo Endosurgery, Inc. | 8/10/2021     | LTI          | Orbera365 Intra gastric Balloon System | Deflation Problem   | Failure of Implant  |
| 30067221<br>12-2021-00064 | 4/10/2021  | Malfunction | Apollo Endosurgery, Inc. | 8/4/2021      | LTI          | Orbera365 Intra gastric Balloon System | Patient-Device Incompatibility  | Failure of Implant; Vomiting; Obstruction/ Occlusion      |
| 30067221<br>12-2021-00066 | NR         | Injury      | Apollo Endosurgery, Inc. | 7/29/2021     | LTI          | Orbera Intra gastric Balloon System    | Insufficient Information  | Failure of Implant; Perforation; Perforation of Esophagus |

| Report Number             | Event Date | Event Type  | Manufacturer                       | Date Received | Product Code | Brand Name                             | Device Problem                 | Patient Problem                                  |
|---------------------------|------------|-------------|------------------------------------|---------------|--------------|--|--------------------------------|--|
| 30067221<br>12-2021-00065 | NR         | Death       | Apollo Endosurgery, Inc.           | 7/29/2021     | LTI          | Orbera® Intra-gastric Balloon System   | Insufficient Information       | Insufficient Information                         |
| 30067221<br>12-2021-00068 | 6/24/2021  | Malfunction | Apollo Endosurgery, Inc.           | 7/22/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Fluid/ Blood Leak              | Failure of Implant; Vomiting                     |
| 30067221<br>12-2021-00038 | 5/24/2021  | Malfunction | Apollo Endosurgery, Inc.           | 7/13/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Mechanical Problem             | Failure of Implant; Nausea; Vomiting             |
| 30067221<br>12-2021-00055 | 5/3/2021   | Malfunction | Apollo Endosurgery, Inc.           | 7/6/2021      | LTI          | Orbera365 Intra-gastric Balloon System | Inflation Problem              | Nausea; Vomiting                                 |
| 30067221<br>12-2021-00056 | 6/2/2021   | Malfunction | Apollo Endosurgery, Inc.           | 6/29/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Fluid/ Blood Leak              | Failure of Implant                               |
| 30067221<br>12-2021-00053 | 5/21/2021  | Malfunction | Apollo Endosurgery, Inc.           | 6/22/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Fluid/ Blood Leak              | Failure of Implant                               |
| 30067221<br>12-2021-00050 | 4/15/2021  | Malfunction | Apollo Endosurgery, Inc.           | 6/22/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Fluid/ Blood Leak              | Failure of Implant                               |
| 30067221<br>12-2021-00062 | NR         | Injury      | Apollo Endosurgery Costa Rica, Srl | 6/22/2021     | LTI          | Intra-gastric Balloon System           | Patient-Device Incompatibility | Nausea; Pain; Abdominal Distention; Pancreatitis |
| 30067221<br>12-2021-00046 | NR         | Malfunction | Apollo Endosurgery, Inc.           | 6/14/2021     | LTI          | Orbera Intra-gastric Balloon System    | Patient-Device Incompatibility | Failure of Implant                               |

| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                     | Device Problem  | Patient Problem  |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|--|---|--|
| 30067221<br>12-2021-<br>00049 | 1/29/2021  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/10/2021     | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Patient-Device<br>Incompatibility                               | Failure of Implant;<br>Pain; Vomiting                  |
| 30067221<br>12-2021-<br>00047 | 3/30/2021  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/8/2021      | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Fluid/ Blood Leak   | Failure of Implant                                     |
| 30067221<br>12-2021-<br>00048 | 5/14/2021  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/8/2021      | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Patient-Device<br>Incompatibility                               | Failure of Implant;<br>Pain; Vomiting                  |
| 30067221<br>12-2021-<br>00051 | 4/14/2021  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/2/2021      | LTI          | Orbera<br>Intragastric<br>Balloon<br>System    | Patient-Device<br>Incompatibility                               | Vomiting   |
| 30067221<br>12-2021-<br>00042 | 5/10/2021  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/2/2021      | LTI          | Orbera365<br>Intragastric<br>Balloon<br>System | Fluid/ Blood Leak   | Abdominal Pain;<br>Failure of Implant;<br>Implant Pain |
| 30067221<br>12-2021-<br>00039 | 1/25/2021  | Injury      | Apollo<br>Endosurgery,<br>Inc. | 6/1/2021      | LTI          | Orbera<br>Intragastric<br>Balloon<br>System    | Adverse Event<br>Without<br>Identified Device<br>or Use Problem | Abdominal Pain;<br>Vomiting                            |
| 30067221<br>12-2021-<br>00032 | 2/5/2021   | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 5/24/2021     | LTI          | Orbera<br>Intragastric<br>Balloon<br>System    | Fluid/ Blood Leak   | Failure of Implant                                     |
| 30067221<br>12-2021-<br>00037 | 4/16/2021  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 5/19/2021     | OCW          | Orbera<br>Intragastric<br>Balloon<br>System    | Patient-Device<br>Incompatibility                               | Eructate; Flatus;<br>Abdominal<br>Distention           |
| 30067221<br>12-2021-<br>00041 | 4/1/2021   | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 5/12/2021     | LTI          | Orbera 365<br>Intragastric<br>Balloon          | Fluid/ Blood Leak   | Failure of Implant                                     |



| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2021-00040 | 3/4/2021   | Malfunction | Apollo Endosurgery, Inc. | 5/11/2021     | LTI          | Orbera Intra-gastric Balloon System    | Fluid/ Blood Leak; Inflation Problem                   | Failure of Implant  |
| 30067221<br>12-2021-00031 | 4/7/2021   | Injury      | Apollo Endosurgery, Inc. | 5/7/2021      | OCZ          | BIB Removal Tool Wire Grasper          | Mechanical Problem; Device Damaged by Another Device   | No Clinical Signs, Symptoms or Conditions                                     |
| 30067221<br>12-2021-00029 | 3/15/2021  | Malfunction | Apollo Endosurgery, Inc. | 4/21/2021     | LTI          | Orbera365 Intra-gastric Balloon System | Deflation Problem                                      | Failure of Implant  |
| 30067221<br>12-2021-00024 | 2/25/2021  | Injury      | Apollo Endosurgery, Inc. | 4/16/2021     | LTI          | Orbera 365 Intra-gastric Balloon       | Adverse Event Without Identified Device or Use Problem | Eructate; Pyrosis/ Heartburn; Nausea; Vertigo; Vomiting; Abdominal Distention |
| 30067221<br>12-2021-00021 | 3/19/2021  | Malfunction | Apollo Endosurgery, Inc. | 4/16/2021     | LTI          | Orbera Intra-gastric Balloon System    | Inflation Problem; Air/ Gas in Device                  | Diarrhea; Nausea; Vomiting; Abdominal Distention                              |
| 30067221<br>12-2021-00026 | 3/11/2021  | Injury      | Apollo Endosurgery, Inc. | 4/8/2021      | LTI          | Orbera365 Intra-gastric Balloon System | Difficult to Remove                                    | Obstruction/ Occlusion  |
| 30067221<br>12-2021-00022 | 3/4/2021   | Malfunction | Apollo Endosurgery, Inc. | 4/1/2021      | LTI          | Orbera Intra-gastric Balloon System    | Mechanical Problem                                     | Failure of Implant  |
| MW5100<br>421             | 2/23/2021  | Malfunction | Apollo Endosurgery Inc   | 3/30/2021     | LTI          | Orbera Balloon                         | Adverse Event Without Identified Device or Use Problem | Nausea; Pain; Vomiting; Pressure Sores; Sleep Dysfunction; Weight Changes     |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|--|---|
| 30067221<br>12-2021-00018 | 11/17/2020 | Malfunction | Apollo Endosurgery       | 3/22/2021     | LTI          | Orbera Intra gastric Balloon System    | Patient-Device Incompatibility                         | Failure of Implant  |
| 30067221<br>12-2021-00019 | 2/18/2021  | Malfunction | Apollo Endosurgery, Inc. | 3/19/2021     | LTI          | Orbera 365 Intra gastric Balloon       | Patient-Device Incompatibility                         | Failure of Implant  |
| 30067221<br>12-2021-00014 | NR         | Death       | Apollo Endosurgery, Inc. | 3/19/2021     | LTI          | Orbera Intra gastric Balloon System    | Insufficient Information                               | Perforation; Insufficient Information                       |
| 30067221<br>12-2021-00017 | 12/22/2020 | Malfunction | Apollo Endosurgery, Inc. | 3/11/2021     | LTI          | Orbera365 Intra gastric Balloon System | Patient-Device Incompatibility                         | Failure of Implant; Vomiting; Ulcer; Obstruction/ Occlusion |
| 30067221<br>12-2021-00013 | 1/25/2021  | Malfunction | Apollo Endosurgery, Inc. | 3/10/2021     | LTI          | Orbera Intra gastric Balloon System    | Fluid/ Blood Leak                                      | Failure of Implant  |
| 30067221<br>12-2021-00016 | 12/10/2020 | Malfunction | Apollo Endosurgery       | 3/10/2021     | LTI          | Orbera365 Intra gastric Balloon System | Patient-Device Incompatibility                         | Failure of Implant; Vomiting; Obstruction/ Occlusion        |
| 30067221<br>12-2021-00012 | 1/23/2021  | Malfunction | Apollo Endosurgery, Inc. | 3/5/2021      | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak                                      | Failure of Implant  |
| 30067221<br>12-2021-00010 | 10/1/2019  | Injury      | Apollo Endosurgery       | 3/1/2021      | LTI          | Orbera Intra gastric Balloon           | Adverse Event Without Identified Device or Use Problem | Dehydration; Vomiting; Obstruction/ Occlusion               |
| 30067221<br>12-2021-00011 | 1/23/2021  | Injury      | Apollo Endosurgery, Inc. | 2/23/2021     | OCW          | Orbera Intra gastric Balloon           | Fluid/ Blood Leak                                      | Failure of Implant; Stomach Ulceration                      |

| Report Number             | Event Date | Event Type  | Manufacturer       | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem  |
|---------------------------|------------|-------------|--------------------|---------------|--------------|--|--|--|
| 30067221<br>12-2021-00009 | 1/16/2021  | Malfunction | Apollo Endosurgery | 2/16/2021     | LTI          | Orbera365 Intra gastric Balloon System | Human-Device Interface Problem                         | Failure of Implant   |
| 30067221<br>12-2021-00008 | 1/19/2021  | Malfunction | Apollo Endosurgery | 2/16/2021     | LTI          | Orbera Intra gastric Balloon System    | Fluid/ Blood Leak                                      | Failure of Implant   |
| 30067221<br>12-2021-00006 | 1/8/2021   | Injury      | Apollo Endosurgery | 2/3/2021      | LTI          | Orbera Intra gastric Balloon System    | Insufficient Information                               | Perforation of Esophagus                                     |
| 30067221<br>12-2021-00005 | 1/15/2021  | Injury      | Apollo Endosurgery | 2/3/2021      | LTI          | Orbera Intra gastric Balloon System    | Adverse Event Without Identified Device or Use Problem | Dehydration; Nausea; Vomiting; Malaise                       |
| 30067221<br>12-2021-00007 | NR         | Injury      | Apollo Endosurgery | 1/28/2021     | LTI          | Orbera Intra gastric Balloon System    | Patient-Device Incompatibility                         | Abdominal Pain; Nausea; Vomiting; Constipation; Pancreatitis |
| 30067221<br>12-2021-00004 | 1/4/2021   | Malfunction | Apollo Endosurgery | 1/27/2021     | LTI          | Orbera Intra gastric Balloon System    | Fluid/ Blood Leak                                      | Obstruction/ Occlusion                                       |
| 30067221<br>12-2021-00003 | 12/24/2020 | Malfunction | Apollo Endosurgery | 1/20/2021     | LTI          | Orbera Intra gastric Balloon System    | Patient-Device Incompatibility                         | Vomiting   |
| 30067221<br>12-2021-00001 | 9/28/2020  | Malfunction | Apollo Endosurgery | 1/11/2021     | LTI          | Orbera Intra gastric Balloon System    | Patient-Device Incompatibility                         | Failure of Implant; Nausea; Pain; Vomiting                   |
| 30067221<br>12-2021-00002 | 11/15/2020 | Malfunction | Apollo Endosurgery | 1/11/2021     | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak                                      | Failure of Implant   |

| Report Number             | Event Date | Event Type  | Manufacturer       | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem  |
|---------------------------|------------|-------------|--------------------|---------------|--------------|--|--|--|
| 30067221<br>12-2020-00134 | 11/25/2020 | Injury      | Apollo Endosurgery | 1/8/2021      | LTI          | Orbera365 Intra-gastric Balloon System | Patient-Device Incompatibility                         | Perforation; Sepsis; Vomiting  |
| 30067221<br>12-2020-00133 | NR         | Malfunction | Apollo Endosurgery | 1/5/2021      | LTI          | Orbera365 Intra-gastric Balloon System | Fluid/ Blood Leak; Migration or Expulsion of Device    | Abdominal Pain   |
| 30067221<br>12-2020-00132 | 9/12/2020  | Injury      | Apollo Endosurgery | 12/28/2020    | LTI          | Orbera365 Intra-gastric Balloon System | Adverse Event Without Identified Device or Use Problem | Abdominal Pain; Perforation; Vomiting; Abdominal Cramps; Hematemesis |
| 30067221<br>12-2020-00130 | 10/31/2020 | Malfunction | Apollo Endosurgery | 12/21/2020    | LTI          | Orbera Intra-gastric Balloon System    | Fluid/ Blood Leak                                      | Failure of Implant   |
| 30067221<br>12-2020-00128 | 10/24/2020 | Injury      | Apollo Endosurgery | 12/17/2020    | LTI          | Orbera Intra-gastric Balloon System    | Adverse Event Without Identified Device or Use Problem | Pancreatitis   |
| 30067221<br>12-2020-00123 | 11/14/2020 | Death       | Apollo Endosurgery | 12/16/2020    | LTI          | Orbera Intra-gastric Balloon System    | Adverse Event Without Identified Device or Use Problem | Nausea; Vomiting; Bowel Perforation                                  |
| 30067221<br>12-2020-00124 | 10/29/2020 | Malfunction | Apollo Endosurgery | 12/15/2020    | LTI          | Orbera Intra-gastric Balloon System    | Patient-Device Incompatibility                         | Pain   |
| 30067221<br>12-2020-00121 | 10/10/2020 | Malfunction | Apollo Endosurgery | 12/8/2020     | LTI          | Orbera Intra-gastric Balloon System    | Patient-Device Incompatibility                         | Abdominal Pain; Vomiting   |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                                     | Device Problem                                      | Patient Problem                                      |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|--|---|--|
| 30067221<br>12-2020-00120 | 10/25/2020 | Malfunction | Apollo Endosurgery       | 12/1/2020     | LTI          | BIB System Intra-gastric Balloon               | Fluid/ Blood Leak; Migration or Expulsion of Device | Abdominal Pain; Failure of Implant; Nausea; Vomiting |
| 30067221<br>12-2020-00115 | 10/9/2020  | Injury      | Apollo Endosurgery       | 11/23/2020    | LTI          | Orbera® Intra-gastric Balloon System           | Patient-Device Incompatibility                      | Vomiting   |
| 30067221<br>12-2020-00118 | 10/28/2020 | Malfunction | Apollo Endosurgery       | 11/23/2020    | LTI          | Orbera Intra-gastric Balloon System            | Patient-Device Incompatibility                      | Abdominal Pain                                       |
| 30067221<br>12-2020-00114 | 7/29/2020  | Malfunction | Apollo Endosurgery       | 11/20/2020    | LTI          | Orbera Intra-gastric Balloon System            | Difficult to Remove; Short Fill                     | Failure of Implant                                   |
| MW5097<br>989             | 9/28/2020  | Injury      | Apollo Endosurgery, Inc. | 11/19/2020    | LTI          | Orbera   | Failure to Deflate                                  | Abdominal Pain; Nausea; Vomiting                     |
| 30067221<br>12-2020-00110 | 9/25/2020  | Malfunction | Apollo Endosurgery       | 11/17/2020    | LTI          | Orbera Intra-gastric Balloon System            | Patient-Device Incompatibility                      | Failure of Implant                                   |
| 30067221<br>12-2020-00111 | 9/30/2020  | Malfunction | Apollo Endosurgery       | 11/17/2020    | LTI          | Orbera Intra-gastric Balloon System Directions | Patient-Device Incompatibility                      | Failure of Implant                                   |
| 30067221<br>12-2020-00107 | 7/19/2020  | Death       | Apollo Endosurgery       | 11/10/2020    | LTI          | BIB System Intra-gastric Balloon               | Insufficient Information                            | Insufficient Information                             |
| 30067221<br>12-2020-00112 | 10/1/2019  | Malfunction | Apollo Endosurgery       | 11/9/2020     | LTI          | Orbera365 Intra-gastric Balloon System         | Fluid/ Blood Leak                                   | No Clinical Signs, Symptoms or Conditions            |

| Report Number             | Event Date | Event Type  | Manufacturer       | Date Received | Product Code | Brand Name                             | Device Problem                                      | Patient Problem                                  |
|---------------------------|------------|-------------|--------------------|---------------|--------------|--|---|--|
| 30067221<br>12-2020-00109 | 9/29/2020  | Malfunction | Apollo Endosurgery | 11/9/2020     | LTI          | Orbera Intra gastric Balloon System    | Patient-Device Incompatibility                      | Abdominal Pain; Vomiting                         |
| 30067221<br>12-2020-00105 | 9/28/2020  | Malfunction | Apollo Endosurgery | 11/4/2020     | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak                                   | No Clinical Signs, Symptoms or Conditions        |
| 30067221<br>12-2020-00104 | 9/12/2020  | Malfunction | Apollo Endosurgery | 11/4/2020     | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak                                   | No Clinical Signs, Symptoms or Conditions        |
| 30067221<br>12-2020-00106 | NR         | Injury      | Apollo Endosurgery | 10/30/2020    | LTI          | Orbera Intra gastric Balloon System    | Insufficient Information                            | Nausea; Vomiting                                 |
| 30067221<br>12-2020-00101 | 5/6/2020   | Malfunction | Apollo Endosurgery | 10/26/2020    | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak                                   | No Clinical Signs, Symptoms or Conditions        |
| 30067221<br>12-2020-00102 | 4/29/2020  | Malfunction | Apollo Endosurgery | 10/26/2020    | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak; Migration or Expulsion of Device | No Clinical Signs, Symptoms or Conditions        |
| 30067221<br>12-2020-00100 | 5/13/2020  | Malfunction | Apollo Endosurgery | 10/26/2020    | LTI          | Orbera365 Intra gastric Balloon System | Fluid/ Blood Leak; Migration or Expulsion of Device | No Clinical Signs, Symptoms or Conditions        |
| 30067221<br>12-2020-00099 | 9/15/2020  | Malfunction | Apollo Endosurgery | 10/26/2020    | LTI          | BIB System Intra gastric Balloon       | Patient-Device Incompatibility                      | No Clinical Signs, Symptoms or Conditions        |
| 30067221<br>12-2020-00098 | 8/28/2020  | Injury      | Apollo Endosurgery | 10/23/2020    | LTI          | Orbera Intra gastric Balloon System    | Fluid/ Blood Leak; Insufficient Information         | Obstruction/ Occlusion; Insufficient Information |

| Report Number             | Event Date | Event Type  | Manufacturer       | Date Received | Product Code | Brand Name                              | Device Problem   | Patient Problem                                |
|---------------------------|------------|-------------|--------------------|---------------|--------------|---|--|--|
| 30067221<br>12-2020-00108 | 9/24/2020  | Malfunction | Apollo Endosurgery | 10/22/2020    | LTI          | Orbera365 Intra-gastric Balloon System  | Migration or Expulsion of Device                       | No Clinical Signs, Symptoms or Conditions      |
| 30067221<br>12-2020-00097 | 7/24/2020  | Malfunction | Apollo Endosurgery | 10/21/2020    | LTI          | Orbera Intra-gastric Balloon System     | Patient-Device Incompatibility                         | Abdominal Pain; Vomiting                       |
| 30067221<br>12-2020-00092 | 5/6/2020   | Malfunction | Apollo Endosurgery | 10/13/2020    | LTI          | Orbera Intra-gastric Balloon System     | Patient-Device Incompatibility                         | Abdominal Pain; Vomiting; Abdominal Distention |
| 30067221<br>12-2020-00103 | 8/29/2020  | Malfunction | Apollo Endosurgery | 10/13/2020    | LTI          | Orbera365 Intra-gastric Balloon System  | Fluid/ Blood Leak                                      | No Clinical Signs, Symptoms or Conditions      |
| 30067221<br>12-2020-00090 | 6/22/2020  | Injury      | Apollo Endosurgery | 10/12/2020    | LTI          | BIB System Intra-gastric Balloon        | Adverse Event Without Identified Device or Use Problem | Abdominal Pain; Perforation                    |
| 30067221<br>12-2020-00089 | 9/8/2020   | Malfunction | Apollo Endosurgery | 10/6/2020     | LTI          | Orbera 365 Intra-gastric Balloon System | Patient-Device Incompatibility                         | No Clinical Signs, Symptoms or Conditions      |
| 30067221<br>12-2020-00088 | NR         | Injury      | Apollo Endosurgery | 10/6/2020     | LTI          | Orbera Intra-gastric Balloon System     | Adverse Event Without Identified Device or Use Problem | Vomiting                                       |
| 30067221<br>12-2020-00091 | NR         | Injury      | Apollo Endosurgery | 10/6/2020     | LTI          | Orbera Intra-gastric Balloon System     | Patient-Device Incompatibility                         | Pancreatitis                                   |
| 30067221<br>12-2020-00087 | 8/26/2020  | Malfunction | Apollo Endosurgery | 10/5/2020     | LTI          | Orbera 365 Intra-gastric Balloon System | Fluid/ Blood Leak                                      | No Clinical Signs, Symptoms or Conditions      |

| Report Number             | Event Date | Event Type  | Manufacturer       | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------|---------------|--------------|-------------------------------------|--|---|
| 30067221<br>12-2020-00086 | 8/29/2020  | Malfunction | Apollo Endosurgery | 9/30/2020     | LTI          | BIB System Intra gastric Balloon    | Patient-Device Incompatibility                         | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available                     |
| 30057221<br>12-2020-00084 | 6/26/2020  | Malfunction | Apollo Endosurgery | 9/18/2020     | LTI          | Orbera Intra gastric Balloon System | Fluid/ Blood Leak                                      | Failure of Implant; Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2020-00082 | 6/16/2020  | Malfunction | Apollo Endosurgery | 9/17/2020     | LTI          | Orbera Intra gastric Balloon System | Fluid/ Blood Leak                                      | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available                     |
| 30067221<br>12-2020-00083 | 7/7/2020   | Injury      | Apollo Endosurgery | 9/17/2020     | LTI          | Orbera Intra gastric Balloon System | Adverse Event Without Identified Device or Use Problem | Vomiting; Increased Sensitivity   |
| 30067221<br>12-2020-00077 | 1/12/2020  | Malfunction | Apollo Endosurgery | 9/16/2020     | LTI          | Orbera Intra gastric Balloon System | Patient-Device Incompatibility                         | Pain; Vomiting  |
| 30067221<br>12-2020-00080 | 3/7/2020   | Malfunction | Apollo Endosurgery | 9/16/2020     | LTI          | Orbera Intra gastric Balloon System | Patient-Device Incompatibility                         | Pain; Vomiting  |
| 30067221<br>12-2020-00076 | 3/9/2020   | Malfunction | Apollo Endosurgery | 9/16/2020     | LTI          | Orbera Intra gastric Balloon System | Off-Label Use; Patient-Device Incompatibility          | Pain; Vomiting  |
| 30067221<br>12-2020-00081 | 3/27/2020  | Malfunction | Apollo Endosurgery | 9/16/2020     | LTI          | Orbera Intra gastric Balloon System | Patient-Device Incompatibility                         | Pain; Vomiting  |



| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|-------------------------------------|--|---|
| 30067221<br>12-2020-00078 | 6/17/2020  | Malfunction | Apollo Endosurgery       | 9/16/2020     | LTI          | Orbera Intra-gastric Balloon System | Patient-Device Incompatibility                         | Pain  |
| 30067221<br>12-2020-00079 | 7/25/2020  | Malfunction | Apollo Endosurgery       | 9/16/2020     | LTI          | Orbera Intra-gastric Balloon System | Patient-Device Incompatibility                         | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available                     |
| 30067221<br>12-2020-00074 | 3/26/2020  | Malfunction | Apollo Endosurgery       | 9/8/2020      | LTI          | Orbera Intra-gastric Balloon System | Fluid/ Blood Leak                                      | Failure of Implant; Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2020-00075 | 4/10/2020  | Injury      | Apollo Endosurgery       | 9/8/2020      | LTI          | Orbera Intra-gastric Balloon System | Adverse Event Without Identified Device or Use Problem | Pain; Abdominal Distention  |
| 30067221<br>12-2020-00073 | 5/23/2020  | Malfunction | Apollo Endosurgery       | 9/8/2020      | LTI          | Orbera Intra-gastric Balloon System | Fluid/ Blood Leak                                      | Failure of Implant; No Known Impact Or Consequence To Patient                                 |
| 30067221<br>12-2020-00072 | NR         | Injury      | Apollo Endosurgery       | 9/8/2020      | LTI          | Orbera Intra-gastric Balloon System | Adverse Event Without Identified Device or Use Problem | Nausea; Vomiting  |
| 30067221<br>12-2020-00071 | 2/3/2020   | Death       | Apollo Endosurgery       | 9/2/2020      | LTI          | Orbera Intra-gastric Balloon System | Adverse Event Without Identified Device or Use Problem | Sudden Cardiac Death  |
| MW5096<br>356             | 8/28/2020  | Injury      | Apollo Endosurgery Inc   | 8/31/2020     | LTI          | Orbera                              | Appropriate Term/ Code Not Available                   | Obstruction/ Occlusion; Foreign Body In Patient   |
| 30067221<br>12-2020-00065 | 7/9/2020   | Malfunction | Apollo Endosurgery, Inc. | 8/18/2020     | LTI          | Orbera365 Intra-gastric Balloon     | Patient-Device Incompatibility                         | Patient Problem/ Medical Problem  |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                              | Device Problem                                      | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|---|---|
| 30067221<br>12-2020-00070 | NR         | Malfunction | Apollo Endosurgery, Inc. | 8/6/2020      | LTI          | Orbera365 Intra gastric Balloon System  | Patient-Device Incompatibility                      | No Code Available                         |
| 30067221<br>12-2020-00069 | 1/31/2020  | Malfunction | Apollo Endosurgery, Inc. | 8/4/2020      | LTI          | Orbera Intra gastric Balloon System     | Patient-Device Incompatibility                      | Nausea; Vomiting; Abdominal Distention    |
| 30067221<br>12-2020-00064 | 6/18/2020  | Malfunction | Apollo Endosurgery, Inc. | 7/20/2020     | LTI          | Orbera 365 Intra gastric Balloon System | Patient-Device Incompatibility                      | Patient Problem/ Medical Problem          |
| 30067221<br>12-2020-00062 | 6/24/2020  | Malfunction | Apollo Endosurgery, Inc. | 7/20/2020     | LTI          | BIB System Intra gastric Balloon        | Patient-Device Incompatibility                      | Patient Problem/ Medical Problem          |
| 30067221<br>12-2020-00063 | 6/24/2020  | Malfunction | Apollo Endosurgery, Inc. | 7/20/2020     | LTI          | Orbera Intra gastric Balloon System     | Patient-Device Incompatibility                      | Pain; Vomiting                            |
| 30067221<br>12-2020-00061 | 5/2/2020   | Malfunction | Apollo Endosurgery, Inc. | 7/20/2020     | LTI          | Orbera 365 Intra gastric Balloon System | Fluid/ Blood Leak                                   | Nausea; Vomiting                          |
| 30067221<br>12-2020-00059 | 5/11/2020  | Malfunction | Apollo Endosurgery, Inc. | 6/10/2020     | LTI          | Orbera Intra gastric Balloon            | Fluid/ Blood Leak                                   | No Known Impact Or Consequence To Patient |
| 30067221<br>12-2020-00058 | 5/6/2020   | Malfunction | Apollo Endosurgery, Inc. | 6/8/2020      | LTI          | Orbera Intra gastric Balloon System     | Fluid/ Blood Leak                                   | Patient Problem/ Medical Problem          |
| 30067221<br>12-2020-00057 | 5/2/2020   | Malfunction | Apollo Endosurgery, Inc. | 6/4/2020      | LTI          | Orbera Intra gastric Balloon System     | Fluid/ Blood Leak; Migration or Expulsion of Device | Nausea; Vomiting                          |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                              | Device Problem   | Patient Problem                           |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|--|---|
| 30067221<br>12-2020-00056 | 4/30/2020  | Injury      | Apollo Endosurgery, Inc. | 6/2/2020      | LTI          | BIB System Intra-gastric Balloon        | Adverse Event Without Identified Device or Use Problem | Nausea; Pain; Obstruction/ Occlusion      |
| 30067221<br>12-2020-00055 | 2/4/2020   | Injury      | Apollo Endosurgery, Inc. | 5/30/2020     | LTI          | Orbera Intra-gastric Balloon System     | Adverse Event Without Identified Device or Use Problem | Internal Organ Perforation                |
| 30067221<br>12-2020-00053 | 4/9/2020   | Malfunction | Apollo Endosurgery, Inc. | 5/22/2020     | LTI          | Orbera 365 Intra-gastric Balloon System | Fluid/ Blood Leak                                      | Failure of Implant                        |
| 30067221<br>12-2020-00049 | 1/29/2020  | Injury      | Apollo Endosurgery, Inc. | 4/13/2020     | LTI          | Orbera 365 Intra-gastric Balloon System | Adverse Event Without Identified Device or Use Problem | Hypersensitivity/ Allergic reaction; Pain |
| 30067221<br>12-2020-00050 | 3/6/2020   | Malfunction | Apollo Endosurgery, Inc. | 4/13/2020     | LTI          | Orbera 365 Intra-gastric Balloon System | Fluid/ Blood Leak                                      | Failure of Implant                        |
| 30067221<br>12-2020-00047 | 3/11/2020  | Injury      | Apollo Endosurgery, Inc. | 4/9/2020      | LTI          | Orbera Intra-gastric Balloon System     | Adverse Event Without Identified Device or Use Problem | Syncope; Dyspnea; Pain; No Code Available |
| 30067221<br>12-2020-00048 | 3/8/2020   | Injury      | Apollo Endosurgery, Inc. | 4/7/2020      | LTI          | Orbera Intra-gastric Balloon System     | Adverse Event Without Identified Device or Use Problem | Pain; Vomiting                            |
| 30067221<br>12-2020-00045 | 1/12/2020  | Injury      | Apollo Endosurgery, Inc. | 4/7/2020      | LTI          | Orbera Intra-gastric Balloon System     | Patient-Device Incompatibility                         | Pain; Vomiting                            |
| 30067221<br>12-2020-00043 | 1/28/2020  | Malfunction | Apollo Endosurgery, Inc. | 4/7/2020      | LTI          | Orbera 365 Intra-gastric Balloon System | Fluid/ Blood Leak                                      | Nausea; Vomiting                          |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                              | Device Problem   | Patient Problem                                       |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|--|---|
| 30067221<br>12-2020-00039 | 2/19/2020  | Malfunction | Apollo Endosurgery, Inc. | 3/17/2020     | LTI          | Orbera 365 Intra gastric Balloon System | Fluid/ Blood Leak; Migration or Expulsion of Device    | Weight Changes; Patient Problem/ Medical Problem      |
| 30067221<br>12-2020-00042 | 2/19/2020  | Injury      | Apollo Endosurgery, Inc. | 3/17/2020     | LTI          | Orbera Intra gastric Balloon System     | Adverse Event Without Identified Device or Use Problem | Nausea; Pain; Vomiting; Abdominal Cramps              |
| 30067221<br>12-2020-00044 | 2/20/2020  | Malfunction | Apollo Endosurgery, Inc. | 3/17/2020     | LTI          | Orbera 365 Intra gastric Balloon System | Fluid/ Blood Leak; Migration or Expulsion of Device    | Malaise   |
| 30067221<br>12-2020-00038 | NR         | Injury      | Apollo Endosurgery, Inc. | 3/9/2020      | LTI          | Orbera® Intra gastric Balloon System    | Adverse Event Without Identified Device or Use Problem | Dehydration; Nausea; Pain; Vomiting; Abdominal Cramps |
| 30067221<br>12-2020-00035 | 12/29/2019 | Malfunction | Apollo Endosurgery, Inc. | 3/9/2020      | LTI          | Orbera® Intra gastric Balloon System    | Fluid/ Blood Leak                                      | Failure of Implant                                    |
| 30067221<br>12-2020-00037 | 1/5/2020   | Malfunction | Apollo Endosurgery, Inc. | 3/9/2020      | LTI          | Orbera® Intra gastric Balloon System    | Patient-Device Incompatibility                         | Nausea; Pain; Vomiting                                |
| 30067221<br>12-2020-00040 | 2/3/2020   | Malfunction | Apollo Endosurgery, Inc. | 3/9/2020      | LTI          | Orbera® Intra gastric Balloon System    | Patient-Device Incompatibility                         | Pain; Vomiting; Abdominal Distention                  |
| 30067221<br>12-2020-00036 | 2/13/2020  | Malfunction | Apollo Endosurgery, Inc. | 3/9/2020      | LTI          | Bib; System Intra gastric Balloon       | Fluid/ Blood Leak                                      | Failure of Implant                                    |
| 30067221<br>12-2020-00034 | NR         | Malfunction | Apollo Endosurgery, Inc. | 3/9/2020      | LTI          | Orbera365® Intra gastric Balloon        | Fluid/ Blood Leak                                      | Failure of Implant                                    |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                              | Device Problem   | Patient Problem                                    |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|--|--|
| 30067221<br>12-2020-00033 | NR         | Malfunction | Apollo Endosurgery, Inc. | 3/9/2020      | LTI          | Orbera365® Intra gastric Balloon        | Fluid/ Blood Leak  | Failure of Implant                                 |
| 30067221<br>12-2020-00030 | 12/2/2019  | Malfunction | Apollo Endosurgery, Inc. | 3/3/2020      | LTI          | Orbera 365 Intra gastric Balloon System | Fluid/ Blood Leak  | Failure of Implant                                 |
| 30067221<br>12-2020-00031 | 12/15/2019 | Malfunction | Apollo Endosurgery, Inc. | 3/3/2020      | LTI          | Orbera Intra gastric Balloon System     | Fluid/ Blood Leak  | Failure of Implant                                 |
| 30067221<br>12-2020-00032 | 1/27/2020  | Injury      | Apollo Endosurgery, Inc. | 3/3/2020      | LTI          | Orbera 365 Intra gastric Balloon System | Adverse Event Without Identified Device or Use Problem       | Nausea; Vomiting; Patient Problem/ Medical Problem |
| 30067221<br>12-2020-00029 | 12/23/2019 | Malfunction | Apollo Endosurgery, Inc. | 3/3/2020      | LTI          | Orbera® Intra gastric Balloon System    | Patient-Device Incompatibility                               | Pain; Vomiting                                     |
| 30067221<br>12-2020-00028 | 11/5/2019  | Malfunction | Apollo Endosurgery, Inc. | 2/24/2020     | LTI          | Orbera Intra gastric Balloon System     | Fluid/ Blood Leak  | Failure of Implant                                 |
| 30067221<br>12-2020-00025 | 11/5/2019  | Malfunction | Apollo Endosurgery, Inc. | 2/24/2020     | LTI          | Orbera Intra gastric Balloon System     | Fluid/ Blood Leak  | Failure of Implant; Nausea; Vomiting               |
| 30067221<br>12-2020-00023 | 12/16/2019 | Malfunction | Apollo Endosurgery, Inc. | 2/24/2020     | LTI          | Orbera Intra gastric Balloon System     | Fungus in Device Environment; Patient-Device Incompatibility | Pain; Vomiting                                     |
| 30067221<br>12-2020-00027 | 12/17/2019 | Malfunction | Apollo Endosurgery, Inc. | 2/24/2020     | LTI          | Orbera Intra gastric Balloon System     | Fluid/ Blood Leak  | Failure of Implant                                 |

| Report Number             | Event Date | Event Type  | Manufacturer                          | Date Received | Product Code | Brand Name                             | Device Problem   | Patient Problem  |
|---------------------------|------------|-------------|---------------------------------------|---------------|--------------|--|--|--|
| 30067221<br>12-2020-00024 | 12/27/2019 | Malfunction | Apollo Endosurgery, Inc.              | 2/24/2020     | LTI          | Orbera Intragastric Balloon System     | Fluid/ Blood Leak                                      | Failure of Implant   |
| 30067221<br>12-2020-00018 | 1/25/2020  | Injury      | Apollo Endosurgery, Inc.              | 2/24/2020     | LTI          | Orbera® Intragastric Balloon System    | Patient-Device Incompatibility                         | Aspiration/ Inhalation; Ischemia; Necrosis; Pain; Vomiting |
| 30067221<br>12-2020-00022 | 1/28/2020  | Injury      | Apollo Endosurgery, Inc.              | 2/24/2020     | LTI          | Bib; System Intragastric Balloon       | Detachment of Device or Device Component               | Perforation of Esophagus                                   |
| 30067221<br>12-2020-00026 | 12/7/2019  | Malfunction | Apollo Endosurgery, Inc.              | 2/24/2020     | LTI          | Orbera® Intragastric Balloon System    | Patient-Device Incompatibility                         | Pain; Vomiting   |
| MW5093<br>139             | 2/4/2020   | Injury      | Appollo Endosurgery Costa Rica S.R.L. | 2/19/2020     | LTI          | Orbera                                 | Adverse Event Without Identified Device or Use Problem | Perforation  |
| 30067221<br>12-2020-00019 | 1/24/2020  | Injury      | Apollo Endosurgery, Inc.              | 2/10/2020     | LTI          | Orbera365 Intragastric Balloon System  | Patient-Device Incompatibility                         | Nausea; Pain; Vomiting                                     |
| 30067221<br>12-2020-00015 | 12/14/2019 | Malfunction | Apollo Endosurgery, Inc.              | 2/3/2020      | LTI          | Orbera365® Intragastric Balloon System | Fluid/ Blood Leak                                      | Failure of Implant   |
| 30067221<br>12-2020-00016 | 12/15/2019 | Malfunction | Apollo Endosurgery, Inc.              | 2/3/2020      | LTI          | Orbera365 Intragastric Balloon System  | Fluid/ Blood Leak                                      | Failure of Implant   |
| 30067221<br>12-2020-00013 | 1/3/2020   | Malfunction | Apollo Endosurgery, Inc.              | 2/3/2020      | LTI          | Orbera365® Intragastric Balloon System | Fluid/ Blood Leak                                      | Failure of Implant   |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                              | Device Problem   | Patient Problem                                     |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|--|---|
| 30067221<br>12-2020-00014 | 1/7/2020   | Malfunction | Apollo Endosurgery, Inc. | 2/3/2020      | LTI          | Orbera365® Intra gastric Balloon System | Fluid/ Blood Leak  | Failure of Implant                                  |
| 30067221<br>12-2020-00009 | 1/17/2020  | Malfunction | Apollo Endosurgery, Inc. | 1/31/2020     | LTI          | Orbera365® Intra gastric Balloon System | Fluid/ Blood Leak; Migration or Expulsion of Device                                    | Failure of Implant                                  |
| 30067221<br>12-2020-00012 | 10/9/2019  | Injury      | Apollo Endosurgery, Inc. | 1/27/2020     | LTI          | Orbera365 Intra gastric Balloon System  | Patient-Device Incompatibility   | Dehydration; Vomiting                               |
| 30067221<br>12-2020-00010 | 11/14/2019 | Malfunction | Apollo Endosurgery, Inc. | 1/22/2020     | LTI          | Orbera365® Intra gastric Balloon System | Patient-Device Incompatibility   | Inflammation; Vomiting                              |
| 30067221<br>12-2019-00195 | 10/6/2019  | Injury      | Apollo Endosurgery, Inc. | 1/13/2020     | LTI          | Orbera Intra gastric Balloon System Us  | Patient-Device Incompatibility; Adverse Event Without Identified Device or Use Problem | Gastritis; Inflammation; Nausea; Vomiting; Weakness |
| 30067221<br>12-2020-00006 | 11/30/2019 | Malfunction | Apollo Endosurgery, Inc. | 1/13/2020     | LTI          | Orbera Intra gastric Balloon System     | Patient-Device Incompatibility   | Pain; Vomiting                                      |
| 30067221<br>12-2020-00005 | 12/13/2019 | Injury      | Apollo Endosurgery, Inc. | 1/13/2020     | LTI          | Orbera Intra gastric Balloon System     | Patient-Device Incompatibility   | Vomiting  |
| 30067221<br>12-2020-00003 | 12/2/2019  | Malfunction | Apollo Endosurgery, Inc. | 1/9/2020      | LTI          | Orbera365® Intra gastric Balloon System | Fluid/ Blood Leak  | No Consequences Or Impact To Patient                |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                              | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|--|---|
| 30067221<br>12-2020-00001 | NR         | Malfunction | Apollo Endosurgery, Inc. | 1/6/2020      | LTI          | Orbera365 Intra gastric Balloon System  | Patient-Device Incompatibility                         | Pain; Vomiting  |
| 30067221<br>12-2019-00205 | 11/26/2019 | Malfunction | Apollo Endosurgery, Inc. | 1/3/2020      | LTI          | Orbera® Intra gastric Balloon System    | Fluid/ Blood Leak; Leak/ Splash                        | Failure of Implant  |
| 30067221<br>12-2019-00182 | 9/21/2019  | Malfunction | Apollo Endosurgery, Inc. | 12/16/2019    | LTI          | Orbera® Intra gastric Balloon System    | Patient-Device Incompatibility                         | Pain; Vomiting  |
| 30067221<br>12-2019-00185 | 10/5/2019  | Malfunction | Apollo Endosurgery, Inc. | 12/16/2019    | LTI          | Orbera® Intra gastric Balloon System    | Fluid/ Blood Leak; Leak/ Splash                        | Failure of Implant  |
| 30067221<br>12-2019-00198 | 11/12/2019 | Malfunction | Apollo Endosurgery, Inc. | 12/16/2019    | LTI          | BIB System Intra gastric Balloon        | Material Split, Cut or Torn                            | No Code Available   |
| 30067221<br>12-2019-00203 | 11/24/2019 | Malfunction | Apollo Endosurgery, Inc. | 12/16/2019    | LTI          | Orbera 365 Intra gastric Balloon System | Leak/ Splash   | Discharge; Patient Problem/ Medical Problem                     |
| 30067221<br>12-2019-00178 | NR         | Injury      | Apollo Endosurgery, Inc. | 12/16/2019    | LTI          | Orbera® Intra gastric Balloon System    | Adverse Event Without Identified Device or Use Problem | Dehydration; Nausea; Vomiting; Patient Problem/ Medical Problem |
| 30067221<br>12-2019-00174 | NR         | Malfunction | Apollo Endosurgery, Inc. | 12/16/2019    | LTI          | Orbera365 Intra gastric Balloon System  | Fluid/ Blood Leak; Leak/ Splash                        | No Code Available   |
| 30067221<br>12-2019-00190 | 11/7/2019  | Malfunction | Apollo Endosurgery, Inc. | 12/10/2019    | LTI          | BIB System Intra gastric Balloon        | Detachment of Device or Device Component               | No Code Available   |



| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                            | Device Problem                                 | Patient Problem                     |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|---------------------------------------|--|-------------------------------------|
| 30067221<br>12-2019-<br>00191 | 11/7/2019  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 12/10/2019    | LTI          | BIB System<br>Intragastric<br>Balloon | Detachment of<br>Device or Device<br>Component | No Code Available                   |
| 30067221<br>12-2019-<br>00121 | 4/19/2019  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 5/22/2019     | LTI          | BIB System<br>Intragastric<br>Balloon | Patient-Device<br>Incompatibility              | Vomiting                            |
| 30067221<br>12-2019-<br>00062 | 9/3/2018   | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 2/22/2019     | LTI          | BIB System<br>Intragastric<br>Balloon | Patient-Device<br>Incompatibility              | Patient Problem/<br>Medical Problem |

Table M3. Reports on Endoscopic Sleeve Gastroplasty Devices from the FDA MAUDE Database<sup>2</sup>

Note: Key information is presented in this table, but detailed event descriptions are available at the [FDA MAUDE database](#).

| Report Number                 | Event Date     | Event Type  | Manufacturer                                | Date Received | Product Code | Brand Name                                   | Device Problem   | Patient Problem  |
|-------------------------------|----------------|-------------|---|---------------|--------------|--|--|--|
| 30067221<br>12-2021-<br>00059 | NR             | Malfunction | Apollo<br>Endosurgery<br>Costa Rica,<br>Srl | 6/18/2021     | OCW          | Endoscopic<br>Suturing<br>System             | Adverse Event<br>Without Identified<br>Device or Use<br>Problem; Insufficient<br>Information             | Hematoma;<br>Nausea; Vomiting  |
| 30067221<br>12-2021-<br>00060 | NR             | Injury      | Apollo<br>Endosurgery<br>Costa Rica,<br>Srl | 6/18/2021     | OCW          | Endoscopic<br>Suturing<br>System             | Adverse Event<br>Without Identified<br>Device or Use<br>Problem; Insufficient<br>Information             | Abdominal Pain;<br>Hematoma;<br>Perforation of<br>Esophagus;<br>Pancreatitis;<br>Stomach<br>Ulceration |
| 30067221<br>12-2022-<br>00033 | NR             | Malfunction | Apollo<br>Endosurgery,<br>Inc.              | 3/28/2022     | OCW          | ESS<br>Endoscopic<br>Suturing<br>System      | Insufficient<br>Information  | Hemorrhage/<br>Bleeding;<br>Myocardial<br>Infarction   |
| 30050998<br>03-2024-<br>00031 | 12/11/202<br>3 | Malfunction | Apollo<br>Endosurgery,<br>Inc               | 1/10/2024     | OCW          | Overstitch<br>Endoscopic<br>Suture<br>System | Break; Device-<br>Device<br>Incompatibility;<br>Difficult to Open or<br>Close; Material<br>Twisted/ Bent | No Clinical Signs,<br>Symptoms or<br>Conditions  |
| 30050998<br>03-2023-<br>07079 | 12/14/202<br>3 | Malfunction | Apollo<br>Endosurgery<br>Inc.,              | 1/9/2024      | OCW          | Overstitch<br>Endoscopic<br>Suture<br>System | Failure to Align   | No Clinical Signs,<br>Symptoms or<br>Conditions  |
| 30050998<br>03-2023-<br>07021 | 11/23/202<br>3 | Malfunction | Apollo<br>Endosurgery                       | 1/3/2024      | OCW          | Overstitch<br>Endoscopic<br>Suture<br>System | Failure to Align;<br>Difficult to Open or<br>Close   | No Clinical Signs,<br>Symptoms or<br>Conditions  |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                            | Device Problem                                      | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---------------------------------------|---|---|
| 30067221<br>12-2023-00159 | 7/3/2023   | Malfunction | Apollo Endosurgery, Inc. | 8/25/2023     | OCW          | Overstitch Endoscopic Suture System   | No Apparent Adverse Event; Insufficient Information | Insufficient Information; Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 17224403                  | 2/8/2023   | Malfunction | Apollo Endosurgery, Inc. | 6/29/2023     | HCF          | Overstitch Endoscopic Suturing System | Mechanical Problem; Failure to Cut                  | Tissue Breakdown  |
| 17001751                  | 5/10/2023  | Malfunction | Apollo Endosurgery, Inc. | 5/25/2023     | HCF          | Overstitch Endoscopic Suturing System | Activation, Positioning or Separation Problem       | No Clinical Signs, Symptoms or Conditions   |
| 14134966                  | 3/15/2022  | Malfunction | Apollo Endosurgery, Inc. | 4/18/2022     | OCW          | Overstitch Endoscopic Suturing System | Failure to Align                                    | No Clinical Signs, Symptoms or Conditions   |
| 13652741                  | 1/11/2022  | Malfunction | Apollo Endosurgery, Inc. | 3/2/2022      | OCW          | Overstitch Endoscopic Suturing System | Material Twisted/ Bent                              | Insufficient Information  |
| 30067221<br>12-2021-00102 | NR         | Malfunction | Apollo Endosurgery, Inc. | 10/26/2021    | OCW          | Overstitch Endoscopic Suturing System | Insufficient Information                            | Dysphagia/ Odynophagia  |
| 30067221<br>12-2021-00100 | NR         | Malfunction | Apollo Endosurgery, Inc. | 10/22/2021    | OCW          | Overstitch Endoscopic Suturing System | Insufficient Information                            | Abdominal Pain; Nausea; Vomiting  |
| 30067221<br>12-2021-00101 | NR         | Malfunction | Apollo Endosurgery, Inc. | 10/22/2021    | OCW          | Overstitch Endoscopic Suturing System | Mechanical Problem                                  | Perforation   |

| Report Number                 | Event Date     | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                     | Device Problem                                  | Patient Problem  |
|-------------------------------|----------------|-------------|--------------------------------|---------------|--------------|--|---|--|
| 30067221<br>12-2021-<br>00087 | NR             | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 9/20/2021     | OCW          | Overstitch<br>Endoscopic<br>Suturing<br>System | Insufficient<br>Information                     | Abdominal Pain;<br>Fever; Vomiting;<br>Gastrointestinal<br>Hemorrhage                    |
| 30067221<br>12-2021-<br>00071 | 6/8/2021       | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 8/4/2021      | OCW          | Overstitch<br>Endoscopic<br>Suturing<br>System | Appropriate Term/<br>Code Not Available         | Pneumonia;<br>Respiratory Tract<br>Infection   |
| 30067221<br>12-2021-<br>00057 | 6/5/2021       | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 7/6/2021      | OCW          | Overstitch<br>Endoscopic<br>Suturing<br>System | Insufficient<br>Information                     | Failure of Implant   |
| 30067221<br>12-2021-<br>00044 | 5/12/2021      | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/29/2021     | OCW          | Overstitch<br>Endoscopic<br>Suturing<br>System | Mechanical Problem;<br>Mechanical Jam           | Appropriate<br>Clinical Signs,<br>Symptoms,<br>Conditions Term/<br>Code Not<br>Available |
| 30067221<br>12-2021-<br>00052 | 5/25/2021      | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/23/2021     | LTI          | Overstitch<br>Endoscopic<br>Suturing<br>System | Mechanical Jam                                  | Perforation  |
| 30067221<br>12-2021-<br>00045 | NR             | Injury      | Apollo<br>Endosurgery,<br>Inc. | 6/18/2021     | OCW          | Overstitch<br>Endoscopic<br>Suturing<br>System | Insufficient<br>Information                     | Bowel Perforation  |
| 30067221<br>12-2021-<br>00058 | 4/13/2021      | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 6/14/2021     | OCW          | Overstitch<br>Endoscopic<br>Suturing<br>System | Entrapment of<br>Device; Difficult to<br>Remove | Device Embedded<br>In Tissue or Plaque   |
| 30067221<br>12-2020-<br>00113 | 10/19/202<br>0 | Injury      | Apollo<br>Endosurgery          | 11/19/2020    | OCW          | Overstitch<br>Endoscopic<br>Suturing<br>System | Entrapment of<br>Device; Difficult to<br>Remove | Device Embedded<br>In Tissue or Plaque   |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name  | Device Problem  | Patient Problem                                   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|---|---|
| 30067221<br>12-2019-00172 | 9/16/2019  | Malfunction | Apollo Endosurgery, Inc. | 1/3/2020      | OCW          | Overstitch Endoscopic Suturing System                                   | Entrapment of Device; Detachment of Device or Device Component; Mechanical Jam; Positioning Problem | Suture Abrasion; Patient Problem/ Medical Problem |
| 30067221<br>12-2019-00197 | 10/1/2019  | Malfunction | Apollo Endosurgery, Inc. | 12/16/2019    | OCW          | Overstitch Endoscopic Suturing System                                   | Entrapment of Device; Difficult to Remove; Mechanical Jam   | Suture Abrasion; Patient Problem/ Medical Problem |
| 12165605                  | 7/9/2021   | Malfunction | Apollo Endosurgery, Inc. | 7/14/2021     | OCW          | Overstitch Endoscopic Suturing System Needle Driver And Anchor Exchange | Physical Resistance/ Sticking   | Insufficient Information                          |
| 30050998<br>03-2023-07140 | 12/19/2023 | Malfunction | Apollo Endosurgery       | 1/10/2024     | OCW          | Overstitch Sx Endoscopic Suture System                                  | Entrapment of Device; Detachment of Device or Device Component; Device-Device Incompatibility       | No Clinical Signs, Symptoms or Conditions         |
| 30050998<br>03-2023-07020 | 12/6/2023  | Malfunction | Apollo Endosurgery, Inc  | 1/4/2024      | OCW          | Overstitch Sx Endoscopic Suture System                                  | Entrapment of Device  | No Clinical Signs, Symptoms or Conditions         |
| 30050998<br>03-2023-06588 | 10/3/2023  | Malfunction | Apollo Endosurgery, Inc. | 12/7/2023     | OCW          | Overstitch Sx Endoscopic Suture System                                  | Device-Device Incompatibility; Difficult to Open or Close   | No Clinical Signs, Symptoms or Conditions         |
| 30067221<br>12-2023-00061 | 3/2/2023   | Malfunction | Apollo Endosurgery, Inc  | 3/31/2023     | OCW          | Overstitch Sx Endoscopic Suturing System                                | Insufficient Information  | Perforation of Esophagus                          |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name  | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---|--|---|
| 30067221<br>12-2023-00129 | 3/1/2023   | Malfunction | Apollo Endosurgery, Inc  | 5/26/2023     | OCW          | Overstitch Sx <sub>2</sub> Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem             | Failure of Implant; Liver Damage/ Dysfunction; Pancreatitis; Unspecified Kidney or Urinary Problem  |
| 30067221<br>12-2023-00217 | 10/11/2023 | Malfunction | Apollo Endosurgery, Inc. | 11/3/2023     | OCW          | Overstitch Endoscopic Suture System                   | Material Integrity Problem   | Insufficient Information; Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2023-00211 | 10/8/2023  | Injury      | Apollo Endosurgery, Inc. | 11/2/2023     | OCW          | Overstitch Endoscopic Suture System                   | Insufficient Information   | Failure of Implant; Pneumonia; Sepsis   |
| 30067221<br>12-2023-00195 | 6/29/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/20/2023     | OCW          | Overstitch Endoscopic Suture System                   | Difficult to Remove; Activation, Positioning or Separation Problem | Insufficient Information; Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2023-00179 | 12/12/2022 | Malfunction | Apollo Endosurgery, Inc. | 9/15/2023     | OCW          | Overstitch Endoscopic Suture System                   | Material Perforation; Mechanical Jam                               | Laceration(s) of Esophagus  |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                          | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|-------------------------------------|--|---|
| 30067221<br>12-2023-00125 | 4/28/2023  | Malfunction | Apollo Endosurgery, Inc. | 6/12/2023     | OCW          | Overstitch Endoscopic Suture System | Mechanical Problem                                     | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2023-00086 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/15/2023     | OCW          | Overstitch Endoscopic Suture System | Adverse Event Without Identified Device or Use Problem | Perforation of Vessels; Bowel Perforation                                 |
| 30067221<br>12-2023-00075 | NR         | Malfunction | Apollo Endosurgery, Inc  | 4/12/2023     | OCW          | Overstitch Endoscopic Suture System | Migration  | Perforation   |
| 30067221<br>12-2023-00076 | NR         | Malfunction | Apollo Endosurgery, Inc  | 4/10/2023     | OCW          | Overstitch Endoscopic Suture System | Migration  | Peritonitis   |
| 30067221<br>12-2023-00072 | NR         | Malfunction | Apollo Endosurgery, Inc  | 3/31/2023     | OCW          | Overstitch Endoscopic Suture System | Migration or Expulsion of Device; Migration            | Syncope/ Fainting; Gastrointestinal Hemorrhage                            |
| 30067221<br>12-2023-00056 | NR         | Malfunction | Apollo Endosurgery, Inc  | 3/15/2023     | OCW          | Overstitch Endoscopic Suture System | No Apparent Adverse Event                              | Gastrointestinal Hemorrhage   |
| 30067221<br>12-2023-00008 | 10/25/2022 | Malfunction | Apollo Endosurgery, Inc  | 2/2/2023      | OCW          | Overstitch Endoscopic Suture System | No Apparent Adverse Event                              | Hemorrhage/ Bleeding; Vomiting  |
| 30067221<br>12-2022-00138 | 9/12/2022  | Death       | Apollo Endosurgery, Inc  | 1/13/2023     | OCW          | Overstitch Endoscopic Suture System | Adverse Event Without Identified Device or Use Problem | Pulmonary Embolism; Cardiac Arrest  |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                            | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---------------------------------------|--|---|
| 30067221<br>12-2023-00206 | 9/24/2023  | Malfunction | Apollo Endosurgery, Inc  | 10/30/2023    | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem | Cardiac Arrest  |
| 30067221<br>12-2023-00197 | NR         | Malfunction | Apollo Endosurgery, Inc  | 10/9/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem | Hemorrhage/<br>Bleeding   |
| 30067221<br>12-2023-00171 | 8/14/2023  | Malfunction | Apollo Endosurgery, Inc. | 9/8/2023      | OCW          | Overstitch Endoscopic Suturing System | Premature Separation                                   | Failure of Implant  |
| 30067221<br>12-2023-00154 | 7/18/2023  | Malfunction | Apollo Endosurgery, Inc. | 8/15/2023     | OCW          | Overstitch Endoscopic Suturing System | Entrapment of Device; Difficult to Remove              | Device Embedded In Tissue or Plaque                             |
| 30067221<br>12-2023-00145 | 2/8/2023   | Malfunction | Apollo Endosurgery, Inc  | 7/21/2023     | OCW          | Overstitch Endoscopic Suturing System | Material Perforation                                   | Laceration(s) of Esophagus                                      |
| 30067221<br>12-2023-00134 | 6/12/2022  | Malfunction | Apollo Endosurgery, Inc  | 6/27/2023     | OCW          | Overstitch Endoscopic Suturing System | Migration  | Abdominal Pain;<br>Diarrhea;<br>Perforation;<br>Vomiting        |
| 30067221<br>12-2023-00138 | NR         | Malfunction | Apollo Endosurgery, Inc  | 6/27/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem | Hemorrhage/<br>Bleeding; Ulcer                                  |
| 30067221<br>12-2023-00135 | NR         | Malfunction | Apollo Endosurgery, Inc  | 6/21/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem | Pulmonary Embolism;<br>Hemorrhage/<br>Bleeding; Nausea;<br>Pain |



| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                            | Device Problem  | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---------------------------------------|---|---|
| 30067221<br>12-2023-00137 | NR         | Malfunction | Apollo Endosurgery, Inc  | 6/21/2023     | OCW          | Overstitch Endoscopic Suturing System | Fluid/ Blood Leak   | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2023-00133 | NR         | Malfunction | Apollo Endosurgery, Inc  | 6/21/2023     | OCW          | Overstitch Endoscopic Suturing System | Migration   | Hemorrhage/ Bleeding; Pain; Ulcer; Pancreatitis                           |
| 30067221<br>12-2023-00113 | 5/5/2023   | Malfunction | Apollo Endosurgery, Inc. | 6/5/2023      | OCW          | Overstitch Endoscopic Suturing System | Difficult to Fold, Unfold or Collapse; Unintended Collision                 | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2023-00110 | 4/19/2023  | Malfunction | Apollo Endosurgery, Inc  | 6/5/2023      | OCW          | Overstitch Endoscopic Suturing System | Insufficient Information; Loosening of Implant Not Related to Bone-Ingrowth | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2023-00127 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/31/2023     | OCW          | Overstitch Endoscopic Suturing System | Insufficient Information  | Hemorrhage/ Bleeding  |
| 30067221<br>12-2023-00123 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/30/2023     | OCW          | Overstitch Endoscopic Suturing System | Appropriate Term/ Code Not Available  | Chest Pain; Fever; Hemorrhage/ Bleeding; Implant Pain                     |
| 30067221<br>12-2023-00122 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/30/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem                      | Hemorrhage/ Bleeding; Liver Damage/ Dysfunction; Pain                     |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                            | Device Problem  | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---------------------------------------|---|---|
| 30067221<br>12-2023-00106 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/23/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem            | Fistula; Hemorrhage/ Bleeding; Perforation; Ulcer                         |
| 30067221<br>12-2023-00090 | 4/24/2023  | Malfunction | Apollo Endosurgery, Inc. | 5/19/2023     | OCW          | Overstitch Endoscopic Suturing System | Mechanical Problem; Mechanical Jam                                | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2023-00100 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/18/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem            | Pulmonary Embolism; Hemorrhage/ Bleeding; Perforation                     |
| 30067221<br>12-2023-00104 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/18/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem            | Hemorrhage/ Bleeding; Pancreatitis  |
| 30067221<br>12-2023-00101 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/18/2023     | OCW          | Overstitch Endoscopic Suturing System | Mechanical Problem  | Hemorrhage/ Bleeding; Perforation of Esophagus                            |
| 30067221<br>12-2023-00121 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/18/2023     | OCW          | Overstitch Endoscopic Suturing System | Insufficient Information  | Fistula; Sepsis   |
| 30067221<br>12-2023-00098 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/17/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem            | Pulmonary Embolism; Dehydration; Thrombosis/ Thrombus                     |
| 30067221<br>12-2023-00099 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/17/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem; Migration | Abdominal Pain; Hemorrhage/ Bleeding; Ulcer                               |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                            | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---------------------------------------|--|---|
| 30067221<br>12-2023-00097 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/16/2023     | OCW          | Overstitch Endoscopic Suturing System | Partial Blockage                                       | Obstruction/ Occlusion;<br>Thrombosis/<br>Thrombus  |
| 30067221<br>12-2023-00096 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/16/2023     | OCW          | Overstitch Endoscopic Suturing System | Leak/ Splash   | Hemorrhage/<br>Bleeding   |
| 30067221<br>12-2023-00095 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/16/2023     | OCW          | Overstitch Endoscopic Suturing System | Migration  | Vomiting;<br>Obstruction/<br>Occlusion  |
| 30067221<br>12-2023-00094 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/16/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem | Abdominal Pain;<br>Perforation of Esophagus   |
| 30067221<br>12-2023-00087 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/15/2023     | OCW          | Overstitch Endoscopic Suturing System | Material Perforation;<br>Migration                     | Abdominal Pain;<br>Perforation of Vessels;<br>Bowel Perforation   |
| 30067221<br>12-2023-00082 | 4/11/2023  | Malfunction | Apollo Endosurgery, Inc. | 5/12/2023     | OCW          | Overstitch Endoscopic Suturing System | No Apparent Adverse Event;<br>Insufficient Information | Insufficient Information;<br>Appropriate Clinical Signs,<br>Symptoms,<br>Conditions Term/<br>Code Not Available |
| 30067221<br>12-2023-00088 | NR         | Malfunction | Apollo Endosurgery, Inc  | 5/12/2023     | OCW          | Overstitch Endoscopic Suturing System | Patient Device Interaction Problem                     | Perforation of Esophagus  |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                            | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---------------------------------------|--|---|
| 30067221<br>12-2023-00039 | NR         | Malfunction | Apollo Endosurgery, Inc  | 2/21/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem | Perforation; Perforation of Vessels; Perforation of Esophagus             |
| 30067221<br>12-2023-00027 | 12/15/2022 | Malfunction | Apollo Endosurgery, Inc  | 2/14/2023     | OCW          | Overstitch Endoscopic Suturing System | Difficult to Insert; Difficult to Advance              | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2023-00004 | 1/4/2023   | Malfunction | Apollo Endosurgery, Inc  | 1/31/2023     | OCW          | Overstitch Endoscopic Suturing System | Material Perforation                                   | Hemorrhage/ Bleeding  |
| 30067221<br>12-2023-00010 | NR         | Malfunction | Apollo Endosurgery, Inc  | 1/25/2023     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem | Post Operative Wound Infection  |
| 30067221<br>12-2022-00088 | 8/12/2022  | Malfunction | Apollo Endosurgery, Inc. | 9/6/2022      | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem | Gastrointestinal Hemorrhage   |
| 30067221<br>12-2022-00065 | 4/25/2022  | Malfunction | Apollo Endosurgery, Inc. | 6/21/2022     | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem | Pain; Malaise; Cramp(s)/ Muscle Spasm(s)                                  |
| 30067221<br>12-2022-00047 | NR         | Malfunction | Apollo Endosurgery, Inc. | 5/27/2022     | OCW          | Overstitch Endoscopic Suturing System | Appropriate Term/ Code Not Available                   | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |

| Report Number             | Event Date | Event Type  | Manufacturer             | Date Received | Product Code | Brand Name                            | Device Problem   | Patient Problem   |
|---------------------------|------------|-------------|--------------------------|---------------|--------------|---------------------------------------|--|---|
| 30067221<br>12-2022-00042 | 4/5/2022   | Malfunction | Apollo Endosurgery, Inc. | 5/7/2022      | OCW          | Overstitch Endoscopic Suturing System | Mechanical Problem   | Appropriate Clinical Signs, Symptoms, Conditions Term/ Code Not Available |
| 30067221<br>12-2022-00030 | 2/16/2022  | Malfunction | Apollo Endosurgery, Inc. | 3/17/2022     | OCW          | Overstitch Endoscopic Suturing System | Mechanical Problem; Detachment of Device or Device Component | Perforation   |
| 30067221<br>12-2022-00027 | 2/10/2022  | Malfunction | Apollo Endosurgery, Inc. | 3/8/2022      | OCW          | Overstitch Endoscopic Suturing System | Mechanical Problem; Mechanical Jam                           | Perforation   |
| 30067221<br>12-2021-00116 | NR         | Malfunction | Apollo Endosurgery, Inc. | 12/8/2021     | OCW          | Overstitch Endoscopic Suturing System | Insufficient Information                                     | Nausea; Vomiting; Gastrointestinal Hemorrhage                             |
| 30067221<br>12-2021-00114 | NR         | Malfunction | Apollo Endosurgery, Inc. | 11/17/2021    | OCW          | Overstitch Endoscopic Suturing System | Adverse Event Without Identified Device or Use Problem       | Hematoma; Genital Bleeding  |
| 30067221<br>12-2020-00007 | 1/9/2020   | Malfunction | Apollo Endosurgery, Inc. | 1/20/2020     | OCW          | Overstitch Endoscopic Suturing System | Difficult to Remove  | Patient Problem/ Medical Problem  |
| 30067221<br>12-2020-00002 | 11/19/2019 | Malfunction | Apollo Endosurgery, Inc. | 1/7/2020      | OCW          | Overstitch Endoscopic Suturing System | Difficult to Remove  | Patient Problem/ Medical Problem  |
| 30067221<br>12-2019-00176 | 9/25/2019  | Malfunction | Apollo Endosurgery, Inc. | 12/16/2019    | OCW          | Overstitch Endoscopic Suturing System | Material Integrity Problem; Material Twisted/ Bent           | No Code Available   |

| Report Number                 | Event Date | Event Type  | Manufacturer                   | Date Received | Product Code | Brand Name                                     | Device Problem              | Patient Problem               |
|-------------------------------|------------|-------------|--------------------------------|---------------|--------------|--|-----------------------------|-------------------------------|
| 30067221<br>12-2019-<br>00193 | 11/6/2019  | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 12/10/2019    | OCW          | Overstitch<br>Endoscopic<br>Suturing<br>System | Insufficient<br>Information | No Code Available             |
| 30067221<br>12-2022-<br>00083 | NR         | Malfunction | Apollo<br>Endosurgery,<br>Inc. | 8/30/2022     | OCW          | Overstitch<br>Endoscopic<br>Suture<br>System   | Insufficient<br>Information | Hemorrhage/<br>Bleeding; Pain |

## Appendix M References

1. US Food and Drug Administration. Medical Device Recalls. 2024; <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfRES/res.cfm>. Accessed February 2023.
2. US Food and Drug Administration. MAUDE - Manufacturer and User Facility Device Experience. 2024; <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfMAUDE/search.CFM>. Accessed February 27, 2023.